

#### US010023921B2

# (12) United States Patent

Barbour et al.

# (54) CORN EVENT TC1507 AND METHODS FOR DETECTION THEREOF

(71) Applicants: Pioneer Hi-Bred International, Inc.,
Johnston, IA (US); Dow Agrosciences
LLC, Indianapolis, IN (US); E. I. Du
Pont De Nemours and Company,
Wilmington, DE (US)

Inventors: Eric Barbour, Johnston, IA (US); James Wayne Bing, Ankeny, IA (US); Guy A. Cardineau, Tempe, AZ (US); Robert F. Cressman, Jr., Wilmington, DE (US); Manju Gupta, Carmel, IN (US); Mary E. Hartnett Locke, Mickleton, NJ (US); David Hondred, Ankeny, IA (US); Joseph W. Keaschall, Clive, IA (US); Michael G. Koziel, Raleigh, NC (US); Terry **EuClaire Meyer**, Urbandale, IA (US); Daniel Moellenbeck, Granger, IA (US); Kenneth Edwin Narva, Carlsbad, CA (US); Wilas Nirunsuksiri, Auburn, WA (US); Stephen W. Ritchie, Omaha, NE (US); Marjorie L. Rudert, Boone, IA (US); Craig D. Sanders, Bear, DE (US); Aihua Shao, Johnston, IA (US); Steven Jeffrey Stelman, San Diego, CA (US); **David S. Stucker**, Johnston, IA (US); Laura Ann Tagliani, Zionsville, IN (US); William M. Van Zante, Urbandale, IA (US)

(73) Assignees: PIONEER HI-BRED
INTERNATIONAL, INC., Johnston,
IA (US); E. I. DU PONT DE
NEMOURS AND COMPANY,
Wilmington, DE (US); DOW
AGROSCIENCES LLC, Indianapolis,
IN (US)

(\*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 703 days.

(21) Appl. No.: 14/547,589

(22) Filed: Nov. 19, 2014

(65) Prior Publication Data

US 2015/0072896 A1 Mar. 12, 2015

## Related U.S. Application Data

- (60) Division of application No. 12/914,271, filed on Oct. 28, 2010, now Pat. No. 8,901,378, which is a division of application No. 12/333,044, filed on Mar. 13, 2009, now Pat. No. 7,989,607, which is a continuation of application No. 11/774,236, filed on Jul. 6, 2007, now Pat. No. 7,514,544, which is a division of application No. 10/837,105, filed on Apr. 30, 2004, now Pat. No. 7,288,643.
- (60) Provisional application No. 60/467,772, filed on May 2, 2003.

## (10) Patent No.: US 10,023,921 B2

(45) **Date of Patent:** Jul. 17, 2018

(51) **Int. Cl.** 

 C12Q 1/68
 (2018.01)

 C12N 15/82
 (2006.01)

 C12Q 1/6895
 (2018.01)

 C12Q 1/686
 (2018.01)

(52) **U.S. Cl.** 

CPC ...... *C12Q 1/6895* (2013.01); *C12N 15/8277* (2013.01); *C12N 15/8286* (2013.01); *C12Q 1/686* (2013.01); *C12Q 2600/13* (2013.01); *C12Q 2600/16* (2013.01); *Y02A 40/162* (2018.01); *Y10S 435/975* (2013.01)

(58) Field of Classification Search

None

See application file for complete search history.

## (56) References Cited

#### U.S. PATENT DOCUMENTS

| 5,474,796    | A          | 12/1995 | Brennan        |
|--------------|------------|---------|----------------|
| 5,877,012    | A          | 3/1999  | Estruch et al. |
| 6,166,302    | A          | 12/2000 | Merlo et al.   |
| 6,646,184    | B2         | 11/2003 | Hohn et al.    |
| 6,713,259    | B2         | 3/2004  | Levine         |
| 6,825,400    | B2         | 11/2004 | Behr et al.    |
| 7,288,643    | B2         | 10/2007 | Barbour et al. |
| 7,417,132    | B2         | 8/2008  | Barbour et al. |
| 7,435,807    | B1         | 10/2008 | Barbour et al. |
| 7,449,564    | B2         | 11/2008 | Barbour et al. |
| 7,514,544    | B2         | 4/2009  | Barbour et al. |
| 7,989,607    | B2         | 8/2011  | Barbour et al. |
| 2001/0053519 | <b>A</b> 1 | 12/2001 | Fodor et al.   |
| 2002/0013960 | <b>A</b> 1 | 1/2002  | Behr et al.    |
|              |            | (Cont   | inued)         |

## FOREIGN PATENT DOCUMENTS

WO WO-2004/011601 A2 2/2004

## OTHER PUBLICATIONS

Genbank Accession No. S94464—polyubiquitin [maize, genomic, 3841 nt] (submitted May 23, 1993, retrieved on Aug. 24, 2017 from http://www.ncbi.nlm.nih.gov/nuccore/S94464).\*

(Continued)

Primary Examiner — Samuel C Woolwine

Assistant Examiner — Olayinka A Oyeyemi

(74) Attorney, Agent, or Firm — Ballard Spahr LLP

## (57) ABSTRACT

The invention provides DNA compositions that relate to transgenic insect resistant maize plants. Also provided are assays for detecting the presence of the maize TC1507 event based on the DNA sequence of the recombinant construct inserted into the maize genome and the DNA sequences flanking the insertion site. Kits and conditions useful in conducting the assays are provided.

## 6 Claims, 1 Drawing Sheet

## (56) References Cited

#### U.S. PATENT DOCUMENTS

2002/0102582 A1 8/2002 Levine 2002/0162136 A1 10/2002 Hohn et al. 2002/0177526 A1 11/2002 Chen et al.

## OTHER PUBLICATIONS

Genbank Accession No. AY562548—Binary Vector pORE\_R5, complete sequence (submitted Mar. 1, 2004, retrieved on Aug. 24, 2017 from http://www.ncbi.nlm.nih.gov/nuccore/AY562548).\*

Anonymous, "EM STD:D28338", available at http://ibis.internal.epo.orgjexamjdbfetch.jsp?id=EM STD%3AD28338, created Jan. 27, 1995 (1 page).

Written Opinion issued by the International Searching Authority dated Feb. 22, 2005 for international application PCT/US2004/013538, filed on Apr. 29, 2004 and published as WO 2004/099447 on Nov. 18, 2004 (Applicant—Dow AgroSciences, LLC) (Inventor—Barbour, et al.) (9 pages).

Anonymous, "EM\_STD:D28338", available at http://ibis.internal.epo.orgjexamjdbfetch.jsp?id=EM\_STD%3AD28338, created Jan. 27, 1995 (1 page).

Extended European Search Report issued by the European Patent Office dated Oct. 2, 2015 for application EP 15168051.9, filed on Apr. 29, 2004 and published as EP 2942402 on Nov. 11, 2015 (Applicant—Dow AgroSciences, LLC // Inventor—Hondred, et al.) (13 pages).

International Search Report issued by the International Searching Authority dated Feb. 22, 2005 for international application PCT/US2004/013538, filed on Apr. 29, 2004 and published as WO 2004/099447 on Nov. 18, 2004 (Applicant—Dow AgroSciences, LLC) (Inventor—Barbour, et al.) (4 pages).

Written Opinion issued by the International Searching Authority dated Feb. 22, 2005 for international application PCT/U52004/013538, filed on Apr. 29, 2004 and published as WO 2004/099447 on Nov. 18, 2004 (Applicant—Dow AgroSciences, LLC) (Inventor—Barbour, et al.) (9 pages).

International Preliminary Report on Patentability issued by the International Searching Authority dated Nov. 4, 2005 for international application PCT/US2004/013538, filed on Apr. 29, 2004 and published as WO 2004/099447 on Nov. 18, 2004 (Applicant—Dow AgroSciences, LLC) (Inventor—Barbour, et al.) (10 pages).

Federal Register Environmental Documents, Mycogen c/o Dow and Pioneer; Availability of Determination of Nonregulated Status for Corn Genetically Engineered for Insect Resistance and Glufosinate Herbicide Tolerance. U.S. Environmental Protection Agency. Aug. 14, 2001; Federal Register. 66(157):42624-5; Printout from http://www.aqbios.com/docroot/decdocs/02122001/pdf (12 pages).

Petition for Determination of Non-Regulated Status B.t. Cry1F Insect-Resistant, Glufosinate-Tolerant Maize Line 1507. May 10, 2000. Printout from http://www.aqbios.com/docroot/decdocs/04-225-004.pdf (94 pages).

Anonymous, "DAS-01507-1 (TC1507)", available at: http://www.agbios.com/dbase.php?action=showprod&data=TC1507

&frmat=long, created on Oct. 17, 2002.

Anonymous, "EM\_STD:D28338", available at http://ibis.internal.epo.orgjexamjdbfetch.jsp?id=EMSTD%3AD28338, created Jan. 27, 1995 (1 page).

Anonymous, "GSN: AAN80062", available at http://ibis.internal.epo.orgjexamjdbfetch.jsp?id=GSN%3AAAN80062, created on Jan. 21, 1987 (1 page).

Anonymous, "ig56e02.b1 WGS-ZmaysF (DH5a methyl filtered) Zea mays genomic clone ig—GSS—NCBI", available at http://www.ncbi.nlm.nih.govjnucgss/BZ616018, created on Sep. 23, 2015 (1 page).

Anonymous, "PUBAQ15TD ZM\_0.6\_1.0 KB Zea mays genomic clone ZMMBTa012C05, genomic s—GSS—NCBI", available at http://www.ncbi.nlm.nih.govjnucgss/BZ677673, created Feb. 5, 2003 (1 page).

Anonymous, "Skov-og Naturstyrelsen", available at: http://www.sns.dk/erhvogadm/hoerinaer/c\_n1\_00/10/indhold.htm, accessed on Jan. 6, 2005.

Cornejo, et al., (1993) "Activity of a maize ubiquitin promoter in transgenic rice", Plant Molecular Biology vol. 23, pp. 567-581. Database EMBL (2002), Mus musculus chromosome 5 clone rp23-426e16 strain C57BL/6J, complete sequence, retrieved from EBI, Database accession No. AC084071.

New England Biolabs 1998/99 Catalog (NEB Catalog) (1 page). Printout from http://www.aqbios.com/docroot/decdocs/02122001/pdf.

Printout from http://www.aqbios.com/docroot/decdocs/04-225-004.pdf.

Printout from http://www.epa.gov/fedrgstr/EPA-IMPACT/2001/Aug./Day-14/20307.htm.

Spencer, et al., (1992) "Segregation of transgenes in maize", Plant Molecular Biology vol. 18, pp. 201-210.

Stratagene ("Gene Characterization Kits" 1988) (1 page).


Extended European Search Report issued by the European Patent Office dated Oct. 2, 2015 for application EP 15168051.9, filed Apr. 29, 2004 and published as EP 2942402 dated Nov. 11, 2015 (Applicant—Dow AgroSciences, LLC // Inventor—Hondred, et al.) (13 pages).

International Search Report issued by the International Searching Authority dated Feb. 22, 2005 for international application PCT/US2004/013538, filed Apr. 29, 2004 and published as WO 2004/099447 dated Nov. 18, 2004 (Applicant—Dow AgroSciences, LLC) (Inventor—Barbour, et al.) (4 pages).

Written Opinion issued by the International Searching Authority dated Feb. 22, 2005 for international application PCT/US2004/013538, filed Apr. 29, 2004 and published as WO 2004/099447 dated Nov. 18, 2004 (Applicant—Dow AgroSciences, LLC) (Inventor—Barbour, et al.) (9 pages).

International Preliminary Report on Patentability issued by the International Searching Authority dated Nov. 4, 2005 for international application PCT/US2004/013538, filed Apr. 29, 2004 and published as WO 2004/099447 dated Nov. 18, 2004 (Applicant—Dow AgroSciences, LLC) (Inventor—Barbour, et al.) (10 pages).

<sup>\*</sup> cited by examiner



## CORN EVENT TC1507 AND METHODS FOR **DETECTION THEREOF**

## CROSS-REFERENCE TO RELATED APPLICATIONS

This application is a divisional of U.S. Ser. No. 12/914, 271, filed Oct. 28, 2010, which is a divisional of U.S. Ser. No. 12/333,044, filed Mar. 13, 2009 (now U.S. Pat. No. 7,989,607, issued Aug. 2, 2011), which is a continuation of  $^{10}$ U.S. Ser. No. 11/774,236, filed Jul. 6, 2007 (now U.S. Pat. No. 7,514,544, issued Apr. 7, 2009), which is a divisional of U.S. Ser. No. 10/837,105, filed Apr. 30, 2004 (now U.S. Pat. No. 7,288,643, issued Oct. 30, 2007), which claims the benefit of U.S. Provisional Application Ser. No. 60/467,772, filed May 2, 2003, the contents of all of which are herein incorporated by reference in their entirety.

## REFERENCE TO SEQUENCE LISTING

The Sequence Listing submitted on Nov. 17, 2014 as a text file named "36446\_0046U2\_Sequence Listing.txt," created on Nov. 17, 2014, and having a size of 57,344 bytes is hereby incorporated by reference pursuant to 37 C.F.R.  $\S 1.52(e)(5)$ .

## FIELD OF THE INVENTION

The present invention relates to the field of plant molecular biology, specifically the invention relates to a DNA 30 construct for conferring insect resistance to a plant. The invention more specifically relates to an insect resistant corn plant TC1507 and to assays for detecting the presence of corn plant TC1507 DNA in a sample and compositions thereof.

## BACKGROUND OF THE INVENTION

This invention relates to the insect resistant corn (Zea mays) plant TC1507, also referred to as maize line TC1507 or maize event TC1507, and to the DNA plant expression construct of corn plant TC1507 and the detection of the transgene/flanking insertion region in corn plant TC1507 and progeny thereof.

Corn is an important crop and is a primary food source in 45 many areas of the world. Damage caused by insect pests is a major factor in the loss of the world's corn crops, despite the use of protective measures such as chemical pesticides. In view of this, insect resistance has been genetically engineered into crops such as corn in order to control insect 50 damage and to reduce the need for traditional chemical pesticides. One group of genes which have been utilized for the production of transgenic insect resistant crops are the delta-endotoxins from *Bacillus thuringiensis* (B.t.). Deltaendotoxins have been successfully expressed in crop plants 55 such as cotton, potatoes, rice, sunflower, as well as corn, and have proven to provide excellent control over insect pests. (Perlak, F. J et al. (1990) *Bio/Technology* 8, 939-943; Perlak, F. J. et al. (1993) *Plant Mol. Biol.* 22: 313-321; Fujimoto H. et al. (1993) Bio/Technology 11: 1151-1155; Tu et al. (2000) 60 ing and selecting an insect resistant monocot crop plant. Nature Biotechnology 18:1101-1104; PCT publication number WO 01/13731; and Bing J W et al. (2000) Efficacy of Cry1F Transgenic Maize, 14<sup>th</sup> Biennial International Plant Resistance to Insects Workshop, Fort Collins, Colo.).

The expression of foreign genes in plants is known to be 65 influenced by their location in the plant genome, perhaps due to chromatin structure (e.g., heterochromatin) or the prox-

imity of transcriptional regulatory elements (e.g., enhancers) close to the integration site (Weising et al., Ann. Rev. Genet 22:421-477, 1988). At the same time the presence of the transgene at different locations in the genome will influence the overall phenotype of the plant in different ways. For this reason, it is often necessary to screen a large number of events in order to identify an event characterized by optimal expression of an introduced gene of interest. For example, it has been observed in plants and in other organisms that there may be a wide variation in levels of expression of an introduced gene among events. There may also be differences in spatial or temporal patterns of expression, for example, differences in the relative expression of a transgene in various plant tissues, that may not correspond to the patterns expected from transcriptional regulatory elements present in the introduced gene construct. For this reason, it is common to produce hundreds to thousands of different events and screen those events for a single event that has desired transgene expression levels and patterns for com-20 mercial purposes. An event that has desired levels or patterns of transgene expression is useful for introgressing the transgene into other genetic backgrounds by sexual outcrossing using conventional breeding methods. Progeny of such crosses maintain the transgene expression characteristics of 25 the original transformant. This strategy is used to ensure reliable gene expression in a number of varieties that are well adapted to local growing conditions.

It would be advantageous to be able to detect the presence of a particular event in order to determine whether progeny of a sexual cross contain a transgene of interest. In addition, a method for detecting a particular event would be helpful for complying with regulations requiring the pre-market approval and labeling of foods derived from recombinant crop plants, for example, or for use in environmental moni-35 toring, monitoring traits in crops in the field, or monitoring products derived from a crop harvest, as well as for use in ensuring compliance of parties subject to regulatory or contractual terms.

It is possible to detect the presence of a transgene by any nucleic acid detection method known in the art including, but not limited to, the polymerase chain reaction (PCR) or DNA hybridization using nucleic acid probes. These detection methods generally focus on frequently used genetic elements, such as promoters, terminators, marker genes, etc., because for many DNA constructs, the coding region is interchangeable. As a result, such methods may not be useful for discriminating between different events, particularly those produced using the same DNA construct or very similar constructs unless the DNA sequence of the flanking DNA adjacent to the inserted heterologous DNA is known. For example, an event-specific PCR assay is described in U.S. Pat. No. 6,395,485 for the detection of elite event GAT-ZM1. Accordingly, it would be desirable to have a simple and discriminative method for the identification of event TC1507.

### SUMMARY OF THE INVENTION

This invention relates preferably to methods for produc-More specifically, a DNA construct is provided that when expressed in plant cells and plants confers resistance to insects. According to one aspect of the invention, a DNA construct, capable of introduction into and replication in a host cell, is provided that when expressed in plant cells and plants confers insect resistance to the plant cells and plants. The DNA construct is comprised of a DNA molecule named

PHI8999A and it includes two transgene expression cassettes. The first expression cassette comprises a DNA molecule which includes the promoter, 5' untranslated exon, and first intron of the maize ubiquitin (Ubi-1) gene (Christensen et al. (1992) Plant Mol. Biol. 18:675-689 and Christensen and Quail (1996) Transgenic Res. 5:213-218) operably connected to a DNA molecule encoding a B.t. δ-endotoxin identified as Cry1F (U.S. Pat. Nos. 5,188,960 and 6,218, 188) operably connected to a DNA molecule comprising a 3' ORF25 transcriptional terminator isolated from Agrobacterium tumefaciens (Barker et al. (1983) Plant Mol. Biol. 2:335-350). The second transgene expression cassette of the DNA construct comprises a DNA molecule of the cauliflower mosaic virus (CaMV) 35S promoter (Odell J. T. et al.  $_{15}$ (1985) Nature 313: 810-812; Mitsuhara et al. (1996) Plant Cell Physiol. 37: 49-59) operably connected to a DNA molecule encoding a phosphinothricin acetyltransferase (PAT) gene (Wohlleben W. et al. (1988) *Gene* 70: 25-37) operably connected to a DNA molecule comprising a 3' transcriptional terminator from (CaMV) 35S (see Mitsuhara et al. (1996) Plant Cell Physiol. 37: 49-59). Plants containing the DNA construct are also provided.

According to another aspect of the invention, compositions and methods are provided for identifying a novel corn plant designated TC1507, which methods are based on primers or probes which specifically recognize the 5' and/or 3' flanking sequence of TC1507. DNA molecules are provided that comprise primer sequences that when utilized in a PCR reaction will produce amplicons unique to the transgenic event TC1507. These molecules may be selected from the group consisting of:

```
(SEQ ID NO: 1)

5'-GTAGTACTATAGATTATATTATTCGTAGAG-3';

(SEQ ID NO: 2)

5'-GCCATACAGAACTCAAAATCTTTTCCGGAG-3';

(SEQ ID NO: 23)

5'-CTTCAAACAAGTGTGACAAA-3';

(SEQ ID NO: 3)

5'-TGTGGTGTTTTGTGGCTCTTAA-3';

(SEQ ID NO: 4)

5'-AGCACCTTTTCATTCTTTCATATAC-3';

(SEQ ID NO: 5)

5'-GACCTCCCCA CAGGCATGAT TGATC-3';
```

and complements thereof. The corn plant and seed comprising these molecules is an aspect of this invention. Further, 50 kits utilizing these primer sequences for the identification of the TC1507 event are provided.

An additional aspect of the invention relates to the specific flanking sequences of TC1507 described herein, which can be used to develop specific identification methods for 55 TC1507 in biological samples. More particularly, the invention relates to the 5' and/or 3' flanking regions of TC1507, SEQ ID NO:21 and SEQ ID NO:22, respectively, which can be used for the development of specific primers and probes. The invention further relates to identification methods for 60 the presence of TC1507 in biological samples based on the use of such specific primers or probes.

According to another aspect of the invention, methods of detecting the presence of DNA corresponding to the corn event TC1507 in a sample are provided. Such methods 65 comprise: (a) contacting the sample comprising DNA with a DNA primer set, that when used in a nucleic acid amplifi-

4

cation reaction with genomic DNA extracted from corn event TC1507 produces an amplicon that is diagnostic for corn event TC1507; (b) performing a nucleic acid amplification reaction, thereby producing the amplicon; and (c) detecting the amplicon.

DNA molecules that comprise the novel transgene/flanking insertion region, SEQ ID NO: 26 and SEQ ID NO: 27 and are homologous or complementary to SEQ ID NO: 26 and SEQ ID NO: 27 are an aspect of this invention.

DNA sequences that comprise the novel transgene/flanking insertion region, SEQ ID NO:26 are an aspect of this invention. DNA sequences that comprise a sufficient length of polynucleotides of transgene insert sequence and a sufficient length of polynucleotides of maize genomic and/or flanking sequence from maize plant TC1507 of SEQ ID NO:26 that are useful as primer sequences for the production of an amplicon product diagnostic for maize plant TC1507 are included.

In addition, DNA sequences that comprise the novel transgene/flanking insertion region, SEQ ID NO:27 are provided. DNA sequences that comprise a sufficient length of polynucleotides of transgene insert sequence and a sufficient length of polynucleotides of maize genomic and/or flanking sequence from maize plant TC1507 of SEQ ID NO:27 that are useful as primer sequences for the production of an amplicon product diagnostic for maize plant TC1507 are included.

According to another aspect of the invention, the DNA sequences that comprise at least 11 or more nucleotides of the transgene portion of the DNA sequence of SEQ ID NO:26 or complements thereof, and a similar length of 5' flanking maize DNA sequence of SEQ ID NO:26 or complements thereof are useful as DNA primers in DNA amplification methods. The amplicons produced using these primers are diagnostic for maize event TC1507. Therefore, the invention also includes the amplicons produced by DNA primers homologous or complementary to SEQ ID NO:26.

According to another aspect of the invention, the DNA sequences that comprise at least 11 or more nucleotides of the transgene portion of the DNA sequence of SEQ ID NO:27 or complements thereof, and a similar length of 3' flanking maize DNA sequence of SEQ ID NO:27 or complements thereof are useful as DNA primers in DNA amplification methods. The amplicons produced using these primers are diagnostic for maize event TC1507. Therefore, the invention also includes the amplicons produced by DNA primers homologous or complementary to SEQ ID NO:27.

More specifically, a pair of DNA molecules comprising a DNA primer set, wherein the DNA molecules are identified as SEQ ID NO: 1 or complements thereof and SEQ ID NO: 2 or complements thereof; SEQ ID NO: 2 or complements thereof and SEQ ID NO: 23 or complements thereof; SEQ ID NO: 3 or complements thereof and SEQ ID NO: 5 or complements thereof; SEQ ID NO: 4 or complements thereof and SEQ ID NO: 5 or complements thereof are aspects of the invention.

Further aspects of the invention include the amplicon comprising the DNA molecules of SEQ ID NO: 1 and SEQ ID NO: 2; the amplicon comprising the DNA molecules of SEQ ID NO: 2 and SEQ ID NO: 23; the amplicon comprising the DNA molecules of SEQ ID NO: 3 and SEQ ID NO: 5; and the amplicon comprising the DNA molecules of SEQ ID NO: 4 and SEQ ID NO: 5.

According to another aspect of the invention, methods of detecting the presence of a DNA molecule corresponding to the TC1507 event in a sample, such methods comprising: (a) contacting the sample comprising DNA extracted from a

corn plant with a DNA probe, molecule that hybridizes under stringent hybridization conditions with DNA extracted from corn event TC1507 and does not hybridize under the stringent hybridization conditions with a control corn plant DNA; (b) subjecting the sample and probe to 5 stringent hybridization conditions; and (c) detecting hybridization of the probe to the DNA. More specifically, a method for detecting the presence of a DNA molecule corresponding to the TC1507 event in a sample, such methods, consisting of (a) contacting the sample comprising DNA extracted from 10 a corn plant with a DNA probe molecule that consists of sequences that are unique to the event, e.g. junction sequences, wherein said DNA probe molecule hybridizes under stringent hybridization conditions with DNA extracted from corn event TC1507 and does not hybridize 15 under the stringent hybridization conditions with a control corn plant DNA; (b) subjecting the sample and probe to stringent hybridization conditions; and (c) detecting hybridization of the probe to the DNA.

In addition, a kit and methods for identifying event 20 TC1507 in a biological sample which detects a TC1507 specific region within SEQ ID NO: 24 are provided.

DNA molecules are provided that comprise at least one junction sequence of TC1507 selected from the group consisting of SEQ ID NO:45, 46, 47, 48, 49, 50, 51, 52, 53, 54, 25 55, 56 and 57 and complements thereof; wherein a junction sequence spans the junction between heterologous DNA inserted into the genome and the DNA from the corn cell flanking the insertion site, i.e. flanking DNA, and is diagnostic for the TC1507 event.

According to another aspect of the invention, methods of producing an insect resistant corn plant that comprise the steps of: (a) sexually crossing a first parental corn line comprising the expression cassettes of the present invention, which confers resistance to insects, and a second parental 35 ordinary skill in the art in the practice of the present corn line that lacks insect resistance, thereby producing a plurality of progeny plants; and (b) selecting a progeny plant that is insect resistant. Such methods may optionally comprise the further step of back-crossing the progeny plant to the second parental corn line to producing a true-breeding 40 corn plant that is insect resistant.

The present invention provides a method of producing a corn plant that is resistant to insects comprising transforming a corn cell with the DNA construct PHI8999A (SEQ ID NO:25), growing the transformed corn cell into a corn plant, 45 selecting the corn plant that shows resistance to insects, and further growing the corn plant into a fertile corn plant. The fertile corn plant can be self pollinated or crossed with compatible corn varieties to produce insect resistant progeny.

The invention further relates to a DNA detection kit for identifying maize event TC1507 in biological samples. Preferably the kit of the invention comprises a first primer which specifically recognizes the 5' or 3' flanking region of TC1507, and a second primer which specifically recognizes 55 a sequence within the foreign DNA of TC1507, or within the flanking DNA, for use in a PCR identification protocol. The invention also relates to a kit for identifying event TC1507 in biological samples, which kit comprises a specific probe having a sequence which corresponds or is complementary 60 to, a sequence having between 80% and 100% sequence identity with a specific region of event TC1507. Preferably the sequence of the probe corresponds to a specific region comprising part of the 5' or 3' flanking region of event TC1507.

The methods and kits encompassed by the present invention can be used for different purposes such as, but not

limited to the following: to identify event TC1507 in plants, plant material or in products such as, but not limited to, food or feed products (fresh or processed) comprising, or derived from plant material; additionally or alternatively, the methods and kits of the present invention can be used to identify transgenic plant material for purposes of segregation between transgenic and non-transgenic material; additionally or alternatively, the methods and kits of the present invention can be used to determine the quality of plant material comprising maize event TC1507. The kits may also contain the reagents and materials necessary for the performance of the detection method.

This invention further relates to the TC1507 corn plant or its parts, including, but not limited to, pollen, ovules, vegetative cells, the nuclei of pollen cells, and the nuclei of egg cells of the corn plant TC1507 and the progeny derived thereof. The corn plant and seed TC1507 from which the DNA primer molecules of the present invention provide a specific amplicon product is an aspect of the invention.

The foregoing and other aspects of the invention will become more apparent from the following detailed description and accompanying drawing.

#### BRIEF DESCRIPTION OF THE DRAWINGS

FIG. 1. Linear map showing the transgenic insert PHI8999A, as well as the sequences flanking the transgenic insert.

#### DETAILED DESCRIPTION

The following definitions and methods are provided to better define the present invention and to guide those of invention. Unless otherwise noted, terms are to be understood according to conventional usage by those of ordinary skill in the relevant art. Definitions of common terms in molecular biology may also be found in Rieger et al., Glossary of Genetics: Classical and Molecular, 5<sup>th</sup> edition, Springer-Verlag; New York, 1991; and Lewin, Genes V, Oxford University Press: New York, 1994. The nomenclature for DNA bases as set forth at 37 CFR 1.822 is used.

As used herein, the term "comprising" means "including but not limited to".

As used herein, the term "corn" means *Zea mays* or maize and includes all plant varieties that can be bred with corn, including wild maize species.

As used herein, the term "TC1507 specific" refers to a 50 nucleotide sequence which is suitable for discriminatively identifying event TC1507 in plants, plant material, or in products such as, but not limited to, food or feed products (fresh or processed) comprising, or derived from plant material.

As used herein, the terms "insect resistant" and "impacting insect pests" refers to effecting changes in insect feeding, growth, and/or behavior at any stage of development, including but not limited to: killing the insect; retarding growth; preventing reproductive capability; and the like.

As used herein, the terms "pesticidal activity" and "insecticidal activity" are used synonymously to refer to activity of an organism or a substance (such as, for example, a protein) that can be measured by numerous parameters including, but not limited to, pest mortality, pest weight loss, pest attrac-65 tion, pest repellency, and other behavioral and physical changes of a pest after feeding on and/or exposure to the organism or substance for an appropriate length of time. For

example "pesticidal proteins" are proteins that display pesticidal activity by themselves or in combination with other proteins.

"Coding sequence" refers to a nucleotide sequence that codes for a specific amino acid sequence. As used herein, the terms "encoding" or "encoded" when used in the context of a specified nucleic acid mean that the nucleic acid comprises the requisite information to guide translation of the nucleotide sequence into a specified protein. The information by which a protein is encoded is specified by the use of codons. 10 A nucleic acid encoding a protein may comprise nontranslated sequences (e.g., introns) within translated regions of the nucleic acid or may lack such intervening nontranslated sequences (e.g., as in cDNA).

specific protein, including regulatory sequences preceding (5' non-coding sequences) and following (3' non-coding sequences) the coding sequence. "Native gene" refers to a gene as found in nature with its own regulatory sequences. "Chimeric gene" refers any gene that is not a native gene, 20 comprising regulatory and coding sequences that are not found together in nature. Accordingly, a chimeric gene may comprise regulatory sequences and coding sequences that are derived from different sources, or regulatory sequences and coding sequences derived from the same source, but 25 arranged in a manner different than that found in nature. "Endogenous gene" refers to a native gene in its natural location in the genome of an organism. "Foreign" refers to material not normally found in the location of interest. Thus "foreign DNA" may comprise both recombinant DNA as 30 well as newly introduced, rearranged DNA of the plant. A "foreign" gene refers to a gene not normally found in the host organism, but that is introduced into the host organism by gene transfer. Foreign genes can comprise native genes inserted into a non-native organism, or chimeric genes. A 35 is understood by those skilled in the art that different "transgene" is a gene that has been introduced into the genome by a transformation procedure. The site in the plant genome where a recombinant DNA has been inserted may be referred to as the "insertion site" or "target site".

As used herein, "insert DNA" refers to the heterologous 40 DNA within the expression cassettes used to transform the plant material while "flanking DNA" can exist of either genomic DNA naturally present in an organism such as a plant, or foreign (heterologous) DNA introduced via the transformation process which is extraneous to the original 45 insert DNA molecule, e.g. fragments associated with the transformation event. A "flanking region" or "flanking sequence" as used herein refers to a sequence of at least 20 base pair, preferably at least 50 base pair, and up to 5000 base pair which is located either immediately upstream of 50 and contiguous with or immediately downstream of and contiguous with the original foreign insert DNA molecule. Transformation procedures leading to random integration of the foreign DNA will result in transformants containing different flanking regions characteristic and unique for each 55 transformant. When recombinant DNA is introduced into a plant through traditional crossing, its flanking regions will generally not be changed. Transformants will also contain unique junctions between a piece of heterologous insert DNA and genomic DNA, or 2 pieces of genomic DNA, or 60 2 pieces of heterologous DNA. A "junction" is a point where 2 specific DNA fragments join. For example, a junction exists where insert DNA joins flanking DNA. A junction point also exists in a transformed organism where 2 DNA fragments join together in a manner that is modified from 65 that found in the native organism. "Junction DNA" refers to DNA that comprises a junction point.

As used herein, "heterologous" in reference to a nucleic acid is a nucleic acid that originates from a foreign species, or, if from the same species, is substantially modified from its native form in composition and/or genomic locus by deliberate human intervention. For example, a promoter operably linked to a heterologous nucleotide sequence can be from a species different from that from which the nucleotide sequence was derived, or, if from the same species, the promoter is not naturally found operably linked to the nucleotide sequence. A heterologous protein may originate from a foreign species, or, if from the same species, is substantially modified from its original form by deliberate human intervention.

"Regulatory sequences" refer to nucleotide sequences "Gene" refers to a nucleic acid fragment that expresses a 15 located upstream (5' non-coding sequences), within, or downstream (3' non-coding sequences) of a coding sequence, and which influence the transcription, RNA processing or stability, or translation of the associated coding sequence. Regulatory sequences may include promoters, translation leader sequences, introns, and polyadenylation recognition sequences.

> "Promoter" refers to a nucleotide sequence capable of controlling the expression of a coding sequence or functional RNA. In general, a coding sequence is located 3' to a promoter sequence. The promoter sequence consists of proximal and more distal upstream elements, the latter elements are often referred to as enhancers. Accordingly, an "enhancer" is a nucleotide sequence that can stimulate promoter activity and may be an innate element of the promoter or a heterologous element inserted to enhance the level or tissue-specificity of a promoter. Promoters may be derived in their entirety from a native gene, or be composed of different elements derived from different promoters found in nature, or even comprise synthetic nucleotide segments. It promoters may direct the expression of a gene in different tissues or cell types, or at different stages of development, or in response to different environmental conditions. Promoters that cause a nucleic acid fragment to be expressed in most cell types at most times are commonly referred to as "constitutive promoters". New promoters of various types useful in plant cells are constantly being discovered; numerous examples may be found in the compilation by Okamuro and Goldberg (1989) Biochemistry of Plants 15:1-82. It is further recognized that since in most cases the exact boundaries of regulatory sequences have not been completely defined, nucleic acid fragments of different lengths may have identical promoter activity.

> The "translation leader sequence" refers to a nucleotide sequence located between the promoter sequence of a gene and the coding sequence. The translation leader sequence is present in the fully processed mRNA upstream of the translation start sequence. The translation leader sequence may affect numerous parameters including, processing of the primary transcript to mRNA, mRNA stability and/or translation efficiency. Examples of translation leader sequences have been described (Turner and Foster (1995) Mol. Biotechnol. 3:225-236).

> The "3' non-coding sequences" refer to nucleotide sequences located downstream of a coding sequence and include polyadenylation recognition sequences and other sequences encoding regulatory signals capable of affecting mRNA processing or gene expression. The polyadenylation signal is usually characterized by affecting the addition of polyadenylic acid tracts to the 3' end of the mRNA precursor. The use of different 3' non-coding sequences is exemplified by Ingelbrecht et al. (1989) Plant Cell 1:671-680.

A "protein" or "polypeptide" is a chain of amino acids arranged in a specific order determined by the coding sequence in a polynucleotide encoding the polypeptide.

A DNA construct is an assembly of DNA molecules linked together that provide one or more expression cas- 5 settes. The DNA construct may be a plasmid that is enabled for self replication in a bacterial cell and contains various endonuclease enzyme restriction sites that are useful for introducing DNA molecules that provide functional genetic elements, i.e., promoters, introns, leaders, coding sequences, 10 3' termination regions, among others; or a DNA construct may be a linear assembly of DNA molecules, such as an expression cassette. The expression cassette contained within a DNA construct comprise the necessary genetic elements to provide transcription of a messenger RNA. The 15 expression cassette can be designed to express in prokaryote cells or eukaryotic cells. Expression cassettes of the present invention are designed to express most preferably in plant cells.

The DNA molecules of the invention are provided in 20 expression cassettes for expression in an organism of interest. The cassette will include 5' and 3' regulatory sequences operably linked to a coding sequence of the invention. "Operably linked" means that the nucleic acid sequences being linked are contiguous and, where necessary to join two 25 protein coding regions, contiguous and in the same reading frame. Operably linked is intended to indicate a functional linkage between a promoter and a second sequence, wherein the promoter sequence initiates and mediates transcription of the DNA sequence corresponding to the second sequence. 30 The cassette may additionally contain at least one additional gene to be cotransformed into the organism. Alternatively, the additional gene(s) can be provided on multiple expression cassettes or multiple DNA constructs.

tion of transcription: a transcriptional and translational initiation region, a coding region, and a transcriptional and translational termination region functional in the organism serving as a host. The transcriptional initiation region (i.e., the promoter) may be native or analogous, or foreign or 40 heterologous to the host organism. Additionally, the promoter may be the natural sequence or alternatively a synthetic sequence. The expression cassettes may additionally contain 5' leader sequences in the expression cassette construct. Such leader sequences can act to enhance translation. 45

It is to be understood that as used herein the term "transgenic" includes any cell, cell line, callus, tissue, plant part, or plant the genotype of which has been altered by the presence of a heterologous nucleic acid including those transgenics initially so altered as well as those created by 50 sexual crosses or asexual propagation from the initial transgenic. The term "transgenic" as used herein does not encompass the alteration of the genome (chromosomal or extrachromosomal) by conventional plant breeding methods or by naturally occurring events such as random cross-fertilization, non-recombinant viral infection, non-recombinant bacterial transformation, non-recombinant transposition, or spontaneous mutation.

A transgenic "event" is produced by transformation of plant cells with a heterologous DNA construct(s), including 60 a nucleic acid expression cassette that comprises a transgene of interest, the regeneration of a population of plants resulting from the insertion of the transgene into the genome of the plant, and selection of a particular plant characterized by insertion into a particular genome location. An event is 65 characterized phenotypically by the expression of the transgene. At the genetic level, an event is part of the genetic

makeup of a plant. The term "event" also refers to progeny produced by a sexual outcross between the transformant and another variety that include the heterologous DNA. Even after repeated back-crossing to a recurrent parent, the inserted DNA and flanking DNA from the transformed parent is present in the progeny of the cross at the same chromosomal location. The term "event" also refers to DNA from the original transformant comprising the inserted DNA and flanking sequence immediately adjacent to the inserted DNA that would be expected to be transferred to a progeny that receives inserted DNA including the transgene of interest as the result of a sexual cross of one parental line that includes the inserted DNA (e.g., the original transformant and progeny resulting from selfing) and a parental line that does not contain the inserted DNA.

An insect resistant TC1507 corn plant can be bred by first sexually crossing a first parental corn plant consisting of a corn plant grown from the transgenic TC1507 corn plant and progeny thereof derived from transformation with the expression cassettes of the present invention that confers insect resistance, and a second parental corn plant that lacks insect resistance, thereby producing a plurality of first progeny plants; and then selecting a first progeny plant that is resistant to insects; and selfing the first progeny plant, thereby producing a plurality of second progeny plants; and then selecting from the second progeny plants an insect resistant plant. These steps can further include the backcrossing of the first insect resistant progeny plant or the second insect resistant progeny plant to the second parental corn plant or a third parental corn plant, thereby producing a corn plant that is resistant to insects.

As used herein, the term "plant" includes reference to whole plants, plant organs (e.g., leaves, stems, roots, etc.), seeds, plant cells, and progeny of same. Parts of transgenic The expression cassette will include in the 5' to 3' direc- 35 plants understood to be within the scope of the invention comprise, for example, plant cells, protoplasts, tissues, callus, embryos as well as flowers, stems, fruits, leaves, and roots originating in transgenic plants or their progeny previously transformed with a DNA molecule of the invention and therefore consisting at least in part of transgenic cells, are also an aspect of the present invention.

As used herein, the term "plant cell" includes, without limitation, seeds, suspension cultures, embryos, meristematic regions, callus tissue, leaves, roots, shoots, gametophytes, sporophytes, pollen, and microspores. The class of plants that can be used in the methods of the invention is generally as broad as the class of higher plants amenable to transformation techniques, including both monocotyledonous and dicotyledonous plants.

"Transformation" refers to the transfer of a nucleic acid fragment into the genome of a host organism, resulting in genetically stable inheritance. Host organisms containing the transformed nucleic acid fragments are referred to as "transgenic" organisms. Examples of methods of plant transformation include *Agrobacterium*-mediated transformation (De Blaere et al. (1987) Meth. Enzymol. 143:277) and particle-accelerated or "gene gun" transformation technology (Klein et al. (1987) Nature (London) 327:70-73; U.S. Pat. No. 4,945,050, incorporated herein by reference). Additional transformation methods are disclosed below.

Thus, isolated polynucleotides of the present invention can be incorporated into recombinant constructs, typically DNA constructs, which are capable of introduction into and replication in a host cell. Such a construct can be a vector that includes a replication system and sequences that are capable of transcription and translation of a polypeptideencoding sequence in a given host cell. A number of vectors

suitable for stable transfection of plant cells or for the establishment of transgenic plants have been described in, e.g., Pouwels et al., (1985; Supp. 1987) Cloning Vectors: A Laboratory Manual, Weissbach and Weissbach (1989) Methods for Plant Molecular Biology, (Academic Press, 5 New York); and Flevin et al., (1990) Plant Molecular Biology Manual, (Kluwer Academic Publishers). Typically, plant expression vectors include, for example, one or more cloned plant genes under the transcriptional control of 5' and 3' regulatory sequences and a dominant selectable marker. 10 Such plant expression vectors also can contain a promoter regulatory region (e.g., a regulatory region controlling inducible or constitutive, environmentally- or developmentally-regulated, or cell- or tissue-specific expression), a transcription initiation start site, a ribosome binding site, an 15 RNA processing signal, a transcription termination site, and/or a polyadenylation signal.

It is also to be understood that two different transgenic plants can also be mated to produce offspring that contain two independently segregating added, exogenous genes. 20 Selfing of appropriate progeny can produce plants that are homozygous for both added, exogenous genes. Back-crossing to a parental plant and out-crossing with a non-transgenic plant are also contemplated, as is vegetative propagation. Descriptions of other breeding methods that are 25 commonly used for different traits and crops can be found in one of several references, e.g., Fehr, in *Breeding Methods for Cultivar Development*, Wilcos J. ed., American Society of Agronomy, Madison Wis. (1987).

A "probe" is an isolated nucleic acid to which is attached a conventional detectable label or reporter molecule, e.g., a radioactive isotope, ligand, chemiluminescent agent, or enzyme. Such a probe is complementary to a strand of a target nucleic acid, in the case of the present invention, to a strand of isolated DNA from corn event TC1507 whether 35 from a corn plant or from a sample that includes DNA from the event. Probes according to the present invention include not only deoxyribonucleic or ribonucleic acids but also polyamides and other probe materials that bind specifically to a target DNA sequence and can be used to detect the 40 presence of that target DNA sequence.

"Primers" are isolated nucleic acids that are annealed to a complementary target DNA strand by nucleic acid hybridization to form a hybrid between the primer and the target DNA strand, then extended along the target DNA strand by 45 a polymerase, e.g., a DNA polymerase. Primer pairs of the present invention refer to their use for amplification of a target nucleic acid sequence, e.g., by the polymerase chain reaction (PCR) or other conventional nucleic-acid amplification methods. "PCR" or "polymerase chain reaction" is a 50 technique used for the amplification of specific DNA segments (see, U.S. Pat. Nos. 4,683,195 and 4,800,159; herein incorporated by reference).

Probes and primers are of sufficient nucleotide length to bind to the target DNA sequence specifically in the hybridization conditions or reaction conditions determined by the operator. This length may be of any length that is of sufficient length to be useful in a detection method of choice. Generally, 11 nucleotides or more in length, preferably 18 nucleotides or more, and more preferably 22 nucleotides or more, are used. Such probes and primers hybridize specifically to a target sequence under high stringency hybridization conditions. Preferably, probes and primers according to the present invention have complete DNA sequence similarity of contiguous nucleotides with the target sequence, 65 although probes differing from the target DNA sequence and that retain the ability to hybridize to target DNA sequences

12

may be designed by conventional methods. Probes can be used as primers, but are generally designed to bind to the target DNA or RNA and not be used in an amplification process.

Specific primers can be used to amplify an integration fragment to produce an amplicon that can be used as a "specific probe" for identifying event TC1507 in biological samples. When the probe is hybridized with the nucleic acids of a biological sample under conditions which allow for the binding of the probe to the sample, this binding can be detected and thus allow for an indication of the presence of event TC1507 in the biological sample. Such identification of a bound probe has been described in the art. The specific probe is preferably a sequence which, under optimized conditions, hybridizes specifically to a region within the 5' or 3' flanking region of the event and preferably also comprises a part of the foreign DNA contiguous therewith. Preferably the specific probe comprises a sequence of at least 80%, preferably between 80 and 85%, more preferably between 85 and 90%, especially preferably between 90 and 95%, and most preferably between 95 and 100% identical (or complementary) to a specific region of the event.

Methods for preparing and using probes and primers are described, for example, in *Molecular Cloning: A Laboratory* Manual, 2<sup>nd</sup> ed., vol. 1-3, ed. Sambrook et al., Cold Spring Harbor Laboratory Press, Cold Spring Harbor, N.Y. 1989 (hereinafter, "Sambrook et al., 1989"); Current Protocols in Molecular Biology, ed. Ausubel et al., Greene Publishing and Wiley-Interscience, New York, 1992 (with periodic updates) (hereinafter, "Ausubel et al., 1992"); and Innis et al., PCR Protocols: A Guide to Methods and Applications, Academic Press: San Diego, 1990. PCR primer pairs can be derived from a known sequence, for example, by using computer programs intended for that purpose such as the PCR primer analysis tool in Vector NTI version 6 (Informax Inc., Bethesda Md.); PrimerSelect (DNASTAR Inc., Madison, Wis.); and Primer (Version 0.5, © 1991, Whitehead Institute for Biomedical Research, Cambridge, Mass.). Additionally, the sequence can be visually scanned and primers manually identified using guidelines known to one of skill in the art.

A "kit" as used herein refers to a set of reagents for the purpose of performing the method of the invention, more particularly, the identification of the event TC1507 in biological samples. The kit of the invention can be used, and its components can be specifically adjusted, for purposes of quality control (e.g. purity of seed lots), detection of event TC1507 in plant material, or material comprising or derived from plant material, such as but not limited to food or feed products. "Plant material" as used herein refers to material which is obtained or derived from a plant.

Primers and probes based on the flanking DNA and insert sequences disclosed herein can be used to confirm (and, if necessary, to correct) the disclosed sequences by conventional methods, e.g., by re-cloning and sequencing such sequences. The nucleic acid probes and primers of the present invention hybridize under stringent conditions to a target DNA sequence. Any conventional nucleic acid hybridization or amplification method can be used to identify the presence of DNA from a transgenic event in a sample. Nucleic acid molecules or fragments thereof are capable of specifically hybridizing to other nucleic acid molecules under certain circumstances. As used herein, two nucleic acid molecules are said to be capable of specifically hybridizing to one another if the two molecules are capable of forming an anti-parallel, double-stranded nucleic acid structure.

A nucleic acid molecule is said to be the "complement" of another nucleic acid molecule if they exhibit complete complementarity. As used herein, molecules are said to exhibit "complete complementarity" when every nucleotide of one of the molecules is complementary to a nucleotide of 5 the other. Two molecules are said to be "minimally complementary" if they can hybridize to one another with sufficient stability to permit them to remain annealed to one another under at least conventional "low-stringency" conditions. Similarly, the molecules are said to be "complementary" if 10 they can hybridize to one another with sufficient stability to permit them to remain annealed to one another under conventional "high-stringency" conditions. Conventional stringency conditions are described by Sambrook et al., 1989, and by Haymes et al., In: Nucleic Acid Hybridization, 15 a Practical Approach, IRL Press, Washington, D.C. (1985), departures from complete complementarity are therefore permissible, as long as such departures do not completely preclude the capacity of the molecules to form a doublestranded structure. In order for a nucleic acid molecule to 20 serve as a primer or probe it need only be sufficiently complementary in sequence to be able to form a stable double-stranded structure under the particular solvent and salt concentrations employed.

In hybridization reactions, specificity is typically the 25 function of post-hybridization washes, the critical factors being the ionic strength and temperature of the final wash solution. The thermal melting point (Tm) is the temperature (under defined ionic strength and pH) at which 50% of a complementary target sequence hybridizes to a perfectly 30 matched probe. For DNA-DNA hybrids, the Tm can be approximated from the equation of Meinkoth and Wahl (1984) Anal. Biochem. 138:267-284: Tm=81.5° C.+16.6  $(\log M)+0.41 (\% GC)-0.61 (\% form)-500/L$ ; where M is the guanosine and cytosine nucleotides in the DNA, % form is the percentage of formamide in the hybridization solution, and L is the length of the hybrid in base pairs. Tm is reduced by about 1° C. for each 1% of mismatching; thus, Tm, hybridization, and/or wash conditions can be adjusted to 40 hybridize to sequences of the desired identity. For example, if sequences with >90% identity are sought, the Tm can be decreased 10° C. Generally, stringent conditions are selected to be about 5° C. lower than the Tm for the specific sequence and its complement at a defined ionic strength and pH. 45 However, severely stringent conditions can utilize a hybridization and/or wash at 1, 2, 3, or 4° C. lower than the Tm; moderately stringent conditions can utilize a hybridization and/or wash at 6, 7, 8, 9, or 10° C. lower than the Tm; low stringency conditions can utilize a hybridization and/or wash 50 at 11, 12, 13, 14, 15, or 20° C. lower than the Tm.

Using the equation, hybridization and wash compositions, and desired Tm, those of ordinary skill will understand that variations in the stringency of hybridization and/or wash solutions are inherently described. If the desired degree of 55 mismatching results in a Tm of less than 45° C. (aqueous solution) or 32° C. (formamide solution), it is preferred to increase the SSC concentration so that a higher temperature can be used. An extensive guide to the hybridization of nucleic acids is found in Tijssen (1993) Laboratory Tech- 60 niques in Biochemistry and Molecular Biology—Hybridization with Nucleic Acid Probes, Part I, Chapter 2 (Elsevier, New York); and Ausubel et al., eds. (1995) Current Protocols in Molecular Biology, Chapter 2 (Greene Publishing (1989) Molecular Cloning: A Laboratory Manual (2d ed., Cold Spring Harbor Laboratory Press, Plainview, N.Y.).

14

As used herein, a substantially homologous sequence is a nucleic acid molecule that will specifically hybridize to the complement of the nucleic acid molecule to which it is being compared under high stringency conditions. Appropriate stringency conditions which promote DNA hybridization, for example, 6x sodium chloride/sodium citrate (SSC) at about 45° C., followed by a wash of 2×SSC at 50° C., are known to those skilled in the art or can be found in *Current* Protocols in Molecular Biology, John Wiley & Sons, N.Y. (1989), 6.3.1-6.3.6. Typically, stringent conditions will be those in which the salt concentration is less than about 1.5 M Na ion, typically about 0.01 to 1.0 M Na ion concentration (or other salts) at pH 7.0 to 8.3 and the temperature is at least about 30° C. for short probes (e.g., 10 to 50 nucleotides) and at least about 60° C. for long probes (e.g., greater than 50 nucleotides). Stringent conditions may also be achieved with the addition of a destabilizing agent such as formamide. Exemplary low stringency conditions include hybridization with a buffer solution of 30 to 35% formamide, 1 M NaCl, 1% SDS (sodium dodecyl sulphate) at 37° C., and a wash in  $1 \times$  to  $2 \times SSC$  ( $20 \times SSC = 3.0$  M NaCl/0.3 M trisodium citrate) at 50 to 55° C. Exemplary moderate stringency conditions include hybridization in 40 to 45% formamide, 1 M NaCl, 1% SDS at 37° C., and a wash in 0.5× to 1×SSC at 55 to 60° C. Exemplary high stringency conditions include hybridization in 50% formamide, 1 M NaCl, 1% SDS at 37° C., and a wash in 0.1×SSC at 60 to 65° C. In a preferred embodiment, a nucleic acid of the present invention will specifically hybridize to one or more of the nucleic acid molecules unique to the TC1507 event or complements thereof or fragments of either under moderately stringent conditions.

Methods of alignment of sequences for comparison are well known in the art. Thus, the determination of percent molarity of monovalent cations, % GC is the percentage of 35 identity between any two sequences can be accomplished using a mathematical algorithm. Non-limiting examples of such mathematical algorithms are the algorithm of Myers and Miller (1988) CABIOS 4:11-17; the local homology algorithm of Smith et al. (1981) Adv. Appl. Math. 2:482; the homology alignment algorithm of Needleman and Wunsch (1970) J. Mol. Biol. 48:443-453; the search-for-similaritymethod of Pearson and Lipman (1988) *Proc. Natl. Acad. Sci.* 85:2444-2448; the algorithm of Karlin and Altschul (1990) Proc. Natl. Acad. Sci. USA 87:2264, modified as in Karlin and Altschul (1993) Proc. Natl. Acad. Sci. USA 90:5873-5877.

> Computer implementations of these mathematical algorithms can be utilized for comparison of sequences to determine sequence identity. Such implementations include, but are not limited to: CLUSTAL in the PC/Gene program (available from Intelligenetics, Mountain View, Calif.); the ALIGN program (Version 2.0); the ALIGN PLUS program (version 3.0, copyright 1997); and GAP, BESTFIT, BLAST, FASTA, and TFASTA in the Wisconsin Genetics Software Package, Version 10 (available from Accelrys, 9685 Scranton Road, San Diego, Calif. 92121, USA). Alignments using these programs can be performed using the default parameters.

The CLUSTAL program is well described by Higgins and Sharp, *Gene* 73: 237-244 (1988); Higgins and Sharp, *CABIOS* 5: 151-153 (1989); Corpet, et al., *Nucleic Acids* Research 16: 10881-90 (1988); Huang, et al., Computer Applications in the Biosciences 8: 155-65 (1992), and Pearson, et al., Methods in Molecular Biology 24: 307-331 and Wiley-Interscience, New York). See Sambrook et al. 65 (1994). The ALIGN and the ALIGN PLUS programs are based on the algorithm of Myers and Miller (1988) supra. The BLAST programs of Altschul et al. (1990) J. Mol. Biol.

215:403 are based on the algorithm of Karlin and Altschul (1990) supra. The BLAST family of programs which can be used for database similarity searches includes: BLASTN for nucleotide query sequences against nucleotide database sequences; BLASTX for nucleotide query sequences against protein database sequences; BLASTP for protein query sequences against protein database sequences; TBLASTN for protein query sequences against nucleotide database sequences; and TBLASTX for nucleotide query sequences against nucleotide database sequences. See, *Current Protocols in Molecular Biology*, Chapter 19, Ausubel, et al., Eds., Greene Publishing and Wiley-Interscience, New York (1995). Alignment may also be performed manually by inspection.

To obtain gapped alignments for comparison purposes, 15 Gapped BLAST (in BLAST 2.0) can be utilized as described in Altschul et al. (1997) *Nucleic Acids Res.* 25:3389. Alternatively, PSI-BLAST (in BLAST 2.0) can be used to perform an iterated search that detects distant relationships between molecules. See Altschul et al. (1997) supra. When 20 utilizing BLAST, Gapped BLAST, PSI-BLAST, the default parameters of the respective programs (e.g., BLASTN for nucleotide sequences, BLASTX for proteins) can be used. See www.ncbi.hlm.nih.gov.

As used herein, "sequence identity" or "identity" in the 25 context of two nucleic acid or polypeptide sequences makes reference to the residues in the two sequences that are the same when aligned for maximum correspondence over a specified comparison window. When percentage of sequence identity is used in reference to proteins it is 30 recognized that residue positions which are not identical often differ by conservative amino acid substitutions, where amino acid residues are substituted for other amino acid residues with similar chemical properties (e.g., charge or hydrophobicity) and therefore do not change the functional 35 properties of the molecule. When sequences differ in conservative substitutions, the percent sequence identity may be adjusted upwards to correct for the conservative nature of the substitution. Sequences that differ by such conservative substitutions are said to have "sequence similarity" or "simi- 40 larity". Means for making this adjustment are well known to those of skill in the art. Typically this involves scoring a conservative substitution as a partial rather than a full mismatch, thereby increasing the percentage sequence identity. Thus, for example, where an identical amino acid is 45 given a score of 1 and a non-conservative substitution is given a score of zero, a conservative substitution is given a score between zero and 1. The scoring of conservative substitutions is calculated, e.g., as implemented in the program PC/GENE (Intelligenetics, Mountain View, Calif.).

As used herein, "percentage of sequence identity" means the value determined by comparing two optimally aligned sequences over a comparison window, wherein the portion of the polynucleotide sequence in the comparison window may comprise additions or deletions (i.e., gaps) as compared 55 to the reference sequence (which does not comprise additions or deletions) for optimal alignment of the two sequences. The percentage is calculated by determining the number of positions at which the identical nucleic acid base or amino acid residue occurs in both sequences to yield the 60 number of matched positions, dividing the number of matched positions by the total number of positions in the window of comparison, and multiplying the result by 100 to yield the percentage of sequence identity.

Regarding the amplification of a target nucleic acid 65 sequence (e.g., by PCR) using a particular amplification primer pair, "stringent conditions" are conditions that permit

**16** 

the primer pair to hybridize only to the target nucleic-acid sequence to which a primer having the corresponding wild-type sequence (or its complement) would bind and preferably to produce a unique amplification product, the amplicon, in a DNA thermal amplification reaction.

The term "specific for (a target sequence)" indicates that a probe or primer hybridizes under stringent hybridization conditions only to the target sequence in a sample comprising the target sequence.

As used herein, "amplified DNA" or "amplicon" refers to the product of nucleic acid amplification of a target nucleic acid sequence that is part of a nucleic acid template. For example, to determine whether a corn plant resulting from a sexual cross contains transgenic event genomic DNA from the corn plant of the present invention, DNA extracted from the corn plant tissue sample may be subjected to a nucleic acid amplification method using a DNA primer pair that includes a first primer derived from flanking sequence adjacent to the insertion site of inserted heterologous DNA, and a second primer derived from the inserted heterologous DNA to produce an amplicon that is diagnostic for the presence of the event DNA. Alternatively, the second primer may be derived from the flanking sequence. The amplicon is of a length and has a sequence that is also diagnostic for the event. The amplicon may range in length from the combined length of the primer pairs plus one nucleotide base pair to any length of amplicon producible by a DNA amplification protocol. Alternatively, primer pairs can be derived from flanking sequence on both sides of the inserted DNA so as to produce an amplicon that includes the entire insert nucleotide sequence of the PHI8999A expression construct, see FIG. 1, approximately 6.2 Kb in size. A member of a primer pair derived from the flanking sequence may be located a distance from the inserted DNA sequence, this distance can range from one nucleotide base pair up to the limits of the amplification reaction, or about twenty thousand nucleotide base pairs. The use of the term "amplicon" specifically excludes primer dimers that may be formed in the DNA thermal amplification reaction.

Nucleic acid amplification can be accomplished by any of the various nucleic acid amplification methods known in the art, including the polymerase chain reaction (PCR). A variety of amplification methods are known in the art and are described, inter alia, in U.S. Pat. Nos. 4,683,195 and 4,683, 202 and in PCR Protocols: A Guide to Methods and Applications, ed. Innis et al., Academic press, San Diego, 1990. PCR amplification methods have been developed to amplify up to 22 Kb of genomic DNA and up to 42 Kb of bacteriophage DNA (Cheng et al., Proc. Natl. Acad. Sci. USA) 91:5695-5699, 1994). These methods as well as other methods known in the art of DNA amplification may be used in the practice of the present invention. It is understood that a number of parameters in a specific PCR protocol may need to be adjusted to specific laboratory conditions and may be slightly modified and yet allow for the collection of similar results. These adjustments will be apparent to a person skilled in the art.

The amplicon produced by these methods may be detected by a plurality of techniques, including, but not limited to, Genetic Bit Analysis (Nikiforov, et al. *Nucleic Acid Res.* 22:4167-4175, 1994) where a DNA oligonucleotide is designed which overlaps both the adjacent flanking DNA sequence and the inserted DNA sequence. The oligonucleotide is immobilized in wells of a microwell plate. Following PCR of the region of interest (using one primer in the inserted sequence and one in the adjacent flanking sequence) a single-stranded PCR product can be hybridized

to the immobilized oligonucleotide and serve as a template for a single base extension reaction using a DNA polymerase and labeled ddNTPs specific for the expected next base. Readout may be fluorescent or ELISA-based. A signal indicates presence of the insert/flanking sequence due to successful amplification, hybridization, and single base extension.

Another detection method is the Pyrosequencing technique as described by Winge (*Innov. Pharma. Tech.* 00: 18-24, 2000). In this method an oligonucleotide is designed 10 that overlaps the adjacent DNA and insert DNA junction. The oligonucleotide is hybridized to a single-stranded PCR product from the region of interest (one primer in the inserted sequence and one in the flanking sequence) and incubated in the presence of a DNA polymerase, ATP, 15 sulfurylase, luciferase, apyrase, adenosine 5' phosphosulfate and luciferin. dNTPs are added individually and the incorporation results in a light signal which is measured. A light signal indicates the presence of the transgene insert/flanking sequence due to successful amplification, hybridization, and 20 single or multi-base extension.

Fluorescence Polarization as described by Chen et al., (*Genome Res.* 9:492-498, 1999) is a method that can be used to detect an amplicon of the present invention. Using this method an oligonucleotide is designed which overlaps the 25 flanking and inserted DNA junction. The oligonucleotide is hybridized to a single-stranded PCR product from the region of interest (one primer in the inserted DNA and one in the flanking DNA sequence) and incubated in the presence of a DNA polymerase and a fluorescent-labeled ddNTP. Single 30 base extension results in incorporation of the ddNTP. Incorporation can be measured as a change in polarization using a fluorometer. A change in polarization indicates the presence of the transgene insert/flanking sequence due to successful amplification, hybridization, and single base extension.

Taqman® (PE Applied Biosystems, Foster City, Calif.) is described as a method of detecting and quantifying the presence of a DNA sequence and is fully understood in the instructions provided by the manufacturer. Briefly, a FRET 40 oligonucleotide probe is designed which overlaps the flanking and insert DNA junction. The FRET probe and PCR primers (one primer in the insert DNA sequence and one in the flanking genomic sequence) are cycled in the presence of a thermostable polymerase and dNTPs. Hybridization of the 45 FRET probe results in cleavage and release of the fluorescent moiety away from the quenching moiety on the FRET probe. A fluorescent signal indicates the presence of the flanking/transgene insert sequence due to successful amplification and hybridization.

Molecular Beacons have been described for use in sequence detection as described in Tyangi et al. (Nature Biotech. 14:303-308, 1996). Briefly, a FRET oligonucleotide probe is designed that overlaps the flanking and insert DNA junction. The unique structure of the FRET probe results in 55 it containing secondary structure that keeps the fluorescent and quenching moieties in close proximity. The FRET probe and PCR primers (one primer in the insert DNA sequence and one in the flanking sequence) are cycled in the presence of a thermostable polymerase and dNTPs. Following suc- 60 cessful PCR amplification, hybridization of the FRET probe to the target sequence results in the removal of the probe secondary structure and spatial separation of the fluorescent and quenching moieties. A fluorescent signal results. A fluorescent signal indicates the presence of the flanking/ 65 transgene insert sequence due to successful amplification and hybridization.

18

A hybridization reaction using a probe specific to a sequence found within the amplicon is yet another method used to detect the amplicon produced by a PCR reaction.

The present invention is further defined in the following Examples. It should be understood that these Examples, while indicating preferred embodiments of the invention, are given by way of illustration only. From the above discussion and these Examples, one skilled in the art can ascertain the essential characteristics of this invention, and without departing from the spirit and scope thereof, can make various changes and modifications of the invention to adapt it to various usages and conditions. Thus, various modifications of the invention, in addition to those shown and described herein, will be apparent to those skilled in the art from the foregoing description. Such modifications are also intended to fall within the scope of the appended claims.

The disclosure of each reference set forth herein is incorporated herein by reference in its entirety.

#### **EXAMPLES**

Example 1. Transformation of Maize by Particle Bombardment and Regeneration of Transgenic Plants Containing the Cry1F Gene

A DNA molecule of 6.2 Kb, designated PHI8999A (see FIG. 1 and SEQ ID NO:25), which includes a first transgene expression cassette comprising the promoter, 5' untranslated exon, and first intron of the maize ubiquitin (Ubi-1) gene (Christensen et al. (1992) Plant Mol. Biol. 18:675-689 and Christensen and Quail (1996) Transgenic Res. 5:213-218) operably connected to a DNA molecule encoding a *Bacillus* thuringiensis δ-endotoxin identified as Cry1F (U.S. Pat. Nos. 5,188,960 and 6,218,188) operably connected to a DNA molecule comprising a 3' ORF25 transcriptional terminator isolated from *Agrobacterium tumefaciens* (Barker et al. (1983) *Plant Mol. Biol.* 2:335-350), and a second transgene expression cassette comprising a DNA molecule of the cauliflower mosaic virus (CaMV) 35S promoter (Odell J. T. et al. (1985) Nature 313: 810-812; Mitsuhara et al. (1996) Plant Cell Physiol. 37:49-59) operably connected to a DNA molecule encoding the selectable marker, phosphinothricin acetyltransferase (PAT) gene (Wohlleben W. et al. (1988) Gene 70:25-37) operably connected to a DNA molecule comprising a 3' transcriptional terminator from (CaMV) 35S (see Mitsuhara et al. (1996) Plant Cell Physiol. 37:49-59) was used to transform maize embryo tissue.

B. t. Cry1F maize plants were obtained by microprojectile bombardment using the Biolistics® PDS-1000He particle 50 gun manufactured by Bio-Rad, Hercules, Calif.; essentially as described by Klein et al. (1987) *Nature*, UK 327(6117): 70-73. Immature embryos isolated from maize ears, harvested soon after pollination were cultured on callus initiation medium for several days. On the day of transformation, microscopic tungsten particles were coated with purified PHI8999A DNA (SEQ ID NO:25) and accelerated into the cultured embryos, where the insert DNA was incorporated into the cell chromosome. Only insert PHI8999A was used during transformation and no additional plasmid DNA was incorporated into the transformant. After bombardment, embryos were transferred to callus initiation medium containing glufosinate as the selection agent. Individual embryos were kept physically separate during culture, and the majority of explants died on the selective medium.

Those embryos that survived and produced healthy, glufosinate-resistant callus tissue were assigned unique identification codes representing putative transformation events,

and continually transferred to fresh selection medium. Plants were regenerated from tissue derived from each unique event and transferred to the greenhouse. Leaf samples were taken for molecular analysis to verify the presence of the transgene by PCR and to confirm expression of the Cry1F 5 protein by ELISA. Plants were then subjected to a whole plant bioassay using European corn borer insects. Positive plants were crossed with inbred lines to obtain seed from the initial transformed plants. A number of lines were evaluated in the field. The TC1507 event was selected from a population of independent transgenic events based on a superior combination of characteristics, including insect resistance and agronomic performance (see Bing J W et al. (2000) Efficacy of Cry1F Transgenic Maize, 14th Biennial International Plant Resistance to Insects Workshop, fort Collins, Colo., herein incorporated by reference).

Example 2. Identification of Nucleotides
Comprising the Flanking Sequence 5' to the
Transgenic Insert DNA in *Bacillus thuringiensis*Cry1F Maize Line TC1507

To identify a DNA fragment that included sequence 5' to the PHI8999A insert in event TC1507, Spe I restriction enzyme fragments from event TC1507 genomic DNA were <sup>25</sup> size selected on agarose gels, purified, and screened by Southern analysis to confirm hybridization to a Cry1F probe. Following confirmation of hybridization and fragment size, the fragments of interest were cloned into a pBluescript II SK (+)<sup>TM</sup> cloning vector to prepare an enriched size selected <sup>30</sup> plasmid based genomic DNA library. A probe homologous to a portion of the Cry1F gene was used to screen the plasmid library for positive clones. A positive clone was identified, purified by additional screening, and confirmed to result in a positive signal when hybridized to the Cry1F 35 probe. Nearly 3 Kb of the Spe I fragment contained in the isolated positive clone was sequenced using a primer walking approach. To initiate the first sequencing run, a primer that binds to a known sequence in the cloning vector DNA was designed to sequence a portion of the DNA of interest.

**20** 

A second sequencing run over the same region using another primer oriented in the reverse direction provided second strand coverage. Primer walking was accomplished by repeatedly using sequence data from previous runs to design new primers that were then used to extend the next round of sequencing further into the DNA of interest until the flanking sequence 5' to the inserted transgenic DNA in maize event TC1507 was obtained. Specific sequence information is provided in Example 4.

Example 3. Confirmation of Flanking Sequence 5' to the B.t. Cry1F Maize Line TC1507 Insert

To confirm the 5' flanking sequence of the B.t. Cry1F maize line TC1507 insert, PCR primer pairs were designed to obtain overlapping PCR products extending from the 5' flanking region into the full-length PHI8999A transgenic insert. PCR products were successfully amplified from B.t. Cry1F maize line TC1507 genomic DNA, isolated, and sequenced for Region 1 through Region 6, shown in Table 1, and confirmed to match the previously determined sequence from the Spe I fragment, described in Example 2. However, the region from bp 2358 to by 2829, immediately adjacent and 5' to the start of the full-length insert was recalcitrant to PCR amplification and appeared to be larger than the sequence obtained from the Spe I clone described above. The use of primer pairs flanking this region and the Advantage®-GC 2 Polymerase Mix (BD Biosciences Clontech, Palo Alto, Calif.) was successful in amplifying PCR products from B.t. Cry1F maize line TC1507 genomic DNA for sequencing. The amplification conditions used to produce amplicons with the Advantage®-GC 2 system are shown in Table 10. The DNA primer pairs used to confirm the sequence in the region from by 2358 to 2829 are those listed in SEQ ID NO:1 and SEQ ID NO:2; and SEQ ID NO:2 and SEQ ID NO:23. Sequence from this region is described in Table 1 (Regions 7a, 7b, 7c, and 8).

Example 4. Event TC1507 5' Flanking Sequence

A Description of Each Region is Provided in Table 1.

Region 1 (SEQ ID NO: 28) Maize genomic (no significant homology)

1 ACTAGTTTCC TAGCCCGCGT CGTGCCCCTA CCCCACCGAC GTTTATGGAA

51 GGTGCCATTC CACGGTTCTT CGTGGCCGCC CCTAAGGATG TAAATGGTCG

101 GTAAAATCCG GTAAATTTCC GGTACCGTTT ACCAGATTTT TCCAGCCGTT

151 TTCGGATTTA TCGGGATATA CAGAAAACGA GACGGAAACG GAATAGGTTT

201 TTTTTCGAAA ACGGTACGGT AAACGGTGAG ACAAACTTAC CGTCCGTTTT

251 CGTATTTCTC GGGAAACTCT GGTATATTCC CGTATTTTCC

301 CCGACCCACG GACCTGCCAA TCAACCATCA GCCAGTCAGC CCATCCCCAC

351 AGCTATGGCC CATGGGGCCA TGTTGGCCAC ATGCCCACGC AACGCAAGGC

401 AGTAAGGCTG GCAGCCTGGC ACGCATTGAC GCATGTGGAC ACACACAGCC

451 GCCGCCTGTT CGTGTTTCTG TGCCGTTGTG CGAGACTGTG ACTGCGAGTG

501 GCGGAGTCGG CGAACGGCGA GGCGTCCCG GAGTCTGGAC TGCGGCTGTG

551 GACAGCGACG CTGTGACGGC GACTCGGCGA AGCCCCAAGC TACCAAGCCC

601 CCAAGTCCCC ATCCATCTCT GCTTCTCTGG TCATCTCCTT CCCCTGGTCG

#### -continued

- Region 2 (SEQ ID NO: 29) Undescribed maize genomic sequence (complement) 670 G CCGAAGCATC ACGAAACGCA CTAAGACCTC

  701 GAAGGAGTCA AACCACTCCT CCGAGGCCTC GGGGGCTACA CCCGGCGGT
  - 801 CCTTGCAAAA AAAGTGCGAC AAAAGCCTCC AAGCGAGTAT TAACACTCAC

751 GCGCTCGCGC GCACCCACCG GAACAAAATG TAACCGAGAA AGGTCGGTCC

- 851 TTTGAGGCTC GGGGGCTAC
- Region 3 (SEQ ID NO: 30) Fragment of maize Huck-1 retrotransposon 870 T GTCGGGGACC ATAATTAGGG GTACCCCCAA
  - 901 GACTCCTAAT CTCAGCTGGT AACCCCCATC AGCACAAAGC TGCAAAGGCC
  - 951 TGATGGGTGC GATTAAGTCA AGGCTCGGTC CACTCAAGGG ACACGATCTC
- 1001 GCCTCGCCCG AGCCCAGCCT CGGGCAAGGG CGGCCGACCC CGAGGATTCA
- 1051 CGTCTCGCCC GAGGGCCCCC TCAAGCGACG GGCACACCTT CGGCTCGCCC
- 1101 GAGGCCCATT CTTCGCCGAG AAGCAACCTT GGCCAGATCG CCACACCGAC
- 1151 CGACCGTATC GCAGGAGCAT TTAATGCGAG GATCGCCTGA CACCTTATCC
- 1201 TGACGCGCC TCTTCAGTCG ACAGAGCCGA AGTGACCGCA ATCACTTCGC
- 1251 CGCTCCACTG ACCGACCTGA CAAGAAGACA GCGCCGCCTG CGTCGCTCCG
- 1301 ACTGCTGTGC CACTCGACAG AGTGAGGCTG ACAGCAGCCA AGTCCGGCCT
- 1351 CGGGCGCCAT AGGAAGCTCC GCCTCGCCCG ACCCTAGGGC TCGGACTCGG
- 1401 CCTCGGCTCC GGAAGACGAC GAACTACGCT TCGCCCGACC CCAGGGCTTG
- 1451 GACTCAGCCT CGGCTCCGGA AGACGACGAA TTCCGCCTCG CCCGACCCCA

1551 GACCCCAGGG CTCGGACTCA GCCTCGGCTC CGGAAGACGA CGAACTCCGC

1501 GGGCTCGGAC TCGGCCTCGG CTCCAGAAGA CGACGAACTC CGCCTCGCCC

- 1601 CTCGCCCGAC CCCAGGGCTC GGACTCAGCC TCGGCCTCAG ACGATGGTCT
- 1651 CCGCCTCGCC CGACCCGGGG CTCGGACTCG A
- Region 4 (SEQ ID NO: 31) Fragment of cry1F gene 1682 CCTTTCTAT CGGACCTTGT
- 1701 CAGATCCTGT CTTCGTCCGA GGAGGCTTTG GCAATCCTCA CTATGTACTC
- 1751 GGTCTTAGGG GAGTGGCCTT TCAACAAACT GGTACGAATC ACACCCGCAC
- 1801 ATTCAGGAAC TCCGGGACCA TTGACTCTCT AGATGAGATA CCACCTCAAG
- 1851 ACAACAGCGG CGCACCTTGG AATGACTACT CCCATGTGCT GAATCATGTT

1951 ACCAATGTTC TCTTGGACGC ATCGTAGCGC TACCCCCACA AACACCATTG

1901 ACCTTTGTGC GCTGGCCAGG TGAGATCTCA GGTTCCGACT CATGGAGAGC

- 2001 ATCCAGAGAG AATCAC
- Region 5 (SEQ ID NO: 32) Fragment of maize chloroplast rpoC2 gene 2017 TCAT TCTTCAAGAA CTGCATATCT TGCCGAGATC
- 2051 CTCATCCCTA AAGGTACTTG ACAATAGTAT TATTGGAGTC GATACACAAC
- 2101 TCACAAAAA TACAAGAAGT CGACTAGGTG GATTGGTCCG AGTGAAGAGA
- 2151 AAAAAAAGCC ATACAGAACT CAAAATCTTT TCCGGAGATA TTCATTTTCC
- 2201 TGAAGAGCG GATAAGATAT TAGGTGGCAG TTTGATACCA CCAGAAAGAG

#### -continued

2251 AAAAAAAAGA TTCTAAGGAA TCAAAAAAAA GGAAAAATTG
GGTTTATGTT

2301 CAACGGAAAA AATTTCTCAA AAGCAAGGAA AAGTATT

Region 6 (SEQ ID NO: 33) Fragment of maize chloroplast or ubiZM1(2) promoter 2338 GTG GCTATTTATC

2351 TATC

Nucleotides 2355-2358 (CGT) connect Region 6 to Region 7a. Region 7a (SEQ ID NO: 34) Fragment of pat gene 2358 GCA GCTGATATGG CCGCGGTTTG TGATATCGTT AACCATTACA

2401 TTGAGACGTC TACAGTGAAC TTTAGGACAG AGCCACAAAC ACCACAAGAG

2451 TGGATTGATG ATCTAGAGAG GTTGCAAGAT AGATACCCTT GGTTGGTTGC

2501 TGAGGTTGAG GGTGTTGTGG CTGGTATTGC TTACGCTGGG CCCTGGAAGG

2551 CTAGGAAC

Region 7b (SEQ ID NO: 35) Fragment of pat gene (complement) 2559 CC TCAACCTCAG CAACCAACCA ATGGTATCTA TCTTGCAACC

2601 TCTCTAGATC ATCAATCCAC TCTTGTGGTG TTTGTGGCTC TGTCCTAAAG

2651 TTCACTGTAG ACGTCTCAAT GTAATGGTTA ACGATATCAC AAACCG

Region 7c (SEQ ID NO: 36) Fragment of cry1F gene (complement) 2697 AGAG

2701 AAGAGGGATC T

Region 8 (SEQ ID NO: 37) Fragment of Polylinker 2712 CGAAGCTTC GGCCGGGCC CATCGATATC CGCGGCCATG

2751 CCTGCAGTGC AGCGTGACCC GGTCGTGCCC CTCTCTAGAG ATAATGAGCA

2801 TTGCATGTCT AAGTTATAAA AAATTACCA

Region 9 (SEQ ID NO: 25) Full-length insert of PHI8999A

# Example 5. Description of the Flanking Sequence 5' to the Insert in Maize Event TC1507

In order to more fully describe the event TC1507 5' flanking sequence, homology searching was done against the GenBank public databases (release 122, 2/01) using the Basic Local Alignment Search Tool (BLAST). The BLAST program performs sequence similarity searching and is particularly useful for identifying homologs to an unknown sequence. In addition to searching the public databases, pairwise alignments were performed using AlignX (Infor-Max Inc., Bethesda, Md.) to look for homology between the maize event TC1507 flanking sequence and the PHI8999A transgenic insert. The results of these homology searches are presented in Table 1. The TC1507 5' flanking sequence is numbered with base 1 being the furthest 5' to the insert and base 2830 at the starting point of the full-length PHI8999A transgenic insert (see FIG. 1). The percent identity values

indicate the percentage of identical matches across the length of the sequences analyzed.

In most cases, similarity searching with the event TC1507 5' flanking sequence resulted in a match to one unique sequence based on a very high percent identity value. Those sequences are identified in Table 1. In addition, there are two regions in the TC1507 5' DNA flanking sequence with high similarity to more than one known sequence. In regions 870-1681 and 2338-2354, the percent identity scores with both sequence fragments are sufficiently high that a single match (homolog) cannot be determined. The two possible homologs for each of these regions are indicated in Table 1.

Highly similar sequences were identified for all but the first 669 base pairs of sequence. Generally, the results of similarity searching indicate high homology with maize genomic sequences 5' to base 1681. The region from base 1682 to the start of the PHI8999A insert at position 2830 contains some fragments associated with the transformation event.

## TABLE 1

| Sequence summary for event TC1507 insert |                                 |         |            |         |                                 |                                  |
|--|---------------------------------|---------|------------|---------|---------------------------------|----------------------------------|
| Region                                   | Location<br>in SEQ<br>ID NO: 24 | Size bp | % Identity | Homolog | Location in homologous sequence | Description                      |
| 1  | 1-669                           | 669     | $N/A^1$    | N/A     | N/A                             | No significant homology detected |

TABLE 1-continued

|        | TABLE 1-continued |         |            |                 |                               |                                       |  |  |
|--------|-------------------|---------|------------|-----------------|-------------------------------|---------------------------------------|--|--|
|        |                   | Sequ    | ence summa | ary for event T | C1507 insert                  |                                       |  |  |
|        | Location in SEQ   |         |            |                 | Location in homologous        |                                       |  |  |
| Region | ID NO: 24         | Size bp | % Identity | Homolog         | sequence                      | Description                           |  |  |
| 2      | 670-869           | 200     | 90.5       | AF123535        | 52432-52632                   | Undescribed                           |  |  |
|        |                   |         |            |                 | (complement)                  | maize genomic sequence                |  |  |
| 3      | 870-1681          | 812     | 89.4       | AF050439        | 1-801                         | Fragment of                           |  |  |
|        |                   |         |            |                 |                               | maize Huck-1 retrotransposon          |  |  |
|        |                   |         | 966        | AF050438        | 1 707                         | 5' LTR <sup>2</sup> Fragment of maiga |  |  |
|        |                   |         | 86.6       | Ar030436        | 1-797                         | Fragment of maize Huck-1 retro-       |  |  |
| 4      | 1682-2016         | 335     | 100.0      | PHI8999A        | 3149-3483                     | transposon 3' LTR<br>Fragment of      |  |  |
| 7      | 1002-2010         | 333     | 100.0      | 111102227       |                               | cry1F gene                            |  |  |
| 5      | 2017-2337         | 321     | 100.0      | X86563          | 29429-29749                   | Fragment of maize chloroplast         |  |  |
|        |                   |         |            |                 |                               | rpoC2 gene                            |  |  |
|        |                   |         |            |                 |                               | (RNA polymerase beta-2 subunit)       |  |  |
| 6      | 2338-2354         | 17      | 100.0      | X86563          | 97643-97659                   | Fragment of                           |  |  |
|        |                   |         |            |                 |                               | maize chloroplast<br>trnI gene (tRNA- |  |  |
|        |                   |         | 82.4       | PHI8999A        | 182-197                       | Ile)                                  |  |  |
|        |                   |         | 02.4       | гптоэээд        | 162-197                       | Fragment of maize ubiZM1(2)           |  |  |
| 7a     | 2358-2558         | 201     | 100.0      | PHI8999A        | 5320-5475                     | promoter<br>Fragment of pat           |  |  |
|        |                   |         | 100.0      |                 |                               | gene                                  |  |  |
| 7b     | 2559-2696         | 138     | 99         | PHI8999A        | 5336-5518<br>(complement)     | Fragment of pat gene                  |  |  |
| 7c     | 2697-2711         | 15      | 100.0      | PHI8999A        | 2544-2558                     | Fragment of                           |  |  |
| 8      | 2712-2829         | 118     | 100.0      | PHI8999A        | (complement)<br>36-153        | cry1F gene<br>Fragment of             |  |  |
|        |                   |         |            |                 |                               | polylinker region                     |  |  |
|        |                   |         |            |                 |                               | (bases 36-80) and ubiZM1(2)           |  |  |
|        |                   |         |            |                 |                               | promoter (bases 81-153)               |  |  |
| 9      | 2830-9015         | 6186    | 100.0      | PHI8999A        | 11-6196                       | Full-length insert                    |  |  |
| 10     | 9016-9565         | 550     | 100.0      | PHI8999A        | 3906-4456                     | of PHI8999A<br>Inverted ORF25         |  |  |
|        |                   |         |            |                 | (complement)                  | terminator                            |  |  |
| 11     | 9566-9693         | 128     | 100.0      | NC_001666       | 121851-121978 (complement) &  | Fragment of maize chloroplast         |  |  |
|        |                   |         |            |                 | 100759-100886                 | rps12 rRNA (23S                       |  |  |
| 12     | 9696-10087        | 392     | 99         | NC_001666       | 17091-17483                   | ribosomal RNA) Fragment of            |  |  |
|        |                   |         |            |                 | (complement)                  | maize chloroplast                     |  |  |
| 13     | 10088-10275       | 188     | 99         | PHI8999A        | 5333-5520                     | genome<br>Fragment of pat             |  |  |
| 14     | 10278-10358       | 81      | 100        | NC_001666       | (complement)<br>137122-137202 | gene<br>Ergament of                   |  |  |
| 14     | 10276-10336       | 61      | 100        | NC_001000       | (complement)                  | Fragment of maize chloroplast         |  |  |
|        |                   |         |            |                 |                               | "ORF241" -<br>hypothetical            |  |  |
|        |                   |         | . •        |                 |                               | protein gene                          |  |  |
| 15     | 10359-10612       | 254     | $N/A^1$    | N/A             | N/A                           | No significant homology               |  |  |
|        | 40045             |         | <b>.</b> 1 | <b>3.</b> T ( ) | <b>-</b> - / -                | detected                              |  |  |
| 16     | 10613-11361       | 749     | $N/A^1$    | N/A             | N/A                           | No description available              |  |  |
|        |                   |         |            |                 |                               |                                       |  |  |

<sup>1</sup>N/A; not applicable

<sup>2</sup>LTR; long terminal repeat

Example 6. Confirmation of the Presence of Regions 1, 2, and 3 in an Unmodified Control Corn Line

PCR analysis was used to determine if Regions 1, 2, and 3 (Table 1) in the 5' flanking region of Event TC1507 are formation to produce maize event TC1507 and thus represents a border with corn genomic DNA. Nine different PCR

analyses were carried out on genomic DNA prepared from <sub>60</sub> TC1507 and the unmodified control corn line Hi-II (see Armstrong (1994) The Maize Handbook, ed. Freeling and Walbot, Springer-Verlag, New York, pp. 663-671, for information on Hi-II) as outlined in Table 2 using the primer sequences shown in Table 3. Two reactions were designed to present in an unmodified control corn line used for trans- 65 amplify DNA within Region 1 of the 5' flanking region from by 25 to 324 (Reaction A—300 bp amplicon); and from by 25 to 480 (Reaction B—456 bp amplicon). The expected

amplicons were present in both the Hi-II unmodified corn line and in maize event TC1507. One PCR primer pair, Reaction C, spanned Region 2 to Region 3 of the 5' flanking region from by 759 to 1182 (424 bp amplicon) and again produced PCR products of the expected size in both Hi-II 5 and TC1507. Reaction D, spanned Region 1 to Region 3 of the 5' flanking region from by 415 to 1182 (768 bp amplicon) and again produced PCR products of the expected size in both Hi-II and TC1507. Reactions E and F were designed as specific primer pairs for the pat gene region of the full-length  $_{10}$ insert of PHI8999A in TC1507 and thus an amplicon in the unmodified Hi-II corn line is not expected. The results indicate that both Reactions E and F are specific for a maize line transformed with a pat gene region and produce the expected amplicon, whereas no amplicon was produced in 15 the unmodified Hi-II corn line. Reaction G was also designed as a primer pair that would produce an amplicon of 366 bp in the maize event TC1507 and no amplicon in the unmodified Hi-II corn line.

28

Reactions H and I were designed as specific primer pairs for TC1507 that would span the end of the transgenic insert into the 5' flanking region. In both Reactions H and I, the reverse primer was located in the ubiquitin promoter region of the full-length PHI8999A insert (Region 9 in Table 1) and the forward primer was located in Region 5, the rpoC2 gene fragment (see Table 1). Reaction H and Reaction I both produced an amplicon in maize line TC1507 and did not produce an amplicon in the unmodified control corn line. These results indicate that both Reactions H and I are specific for the TC1507 event.

The PCR results show that the undescribed sequence (Region 1) is present in the unmodified corn line Hi-II and that Regions 1, 2 and 3, are contiguous in the unmodified corn line Hi-II. The DNA sequences amplified in Reactions A, B, C, and D are not unique to the 5' flanking region of maize event TC1507 but are also present in the unmodified corn line Hi-II.

TABLE 2

| PCR reactions for sequence 5' to the PHI8999A insert in maize event TC1507      |
|---|
| and for regions within the full-length insert of PHI8999A in maize event TC1507 |

| Reaction                                     | PCR<br>Amplicon<br>Location   | Amplicon<br>Size (bp) | Region in TC1507<br>flanking sequence<br>or PHI8999A<br>insert                                  | Amplicon<br>present<br>In Hi-II | Amplicon present in maize line TC1507 |
|--|---|-----------------------|---|---------------------------------|---------------------------------------|
| A  | 25-324 bp<br>in TC1507<br>flanking  | 300                   | Region 1  | Yes                             | Yes                                   |
| В  | sequence<br>25-480 bp<br>in TC1507<br>flanking<br>sequence  | 456                   | Region 1  | Yes                             | Yes                                   |
|  | bp in TC1507<br>flanking<br>sequence  |                       | 3   |                                 |                                       |
| D  | bp in TC1507 5' flanking sequence   | 768                   | Region 1 to<br>Region3  | Yes                             | Yes                                   |
| E<br>Not Unique<br>to TC1507                 | 4750-5794<br>bp<br>in PHI8999A  | 1045                  | Region 9 (in full-<br>length insert of<br>PHI8999A 35S<br>promoter to pat<br>gene)              | No                              | Yes                                   |
| F<br>Not Unique<br>to TC1507                 | 4827-5308<br>bp<br>in PHI8999A  | 482                   | Region 9 (in full-length insert of PHI8999A 35S promoter to pat gene)                           | No                              | Yes                                   |
| Detects cry1F fragment in 5' flanking region | cry1F sequence in 5' flanking sequence and in full-length insert of PHI8999A                      | 366                   | Spans 335 bp cry1F sequence in 5' flanking sequence and same sequence in the full-length insert | No                              | Yes                                   |
| H<br>Unique to<br>TC1507                     | 2158 bp in Region 5 (rpoC2 gene fragment) to 3069 bp in Region 9 (full-length insert of PHI8999A) | 912                   | Region 5 to Region 9 Unique to Insertion Event [SPANS UNIQUE JUNCTION REGIONS]                  | No                              | Yes                                   |

TABLE 2-continued

| PCR reactions for sequence 5' to the PHI8999A insert in maize event TC1507      |
|---|
| and for regions within the full-length insert of PHI8999A in maize event TC1507 |

| Reaction                 | PCR<br>Amplicon<br>Location   | Amplicon<br>Size (bp) | Region in TC1507 flanking sequence or PHI8999A insert                          | Amplicon<br>present<br>In Hi-II | Amplicon present in maize line TC1507 |
|--------------------------|---|-----------------------|--|---------------------------------|---------------------------------------|
| I<br>Unique to<br>TC1507 | 2158 bp in Region 5 (rpoC2 gene fragment) to 3001 bp in Region 9 (full-length insert of PHI8999A) | 844                   | Region 5 to Region 9 Unique to Insertion Event [SPANS UNIQUE JUNCTION REGIONS] | No                              | Yes                                   |

TABLE 3

PCR primers for sequence 5' to the PHI8999A insert in TC1507 and for regions within the full-length insert of PHI8999A in maize event TC1507

| Reaction Am | plicon Size | (bp) Prim | er | Pair | î. | Primer          | Sequences                | 5'   | to   | 3 ' |
|-------------|-------------|-----------|----|------|----|-----------------|--------------------------|------|------|-----|
| A           | 300         |           |    |      |    |                 | CCCCACCGAC               |      |      |     |
| В           | 456         | ~         |    |      |    |                 | CCCCACCGAC               |      |      |     |
| С           | 424         |           |    |      |    |                 | CCACCGGAAC<br>CATTAAATGC |      |      |     |
| D           | 768         | ~         |    |      |    |                 | ACGCATTGAC               |      |      |     |
| E           | 1045        | ~         |    |      |    |                 | ACCTAACAGA<br>GCAACTCAAA |      |      |     |
| F           | 482         | ~         |    |      |    |                 | TTCGTCAACA               |      |      |     |
| G           | 366         | ~         |    |      |    |                 | GACTCGACCI<br>CTTGAAGAAI |      |      |     |
| Н           | 912         | SEQ       | ID | NO:  | 1  | GTAGTAC<br>AGAG | CTATAGATTA               | TAT' | TATI | CGT |
|             |             | SEQ       | ID | NO:  | 2  |                 | CAGAACTCAA               | AAT  | CTTI | TCC |
| I           | 844         | ~         |    |      |    | GGAG            | CAGAACTCAA               |      |      | TCC |
|             |             | SEQ       | ID | NO:  | 23 | CTTCAAA         | ACAAGTGTGA               | CAA  | A    |     |

Example 7. Flanking Sequence 3' to Inserted Transgenic DNA in Maize Event TC1507

Two separate PCR approaches were used to extend the length of the sequence information 3' to the full-length PHI8999A insert in maize event TC1507. In the first approach PCR primer pairs were designed to amplify a product that spanned the junction between the full-length depiction of the inverted ORF25 terminator. A forward primer was located at the end of the full-length PHI8999A insert and a series of reverse primers were located at 100 bp intervals in the inverted sequence. In this manner the length of the inverted fragment present in the maize event TC1507 65 could be determined within a 100 bp region based on the successful PCR reactions. This method indicated the

inverted fragment contained the majority of the ORF25 terminator but no Cry1F sequence. PCR fragments were isolated and sequenced from this region.

In the second approach PCR primers were designed to walk out into the flanking DNA sequence from the inverted ORF25 terminator region as determined in the PCR experiment described above. Genomic DNA isolated from two to insert and the inverted ORF25 terminator, see FIG. 1 for a 60 three individual plants of event TC1507 and an unmodified control corn line was digested with various restriction enzymes and then ligated to adaptors specific for the restriction enzyme used for digestion (Universal Genome Walker<sup>TM</sup> Kit, Clontech Laboratories, Inc. and Devon et al. (1995) Nucleic Acids Res. 23:1644-1645). Primary PCR was carried out using an ORF25 terminator specific primer and a primer homologous to the adaptor sequence ligated onto

the digested DNA. In order to increase the specificity of the reaction a nested secondary PCR was performed again with another ORF25 terminator specific primer and a secondary primer homologous to the adaptor sequence with the secondary primers being internal to the respective primers used 5 in the primary PCR. Products produced by the nested PCR were analyzed by agarose gel electrophoresis and fragments unique to TC1507 DNA samples were isolated and sequenced. Fragments were amplified from both the ORF25 terminator contained within the full-length insert and from the targeted (inverted) ORF25 terminator on the 3' end of the full-length PHI8999A insert. Fragments from the full-length insert were of a predicted size based on the knowledge of the restriction enzyme sites located in the full-length insert. 15 available for Region 16. Fragments produced from the 3' inverted ORF25 terminator appeared as fragments of unexpected size. Sequence analysis of amplified fragments from the 3' inverted ORF25 terminator resulted in flanking DNA sequence of 1043 bp. Resultant sequence from the above series of genome walk- region is not unique to the 3' flanking region of maize event ing experiments was used to design additional primers to walk further out from the insert into the bordering maize genome with a final 3' flanking sequence, of 2346 bp.

In order to describe the TC1507 3' flanking sequence, homology searching was done against the GenBank public 25 databases using the Basic Local Alignment Search Tool (BLAST). The BLAST program performs sequence similarity searching and is particularly useful for identifying homologs to an unknown sequence. In addition to searching

the public databases, alignments were performed using SeqMan 4.05<sup>TM</sup>, Martinez and Needleman-Wunsch alignment algorithms (DNASTAR Inc.) to look for homology between the TC1507 3' flanking sequence and the PHI8999A transgenic insert. The results of these homology searches are presented in Table 1. The percent identity values indicate the percentage of identical matches across the length of the sequences analyzed. The results of similarity searching for the 3' flanking sequence indicate high homology with three regions of maize chloroplast DNA, a 188 bp fragment of the pat gene, and 254 bp of DNA (Region 15, Table 1) with no significant homology. An additional 749 bp (Region 16) beyond Region 15 (see Table 1) was also sequenced. No similarity searching results are

PCR analysis on control and TC1507 genomic DNA determined that the 254 bp sequence (Region 15, fragment of maize chloroplast "ORF241") is present in the maize genome. The DNA sequence of Region 15 in the 3' flanking TC1507 but is also present in the unmodified control corn line. The TC1507 3' flanking sequence is presented in Example 8 and diagrammed in FIG. 1.

Example 8. Sequence of the Region 3' to the End of the Full-Length Insert DNA in Maize Event TC1507

A description of each region is in Table 1.

```
CTCAC TCCGCTTGAT CTTGGCAAAG ATATTTGACG
        CATTTATTAG TATGTGTTAA TTTTCATTTG CAGTGCAGTA TTTTCTATTC
  9051
      GATCTTTATG TAATTCGTTA CAATTAATAA ATATTCAAAT CAGATTATTG
  9151 ACTGTCATTT GTATCAAATC GTGTTTAATG GATATTTTTA TTATAATATT
  9201 GATGATATCT CAATCAAAAC GTAGATAATA ATAATATTTA TTTAATATTT
  9251 TTGCGTCGCA CAGTGAAAAT CTATATGAGA TTACAAAATA CCGACAACAT
  9301 TATTTAAGAA ACATAGACAT TAACCCTGAG ACTGTTGGAC ATCAACGGGT
  9351 AGATTCCTTC ATGCATAGCA CCTCATTCTT GGGGACAAAA GCACGGTTTG
  9401 GCCGTTCCAT TGCTGCACGA ACGAGCTTTG CTATATCCTC GGGTTGGATC
  9451 ATCTCATCAG GTCCAATCAA ATTTGTCCAA GAACTCATGT TAGTCGCAAC
  9501 GAAACCGGGG CATATGTCGG GTATCTCGAG CTCGCGAAAG CTTGGCTGCA
  9551 GGTCGACGGA TCCTT
Region 11 (SEQ ID NO: 39) Fragment of maize chloroplast rps12 rRNA gene (comple-
ment)
  9566 CAACA AAAGGGTACC TGTACCCGAA ACCGACACAG
  9601 GTGGGTAGGT AGAGAATACC TAGGGGCGCG AGACAACTCT CTCTAAGGAA
  9651 CTCGGCAAAA TAGCCCCGTA ACTTCGGGAG AAGGGGTGCC CCC
Nucleotides 9694-9695 (CG) connect Region 11 to Region 12.
Region 12 (SEQ ID NO: 40) Fragment of maize chloroplast genome
  9696 CTAAC
  9701 AATAAACGAA TACGGTTTAT GTATGGATTC CGGTAAAATA CCGGTACTCG
  9751 ATTTCATAAG AGTCGAATAG GAAGTTAAGA TGAGGGTGGT ATCATCATAA
  9801 AAATGGAGTA GTATCCTAAA TTATACTAAT CCACGTATGA TATGTATGCC
  9851 TTTCCTTATC AACCGGAAGT AGTGCAAAAA AAATTCTATA CTGCACTGCT
```

Region 10 (SEQ ID NO: 38) Fragment of ORF25 Terminator (complement)

|      |   |     | _   |
|------|---|-----|-----|
| -con | T | ın: | บอด |

9901 CTCTTTTTAC TGAGAAATGC AAAAAAATAA AAGTGAAGTA AGGGTGCCCC
9951 ATAGATATTT GATCTTGCCT CCTGTCCCCC CCCCCTTTT TTCATCAAAA
10001 ATTTCCATGA AAAAAGAAAA GATGAATTTG TCCATTCATT GAACCCTAGT

10051 TCGGGACTGA CGGGGCTCGA ACCCGCAGCT TCCGCCT

Region 13 (SEQ ID NO: 41) Fragment of pat gene (complement) 10088 GTT CCTAGCCTTC

10101 CAGGGCCCAG CGTAAGCAAT ACCAGCCACA GCACCCTCAA

10201 GTGGTGTTTG TGGCTCTGTC CTAAAGTTCA CTGTAGACGT CTCAATGTAA

10151 CAACCAAGGG TATCTATCTT GCAACCTCTC TAGATCATCA ATCCACTCTT

10251 TGGTTAACGA TATCACAAAC CGCGG

CCTCAGCAAC

Nucleotides 10276-10277 (AA) connect Region 13 to Region 14. Region 14 (SEQ ID NO: 42) Fragment of maize chloroplast ORF241 (complement)

10278 CAC AAGAACGAAA GCACCTTTTC

10301 ATTCTTTCAT ATACTAGGGG TTTTTACTTG GAAAAGACAA TGTTCCATAC

10351 TAAAGGAT

Region 15 (SEQ ID NO: 43) Maize genomic (no significant homology) 10359 AG CTGCAGAAGC CGCCACCGTC TTGAGGACCT TCCGGGGAGC

10401 CAGACCGGTC GAACCGTGCC TCCACTTGCT AAGGAGAAAG

GGAAAATCAG

10451 GGCCAGGACA TACGAAGGAG GAGCCAGAAC GAAGATATCC

TAAGATACTT

10501 ACTCGCTCCG GGCCATGATC AATCATGCCT GTGGGGAGGT CTCTCGCACC

10551 TCGATCCATG AAGGTACCAC CGAGGTCTGC CCCGCCGCCG GCTTCGGTAC

10601 CGTCCTCGCC TT

Region 16 (SEQ ID NO: 44) Maize genomic 10613 GGGCGCCC GAGGCACCCG GGGGATGGAC TGCCCAGGCG

10651 CAGCCACGAC GACCCAAGGA TCACCCTCCT GCGCAGTCGG

CACGAGCAAT

10701 AGTTCTCGGG GAACAGGCAG CTTGGCCTGA CTCCCCGGGG TCACCTCAAC

10751 TACCTCGGCC GAGGGGTCAA GTACCCCCTC AGTCCGCCCC CGCTCTTCGG

10801 ACCGGGACCC CGACGTCCCG GCCCCGGATA CCGACGGCAC CAGCCCGCTC

10851 GGGGGCTGGC TTGACGACCC CTGGCCCAGC CTCAGATCTG GGCTGAGGCC

10901 GAGGCAGGCG GCCATGTCGT CGTCTTCATC ATCGTCTTCA TCATCGTCGT

11001 TGCCGACGAC GAAGCCTTTC CAAGGCATCC CGAGCCCACG TCCGCTCGTG

10951 CGTCATCAGG CGTCTCCGGC GACGGCTCCC TTGGGAGCCC CTCCCTCTCC

11051 GGCCCGAGCC TTCTTTGCGT CCTTCTTCTC CTTCCTCTTC TCCGCGGTGA

11101 CCCTCCGCGC AGCTCGGTCC ACCGCATCCT CCGGGACTGG TGGCAGGGAA

11151 GGCTTGTGAT GCCCTACCTC CTGGAGACAG ACGAAAAGTC TCAGCTATGA

11201 GAACCGAGGG CAATCTGACG CAAGAAGGAA GAAGGAGCGG

ATACTCACCA

#### -continued

11251 GAGACACGCA CCCGCGATCG GGACGCATTA AGGGCTGGGA

35

AAAAGTGCCG

11301 GCCTCTAATT TCGCTACCGT GCCGTCCACC CACCTGTGGA GGTCATCGAT

11351 GGGAAGGGGA A

Example 9. Confirmation of the Presence of Region 15 in the Unmodified Control Corn Line

PCR analysis was used to determine if the undescribed region of sequence on the end of the 3' flanking sequence (Region 15 in Table 1) is present in the unmodified control corn line used for transformation to produce maize event TC1507 and thus represents a border with corn genomic DNA. Successful PCR amplification of Region 15 in both maize line TC1507 and the unmodified Hi-II control corn line revealed that Region 15 was indeed present in corn genomic DNA. Five different PCR analyses were carried out on genomic DNA prepared from TC1507 and the unmodified Hi-II control corn line as outlined in Table 7 below using the primer sequences shown in Table 8. Three reactions were designed to amplify DNA within Region 15 of the 3' flanking region; Reaction L—producing a 175 bp amplicon, 30 Reaction M—producing a 134 bp amplicon, and Reaction N—producing a 107 bp amplicon. The expected amplicons were present in both the unmodified control corn line and in maize line TC1507. Reactions J and K were designed as specific primer pairs for TC1507 that would span the end of <sup>35</sup> the insert into the 3' flanking region. In Reaction J, the forward primer was located in the pat gene fragment on the 3' end of the full-length PHI8999A insert (Region 13 in Table 1) and the reverse primer was located in the undefined 40 Region 15. In Reaction K the forward primer was located in the chloroplast hypothetical protein gene on the 3' end of the full-length insert (Region 14 in Table 1) and the reverse primer was located in the undefined Region 15. Both Reac-

tion J and Reaction K produced an amplicon in maize line TC1507 and did not produce an amplicon in the unmodified control corn line. The results indicate that both Reactions J and K are specific for the TC1507 event.

**36** 

The PCR results indicate that the undescribed sequence (Region 15) of the 3' flanking sequence of TC1507 is also present in genomic DNA isolated from the unmodified Hi-II control corn line. The DNA sequences amplified in Reactions L, M, and N are not unique to the 3' flanking region of TC1507 but are also present in the unmodified control corn line.

TABLE 7

|   | PCR reactions for sequence 3' to the PHI8999A insert |           |   |            |                                |  |  |
|---|--|-----------|---|------------|--------------------------------|--|--|
|   | in maize event TC1507                                |           |   |            |                                |  |  |
| ) | Danatian   | -         | Region in TC1507 3' flanking                    | 1          | Amplicon present in maize line |  |  |
|   | Reaction   | Size (bp) | sequence  | in Control | TC1507                         |  |  |
| 5 | J  | 342       | Region 13 (pat gene fragment) to Region 15      | No         | Yes                            |  |  |
| ) | K  | 252       | Region 14<br>(chloroplast gene) to<br>Region 15 | No         | Yes                            |  |  |
|   | L  | 175       | Region 15                                       | Yes        | Yes                            |  |  |
|   | M  | 134       | Region 15                                       | Yes        | Yes                            |  |  |
|   | $\mathbf{N}$   | 107       | Region 15                                       | Yes        | Yes                            |  |  |

TABLE 8

|              | _   |        | _      | e 3' to the<br>event TC1507   |
|--------------|-----|--------|--------|-------------------------------|
| Reaction Amp |     |        |        | Primer Sequences 5' to 3'     |
| J            | 342 | SEQ ID | NO: 3  | TGTGGTGTTTGTGGCTCTGTCCT<br>AA |
|              |     | SEQ ID | NO: 5  | GACCTCCCCACAGGCATGATTG<br>ATC |
| K            | 252 | SEQ ID | NO: 4  | AGCACCTTTTCATTCTTTCATAT<br>AC |
|              |     | SEQ ID | NO: 5  | GACCTCCCCACAGGCATGATTG<br>ATC |
| L            | 175 | SEQ ID | NO: 16 | AAGCCGCCACCGTCTTGAGGAC<br>CTT |
|              |     | SEQ ID | NO: 5  | GACCTCCCCACAGGCATGATTG<br>ATC |
| M            | 134 | SEQ ID | NO: 17 | GTCGAACCGTGCCTCCACTTGC<br>TAA |

TABLE 8-continued

| PCR primers for sequence 3' to the<br>PHI8999A insert in maize event TC1507 |              |             |        |                               |  |
|---|--------------|-------------|--------|-------------------------------|--|
| Reaction An   | mplicon Size | (bp) Primer | Pair   | Primer Sequences 5' to 3'     |  |
|   |              | SEQ ID      | NO: 5  | GACCTCCCCACAGGCATGATTG<br>ATC |  |
| N   | 107          | SEQ ID      | NO: 18 | AGAAAGGGAAAATCAGGGCCA<br>GGAC |  |
|   |              | SEQ ID      | NO: 5  | GACCTCCCCACAGGCATGATTG<br>ATC |  |

#### Example 10. PCR Primers

DNA event specific primer pairs were used to produce an amplicon diagnostic for TC1507. These event primer pairs include, but are not limited to, SEQ ID NO: 1 and SEQ ID NO: 2; SEQ ID NO: 2 and SEQ ID NO: 23; SEQ ID NO: 20 3 and SEQ ID NO: 5; and SEQ ID NO: 4 and SEQ ID NO: 5. In addition to these primer pairs, any primer pair derived from SEQ ID NO: 26 and SEQ ID NO: 27 that when used in a DNA amplification reaction produces a DNA amplicon diagnostic for TC1507 is an aspect of the present invention. 25 The amplification conditions for this analysis are illustrated in Table 9, however, any modification of these methods that use DNA primers or complements thereof to produce an amplicon DNA molecule diagnostic for TC1507 is within the ordinary skill of the art. The preferred amplification 30 conditions for reactions utilizing the PCR primers identified in SEQ ID NOS: 1, 2, and 23 are illustrated in Table 10. In addition, control primer pairs, which include SEQ ID NOS: 10 and 11; SEQ ID NOS: 10 and 12; SEQ ID NOS: 13 and 14; SEQ ID NOS: 14 and 15; SEQ ID NOS: 5 and 16; SEQ 35 ID NOS: 5 and 17; and SEQ ID NOS: 5 and 18; for amplification of an endogenous corn gene are included as internal standards for the reaction conditions. Also included are primer pairs that will produce an amplicon in transgenic events containing a pat gene (SEQ ID NOS: 6 and 7; SEQ 40 ID NOS: 8 and 9), and a primer pair that will produce an amplicon in transgenic events containing a cry1F gene (SEQ ID NOS: 19 and 20).

The analysis of plant tissue DNA extracts to test for the presence of the TC1507 event should include a positive 45 tissue DNA extract control (a DNA sample known to contain the transgenic sequences). A successful amplification of the positive control demonstrates that the PCR was run under conditions which allow for the amplification of target sequences. A negative, or wild-type, DNA extract control in 50 which the template DNA provided is either genomic DNA prepared from a non-transgenic plant, or is a non-TC1507 transgenic plant, should also be included. Additionally a negative control that contains no template corn DNA extract will be a useful gauge of the reagents and conditions used in 55 the PCR protocol.

Additional DNA primer molecules of sufficient length can be selected from SEQ ID NO: 26 and SEQ ID NO: 27 by those skilled in the art of DNA amplification methods, and conditions optimized for the production of an amplicon that 60 may differ from the methods shown in Table 9 or Table 10 but result in an amplicon diagnostic for event TC1507. The use of these DNA primer sequences with modifications to the methods shown in Table 9 and Table 10 are within the scope of the invention. The amplicon wherein at least one 65 DNA primer molecule of sufficient length derived from SEQ ID NO: 26 and SEQ ID NO: 27 that is diagnostic for event

TC1507 is an aspect of the invention. The amplicon wherein at least one DNA primer of sufficient length derived from any of the genetic elements of PHI8999A that is diagnostic for event TC1507 is an aspect of the invention. The assay for the TC1507 amplicon can be performed by using a Stratagene Robocycler, MJ Engine, Perkin-Elmer 9700, or Eppendorf Mastercycler Gradient thermocycler, or by methods and apparatus known to those skilled in the art.

TABLE 9

| PCR Conditions:<br>Conditions: |                                |  |  |  |  |  |
|--------------------------------|--------------------------------|--|--|--|--|--|
| Kit used:                      | Perkin-Elmer AmpliTAQ Gold kit |  |  |  |  |  |
| Volume                         | Component                      |  |  |  |  |  |
| 5 μl                           | template (10 ng/μl)            |  |  |  |  |  |
| 4 μl                           | 2 μl each primer (10 μM)       |  |  |  |  |  |
| 2 μl                           | 10X PCR Gold Buffer            |  |  |  |  |  |
| 2 μl                           | 25 mM MgCl <sub>2</sub>        |  |  |  |  |  |
| 2 μl                           | 50X dNTP's (10 mM)             |  |  |  |  |  |
| 0.1 μl                         | Amplitaq Gold Polymerase       |  |  |  |  |  |
| 4.9 µl                         | _H <sub>2</sub> O              |  |  |  |  |  |
| 20 μl                          | Total                          |  |  |  |  |  |
|                                | Cycling Parameters             |  |  |  |  |  |
| Gen                            | eAmp ® PCR System 9700         |  |  |  |  |  |
|                                | 9 min 92° C.                   |  |  |  |  |  |
|                                | 30 cycles:                     |  |  |  |  |  |
|                                | 94° C. 30 sec                  |  |  |  |  |  |
|                                | 60° C. 30 sec                  |  |  |  |  |  |
|                                | 72° C. 1 min                   |  |  |  |  |  |
|                                | 7 min 72° C.                   |  |  |  |  |  |
|                                | Hold 4° C.                     |  |  |  |  |  |

## TABLE 10

PCR Conditions used with the Advantage ®-GC 2 Polymerase Mix: Conditions:

|    | Kit used:<br>Volume | Advantage ®-GC 2 Polymerase Mix<br>Component      |
|----|---------------------|---|
| _  | 5 μl<br>5 μl        | template (10 ng/μl)<br>2.5 μl each primer (10 μM) |
| J  | 10 µl               | 5X GC2 Buffer                                     |
|    | 10 μl               | GC melt (1.0M final conc.)                        |
|    | 1.5 µl              | 50X dNTP's (10 mM)                                |
|    | $1.0~\mu l$         | Advantage GC2 Polymerase                          |
| _  | 17.5 μl             | _H <sub>2</sub> O                                 |
| 5_ | 50 μl               | Total   |

## TABLE 10-continued

| Cycling Parameters GeneAmp ® PCR System 9700                               |  |
|--|--|
| 5 min 94° C.<br>35 cycles:   |  |
| 94° C. 1 min<br>60° C. 2 min<br>72° C. 3 min<br>7 min 72° C.<br>Hold 4° C. |  |

## **40**

Having illustrated and described the principles of the present invention, it should be apparent to persons skilled in the art that the invention can be modified in arrangement and detail without departing from such principles. We claim all modifications that are within the spirit and scope of the appended claims.

All publications and published patent documents cited in this specification are incorporated herein by reference to the same extent as if each individual publication or patent application was specifically and individually indicated to be incorporated by reference.

#### SEQUENCE LISTING

```
<160> NUMBER OF SEQ ID NOS: 57
<210> SEQ ID NO 1
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Event specific primer sequence designed for
      TC1507
<400> SEQUENCE: 1
gtagtactat agattatatt attcgtagag
<210> SEQ ID NO 2
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Event specific primer sequence designed for
      TC1507
<400> SEQUENCE: 2
gccatacaga actcaaaatc ttttccggag
<210> SEQ ID NO 3
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Event specific primer sequence designed for
      TC1507
<400> SEQUENCE: 3
tgtggtgttt gtggctctgt cctaa
<210> SEQ ID NO 4
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Event specific primer for TC1507
<400> SEQUENCE: 4
agcacctttt cattctttca tatac
<210> SEQ ID NO 5
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Genomic DNA primer sequence
<400> SEQUENCE: 5
```

```
25
gacctcccca caggcatgat tgatc
<210> SEQ ID NO 6
<211> LENGTH: 25
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer in full length insert, 35S promoter to
      pat gene
<400> SEQUENCE: 6
tagaggacct aacagaactc gccgt
<210> SEQ ID NO 7
<211> LENGTH: 24
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer in full length insert, 35S promoter to
      pat gene
<400> SEQUENCE: 7
                                                                      24
gagctggcaa ctcaaaatcc cttt
<210> SEQ ID NO 8
<211> LENGTH: 25
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer in full length insert, 35S promoter to
      pat gene
<400> SEQUENCE: 8
                                                                      25
aaaatcttcg tcaacatggt ggagc
<210> SEQ ID NO 9
<211> LENGTH: 25
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer in full length insert, 35S promoter to
      pat gene
<400> SEQUENCE: 9
                                                                      25
taatctcaac tggtctcctc tccgg
<210> SEQ ID NO 10
<211> LENGTH: 22
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer - Zea mays genomic DNA
<400> SEQUENCE: 10
                                                                      22
cccctacccc accgacgttt at
<210> SEQ ID NO 11
<211> LENGTH: 22
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer - Zea mays genomic DNA
<400> SEQUENCE: 11
ttgattggca ggtccgtggg tc
                                                                      22
```

```
<210> SEQ ID NO 12
<211> LENGTH: 22
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer - Zea mays genomic DNA
<400> SEQUENCE: 12
                                                                       22
cacaacggca cagaaacacg aa
<210> SEQ ID NO 13
<211> LENGTH: 22
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer - Zea mays genomic DNA
<400> SEQUENCE: 13
                                                                       22
gcgcacccac cggaacaaaa tg
<210> SEQ ID NO 14
<211> LENGTH: 22
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer - Zea mays genomic DNA
<400> SEQUENCE: 14
tcctcgcatt aaatgctcct gc
                                                                       22
<210> SEQ ID NO 15
<211> LENGTH: 22
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer - Zea mays genomic DNA
<400> SEQUENCE: 15
cctggcacgc attgacgcat gt
<210> SEQ ID NO 16
<211> LENGTH: 25
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer - Zea mays genomic DNA
<400> SEQUENCE: 16
                                                                       25
aagccgccac cgtcttgagg acctt
<210> SEQ ID NO 17
<211> LENGTH: 25
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer - Zea mays genomic DNA
<400> SEQUENCE: 17
gtcgaaccgt gcctccactt gctaa
                                                                       25
<210> SEQ ID NO 18
<211> LENGTH: 25
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer - Zea mays genomic DNA
```

| -continuea   |      |
|--|------|
| <400> SEQUENCE: 18   |      |
| agaaagggaa aatcagggcc aggac  | 25   |
| <210> SEQ ID NO 19 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE:   |      |
| <223> OTHER INFORMATION: Cry1F sequence primer   |      |
| <400> SEQUENCE: 19 ggctcggact cgacctttct at  | 22   |
|  |      |
| <210> SEQ ID NO 20 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Cry1F sequence primer                  |      |
| <400> SEQUENCE: 20   |      |
| gcagttcttg aagaatgagt ga   | 22   |
| <210> SEQ ID NO 21 <211> LENGTH: 2829 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: 5' flanking sequence of event TC1507 |      |
| <400> SEQUENCE: 21   |      |
| actagtttcc tagcccgcgt cgtgccccta ccccaccgac gtttatggaa ggtgccattc  | 60   |
| cacggttctt cgtggccgcc cctaaggatg taaatggtcg gtaaaatccg gtaaatttcc  | 120  |
| ggtaccgttt accagatttt tccagccgtt ttcggattta tcgggatata cagaaaacga  | 180  |
| gacggaaacg gaataggttt tttttcgaaa acggtacggt  | 240  |
| cgtccgtttt cgtatttctc gggaaactct ggtatattcc cgtatttgtc ccgtattttc  | 300  |
| ccgacccacg gacctgccaa tcaaccatca gccagtcagc ccatccccac agctatggcc  | 360  |
| catggggcca tgttggccac atgcccacgc aacgcaaggc agtaaggctg gcagcctggc  | 420  |
| acgcattgac gcatgtggac acacacagcc gccgcctgtt cgtgtttctg tgccgttgtg  | 480  |
| cgagactgtg actgcgagtg gcggagtcgg cgaacggcga ggcgtctccg gagtctggac  | 540  |
| tgcggctgtg gacagcgacg ctgtgacggc gactcggcga agccccaagc taccaagccc  | 600  |
| ccaagtcccc atccatctct gcttctctgg tcatctcctt cccctggtcg atctgcaggc  | 660  |
| gccagaccgg ccgaagcatc acgaaacgca ctaagacctc gaaggagtca aaccactcct  | 720  |
| ccgaggcctc gggggctaca cccggcgggt gcgctcgcgc gcacccaccg gaacaaaatg  | 780  |
| taaccgagaa aggtcggtcc ccttgcaaaa aaagtgcgac aaaagcctcc aagcgagtat  | 840  |
| taacactcac tttgaggctc gggggctact gtcggggacc ataattaggg gtacccccaa  | 900  |
| gactcctaat ctcagctggt aacccccatc agcacaaagc tgcaaaggcc tgatgggtgc  | 960  |
| gattaagtca aggctcggtc cactcaaggg acacgatctc gcctcgcccg agcccagcct  | 1020 |
| cgggcaaggg cggccgaccc cgaggattca cgtctcgccc gagggccccc tcaagcgacg  | 1080 |
| ggcacacctt cggctcgccc gaggcccatt cttcgccgag aagcaacctt ggccagatcg  | 1140 |
| ccacaccgac cgaccgtatc gcaggagcat ttaatgcgag gatcgcctga caccttatcc  | 1200 |
|  |      |

2340

2580

-continued tgacgcgcgc tcttcagtcg acagagccga agtgaccgca atcacttcgc cgctccactg 1260 1320 accgacctga caagaagaca gcgccgcctg cgtcgctccg actgctgtgc cactcgacag 1380 agtgaggetg acageageea agteeggeet egggegeeat aggaagetee geetegeeeg 1440 accetaggge teggaetegg ceteggetee ggaagaegae gaaetaeget tegeeegaee 1500 ccagggcttg gactcagcct cggctccgga agacgacgaa ttccgcctcg cccgacccca 1560 gggctcggac tcggcctcgg ctccagaaga cgacgaactc cgcctcgccc gaccccaggg 1620 ctcggactca gcctcggctc cggaagacga cgaactccgc ctcgcccgac cccagggctc ggactcagcc tcggcctcag acgatggtct ccgcctcgcc cgacccgggg ctcggactcg acctttctat cggaccttgt cagatcctgt cttcgtccga ggaggctttg gcaatcctca 1740 1800 ctatgtactc ggtcttaggg gagtggcctt tcaacaaact ggtacgaatc acacccgcac 1860 attcaggaac tccgggacca ttgactctct agatgagata ccacctcaag acaacagcgg

1920 cgcaccttgg aatgactact cccatgtgct gaatcatgtt acctttgtgc gctggccagg tgagatetea ggtteegaet eatggagage accaatgtte tettggaege ategtagege 1980 2040 tacccccaca aacaccattg atccagagag aatcactcat tcttcaagaa ctgcatatct 2100 tgccgagatc ctcatcccta aaggtacttg acaatagtat tattggagtc gatacacaac 2160 tcacaaaaaa tacaagaagt cgactaggtg gattggtccg agtgaagaga aaaaaaagcc

2220 atacagaact caaaatcttt tccggagata ttcattttcc tgaagaggcg gataagatat 2280 taggtggcag tttgatacca ccagaaagag aaaaaaaaga ttctaaggaa tcaaaaaaaa

ggaaaaattg ggtttatgtt caacggaaaa aatttctcaa aagcaaggaa aagtattgtg

ctggtattgc ttacgctggg ccctggaagg ctaggaaccc tcaacctcag caaccaacca

2400 gctatttatc tatccgtgca gctgatatgg ccgcggtttg tgatatcgtt aaccattaca 2460 ttgagacgtc tacagtgaac tttaggacag agccacaaac accacaagag tggattgatg

2520 atctagagag gttgcaagat agataccctt ggttggttgc tgaggttgag ggtgttgtgg

2640 atggtatcta tcttgcaacc tctctagatc atcaatccac tcttgtggtg tttgtggctc 2700 tgtcctaaag ttcactgtag acgtctcaat gtaatggtta acgatatcac aaaccgagag

2760 aagagggate tegaagette ggeegggee categatate egegggeatg eetgeagtge 2820 agegtgaeee ggtegtgeee etetetagag ataatgagea ttgeatgtet aagttataaa

2829 aaattacca

<400> SEQUENCE: 22

60 ctcactccgc ttgatcttgg caaagatatt tgacgcattt attagtatgt gttaattttc atttgcagtg cagtattttc tattcgatct ttatgtaatt cgttacaatt aataaatatt 180 caaatcagat tattgactgt catttgtatc aaatcgtgtt taatggatat ttttattata atattgatga tatctcaatc aaaacgtaga taataataat atttatttaa tattttgcg 240 togoacagtg aaaatotata tgagattaca aaatacogac aacattattt aagaaacata 300 360 gacattaacc ctgagactgt tggacatcaa cgggtagatt ccttcatgca tagcacctca 420 ttcttgggga caaaagcacg gtttggccgt tccattgctg cacgaacgag ctttgctata

<sup>&</sup>lt;210> SEQ ID NO 22

<sup>&</sup>lt;211> LENGTH: 2346

<sup>&</sup>lt;212> TYPE: DNA

<sup>&</sup>lt;213 > ORGANISM: Artificial Sequence

<sup>&</sup>lt;220> FEATURE:

<sup>&</sup>lt;223> OTHER INFORMATION: 3' flanking sequence of event TC1507

-continued

```
teetegggtt ggateatete ateaggteea ateaaatttg teeaagaaet eatgttagte
                                                                     480
                                                                     540
gcaacgaaac cggggcatat gtcgggtatc tcgagctcgc gaaagcttgg ctgcaggtcg
                                                                     600
acggatcctt caacaaaagg gtacctgtac ccgaaaccga cacaggtggg taggtagaga
                                                                     660
atacctaggg gcgcgagaca actctctcta aggaactcgg caaaatagcc ccgtaacttc
                                                                     720
gggagaaggg gtgccccccg ctaacaataa acgaatacgg tttatgtatg gattccggta
                                                                     780
aaataccggt actcgatttc ataagagtcg aataggaagt taagatgagg gtggtatcat
                                                                     840
cataaaaatg gagtagtatc ctaaattata ctaatccacg tatgatatgt atgcctttcc
                                                                     900
ttatcaaccg gaagtagtgc aaaaaaaatt ctatactgca ctgctctctt tttactgaga
                                                                     960
aatgcaaaaa aataaaagtg aagtaagggt gccccataga tatttgatct tgcctcctgt
                                                                    1020
ccccccccc cttttttcat caaaaatttc catgaaaaaa gaaaagatga atttgtccat
                                                                    1080
tcattgaacc ctagttcggg actgacgggg ctcgaacccg cagcttccgc ctgttcctag
                                                                    1140
ccttccaggg cccagcgtaa gcaataccag ccacagcacc ctcaacctca gcaaccaacc
                                                                    1200
aagggtatct atcttgcaac ctctctagat catcaatcca ctcttgtggt gtttgtggct
                                                                    1260
ctgtcctaaa gttcactgta gacgtctcaa tgtaatggtt aacgatatca caaaccgcgg
                                                                    1320
aacacaagaa cgaaagcacc ttttcattct ttcatatact aggggttttt acttggaaaa
                                                                    1380
gacaatgttc catactaaag gatagctgca gaagccgcca ccgtcttgag gaccttccgg
                                                                    1440
ggagccagac cggtcgaacc gtgcctccac ttgctaagga gaaagggaaa atcagggcca
                                                                    1500
ggacatacga aggaggagcc agaacgaaga tatcctaaga tacttactcg ctccgggcca
                                                                    1560
tgatcaatca tgcctgtggg gaggtctctc gcacctcgat ccatgaaggt accaccgagg
                                                                    1620
tetgeecege egeeggette ggtacegtee tegeettggg egeeegagge acceggggga
tggactgccc aggcgcagcc acgacgaccc aaggatcacc ctcctgcgca gtcggcacga
                                                                    1740
gcaatagttc tcggggaaca ggcagcttgg cctgactccc cggggtcacc tcaactacct
                                                                    1800
eggeegaggg gteaagtace eceteagtee geeceegete tteggaeegg gaeecegaeg
teceggeece ggatacegae ggeaceagee egeteggggg etggettgae gaeecetgge
                                                                    1860
                                                                    1920
ccagcctcag atctgggctg aggccgaggc aggcggccat gtcgtcgtct tcatcatcgt
cttcatcatc gtcgtcgtca tcaggcgtct ccggcgacgg ctcccttggg agcccctccc
                                                                    1980
                                                                    2040
teteetgeeg acgaegaage ettteeaagg cateeegage eeaegteege tegtgggeee
                                                                    2100
gageettett tgegteette tteteettee tetteteege ggtgaeeete egegeagete
                                                                    2160
ggtccaccgc atcctccggg actggtggca gggaaggctt gtgatgccct acctcctgga
                                                                    2220
gacagacgaa aagtctcagc tatgagaacc gagggcaatc tgacgcaaga aggaagaagg
                                                                    2280
ageggataet eaceagagae aegeaeeege gategggaeg eattaaggge tgggaaaaag
                                                                    2340
tgccggcctc taatttcgct accgtgccgt ccacccacct gtggaggtca tcgatgggaa
                                                                    2346
ggggaa
```

<400> SEQUENCE: 23

<sup>&</sup>lt;210> SEQ ID NO 23

<sup>&</sup>lt;211> LENGTH: 20

<sup>&</sup>lt;212> TYPE: DNA

<sup>&</sup>lt;213> ORGANISM: Artificial Sequence

<sup>&</sup>lt;220> FEATURE:

<sup>&</sup>lt;223> OTHER INFORMATION: Event specific primer sequence designed for TC1507

-continued

**51** 

<210> SEQ ID NO 24

<211> LENGTH: 11361

<212> TYPE: DNA

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: The sequence represents the transgenic insert in maize line TC1507 as well as the sequence flanking the insertion sites.

<400> SEQUENCE: 24

| actagtttcc tagcccgcgt | cgtgccccta | ccccaccgac | gtttatggaa | ggtgccattc | 60   |
|-----------------------|------------|------------|------------|------------|------|
| cacggttctt cgtggccgcc | cctaaggatg | taaatggtcg | gtaaaatccg | gtaaatttcc | 120  |
| ggtaccgttt accagatttt | tccagccgtt | ttcggattta | tcgggatata | cagaaaacga | 180  |
| gacggaaacg gaataggttt | tttttcgaaa | acggtacggt | aaacggtgag | acaaacttac | 240  |
| cgtccgtttt cgtatttctc | gggaaactct | ggtatattcc | cgtatttgtc | ccgtattttc | 300  |
| ccgacccacg gacctgccaa | tcaaccatca | gccagtcagc | ccatccccac | agctatggcc | 360  |
| catggggcca tgttggccac | atgcccacgc | aacgcaaggc | agtaaggctg | gcagcctggc | 420  |
| acgcattgac gcatgtggac | acacacagcc | gccgcctgtt | cgtgtttctg | tgccgttgtg | 480  |
| cgagactgtg actgcgagtg | gcggagtcgg | cgaacggcga | ggcgtctccg | gagtctggac | 540  |
| tgcggctgtg gacagcgacg | ctgtgacggc | gactcggcga | agccccaagc | taccaagccc | 600  |
| ccaagtcccc atccatctct | gcttctctgg | tcatctcctt | cccctggtcg | atctgcaggc | 660  |
| gccagaccgg ccgaagcatc | acgaaacgca | ctaagacctc | gaaggagtca | aaccactcct | 720  |
| ccgaggcctc gggggctaca | cccggcgggt | gcgctcgcgc | gcacccaccg | gaacaaaatg | 780  |
| taaccgagaa aggtcggtcc | ccttgcaaaa | aaagtgcgac | aaaagcctcc | aagcgagtat | 840  |
| taacactcac tttgaggctc | gggggctact | gtcggggacc | ataattaggg | gtacccccaa | 900  |
| gactcctaat ctcagctggt | aacccccatc | agcacaaagc | tgcaaaggcc | tgatgggtgc | 960  |
| gattaagtca aggctcggtc | cactcaaggg | acacgatctc | gcctcgcccg | agcccagcct | 1020 |
| cgggcaaggg cggccgaccc | cgaggattca | cgtctcgccc | gagggccccc | tcaagcgacg | 1080 |
| ggcacacctt cggctcgccc | gaggcccatt | cttcgccgag | aagcaacctt | ggccagatcg | 1140 |
| ccacaccgac cgaccgtatc | gcaggagcat | ttaatgcgag | gatcgcctga | caccttatcc | 1200 |
| tgacgcgcgc tcttcagtcg | acagagccga | agtgaccgca | atcacttcgc | cgctccactg | 1260 |
| accgacctga caagaagaca | gcgccgcctg | cgtcgctccg | actgctgtgc | cactcgacag | 1320 |
| agtgaggctg acagcagcca | agtccggcct | cgggcgccat | aggaagctcc | gcctcgcccg | 1380 |
| accctagggc tcggactcgg | cctcggctcc | ggaagacgac | gaactacgct | tcgcccgacc | 1440 |
| ccagggcttg gactcagcct | cggctccgga | agacgacgaa | ttccgcctcg | cccgacccca | 1500 |
| gggctcggac tcggcctcgg | ctccagaaga | cgacgaactc | cgcctcgccc | gaccccaggg | 1560 |
| ctcggactca gcctcggctc | cggaagacga | cgaactccgc | ctcgcccgac | cccagggctc | 1620 |
| ggactcagcc tcggcctcag | acgatggtct | ccgcctcgcc | cgacccgggg | ctcggactcg | 1680 |
| acctttctat cggaccttgt | cagatcctgt | cttcgtccga | ggaggctttg | gcaatcctca | 1740 |
| ctatgtactc ggtcttaggg | gagtggcctt | tcaacaaact | ggtacgaatc | acacccgcac | 1800 |
| attcaggaac tccgggacca | ttgactctct | agatgagata | ccacctcaag | acaacagcgg | 1860 |
| cgcaccttgg aatgactact | cccatgtgct | gaatcatgtt | acctttgtgc | gctggccagg | 1920 |
| tgagatctca ggttccgact | catggagagc | accaatgttc | tcttggacgc | atcgtagcgc | 1980 |

| tacccccaca | aacaccattg | atccagagag | aatcactcat | tcttcaagaa | ctgcatatct | 2040 |
|------------|------------|------------|------------|------------|------------|------|
| tgccgagatc | ctcatcccta | aaggtacttg | acaatagtat | tattggagtc | gatacacaac | 2100 |
| tcacaaaaaa | tacaagaagt | cgactaggtg | gattggtccg | agtgaagaga | aaaaaagcc  | 2160 |
| atacagaact | caaaatcttt | tccggagata | ttcattttcc | tgaagaggcg | gataagatat | 2220 |
| taggtggcag | tttgatacca | ccagaaagag | aaaaaaaga  | ttctaaggaa | tcaaaaaaa  | 2280 |
| ggaaaaattg | ggtttatgtt | caacggaaaa | aatttctcaa | aagcaaggaa | aagtattgtg | 2340 |
| gctatttatc | tatccgtgca | gctgatatgg | ccgcggtttg | tgatatcgtt | aaccattaca | 2400 |
| ttgagacgtc | tacagtgaac | tttaggacag | agccacaaac | accacaagag | tggattgatg | 2460 |
| atctagagag | gttgcaagat | agataccctt | ggttggttgc | tgaggttgag | ggtgttgtgg | 2520 |
| ctggtattgc | ttacgctggg | ccctggaagg | ctaggaaccc | tcaacctcag | caaccaacca | 2580 |
| atggtatcta | tcttgcaacc | tctctagatc | atcaatccac | tcttgtggtg | tttgtggctc | 2640 |
| tgtcctaaag | ttcactgtag | acgtctcaat | gtaatggtta | acgatatcac | aaaccgagag | 2700 |
| aagagggatc | tcgaagcttc | ggccggggcc | catcgatatc | cgcgggcatg | cctgcagtgc | 2760 |
| agcgtgaccc | ggtcgtgccc | ctctctagag | ataatgagca | ttgcatgtct | aagttataaa | 2820 |
| aaattaccac | aactggaaga | gcggttaccc | ggaccgaagc | ttcggccggg | gcccatcgat | 2880 |
| atccgcgggc | atgcctgcag | tgcagcgtga | cccggtcgtg | cccctctcta | gagataatga | 2940 |
| gcattgcatg | tctaagttat | aaaaattac  | cacatattt  | ttttgtcaca | cttgtttgaa | 3000 |
| gtgcagttta | tctatcttta | tacatatatt | taaactttac | tctacgaata | atataatcta | 3060 |
| tagtactaca | ataatatcag | tgttttagag | aatcatataa | atgaacagtt | agacatggtc | 3120 |
| taaaggacaa | ttgagtattt | tgacaacagg | actctacagt | tttatcttt  | tagtgtgcat | 3180 |
| gtgttctcct | tttttttgc  | aaatagcttc | acctatataa | tacttcatcc | attttattag | 3240 |
| tacatccatt | tagggtttag | ggttaatggt | ttttatagac | taatttttt  | agtacatcta | 3300 |
| ttttattcta | ttttagcctc | taaattaaga | aaactaaaac | tctattttag | ttttttatt  | 3360 |
| taataattta | gatataaaat | agaataaaat | aaagtgacta | aaaattaaac | aaataccctt | 3420 |
| taagaaatta | aaaaactaa  | ggaaacattt | ttcttgtttc | gagtagataa | tgccagcctg | 3480 |
| ttaaacgccg | tcgacgagtc | taacggacac | caaccagcga | accagcagcg | tcgcgtcggg | 3540 |
| ccaagcgaag | cagacggcac | ggcatctctg | tcgctgcctc | tggacccctc | tcgagagttc | 3600 |
| cgctccaccg | ttggacttgc | tccgctgtcg | gcatccagaa | attgcgtggc | ggagcggcag | 3660 |
| acgtgagccg | gcacggcagg | cggcctcctc | ctcctctcac | ggcacggcag | ctacggggga | 3720 |
| ttcctttccc | accgctcctt | cgctttccct | tcctcgcccg | ccgtaataaa | tagacacccc | 3780 |
| ctccacaccc | tctttcccca | acctcgtgtt | gttcggagcg | cacacacaca | caaccagatc | 3840 |
| tcccccaaat | ccacccgtcg | gcacctccgc | ttcaaggtac | gccgctcgtc | ctcccccc   | 3900 |
| cccctctct  | accttctcta | gatcggcgtt | ccggtccatg | gttagggccc | ggtagttcta | 3960 |
| cttctgttca | tgtttgtgtt | agatccgtgt | ttgtgttaga | tccgtgctgc | tagcgttcgt | 4020 |
| acacggatgc | gacctgtacg | tcagacacgt | tctgattgct | aacttgccag | tgtttctctt | 4080 |
| tggggaatcc | tgggatggct | ctagccgttc | cgcagacggg | atcgatttca | tgatttttt  | 4140 |
| tgtttcgttg | catagggttt | ggtttgccct | tttcctttat | ttcaatatat | gccgtgcact | 4200 |
| tgtttgtcgg | gtcatcttt  | catgctttt  | tttgtcttgg | ttgtgatgat | gtggtctggt | 4260 |
|            |            |            |            | ctacctggtg |            | 4320 |
|            |            |            |            | cgaattgaag |            | 4380 |
|            | 55-5-5     | 5          |            |            |            |      |

| gaaatatcga     | tctaggatag | gtatacatgt     | tgatgcgggt | tttactgatg | catatacaga | 4440 |
|----------------|------------|----------------|------------|------------|------------|------|
| gatgcttttt     | gttcgcttgg | ttgtgatgat     | gtggtgtggt | tgggcggtcg | ttcattcgtt | 4500 |
| ctagatcgga     | gtagaatact | gtttcaaact     | acctggtgta | tttattaatt | ttggaactgt | 4560 |
| atgtgtgtgt     | catacatctt | catagttacg     | agtttaagat | ggatggaaat | atcgatctag | 4620 |
| gataggtata     | catgttgatg | tgggttttac     | tgatgcatat | acatgatggc | atatgcagca | 4680 |
| tctattcata     | tgctctaacc | ttgagtacct     | atctattata | ataaacaagt | atgttttata | 4740 |
| attattttga     | tcttgatata | cttggatgat     | ggcatatgca | gcagctatat | gtggatttt  | 4800 |
| ttagccctgc     | cttcatacgc | tatttatttg     | cttggtactg | tttcttttgt | cgatgctcac | 4860 |
| cctgttgttt     | ggtgttactt | ctgcaggtcg     | actctagagg | atccaacaat | ggagaacaac | 4920 |
| atacagaatc     | agtgcgtccc | ctacaactgc     | ctcaacaatc | ctgaagtaga | gattctcaac | 4980 |
| gaagagaggt     | cgactggcag | attgccgtta     | gacatctccc | tgtcccttac | acgtttcctg | 5040 |
| ttgtctgagt     | ttgttccagg | tgtgggagtt     | gcgtttggcc | tcttcgacct | catctggggc | 5100 |
| ttcatcactc     | catctgattg | gagcctcttt     | cttctccaga | ttgaacagtt | gattgaacaa | 5160 |
| aggattgaga     | ccttggaaag | gaatcgggcc     | atcactaccc | ttcgtggctt | agcagacagc | 5220 |
| tatgagatct     | acattgaagc | actaagagag     | tgggaagcca | atcctaacaa | tgcccaactg | 5280 |
| agagaagatg     | tgcgtatacg | ctttgctaac     | acagatgatg | ctttgatcac | agccatcaac | 5340 |
| aacttcaccc     | ttaccagctt | cgagatccct     | cttctctcgg | tctatgttca | agctgctaac | 5400 |
| ctgcacttgt     | cactactgcg | cgacgctgtg     | tcgtttgggc | aaggttgggg | actggacata | 5460 |
| gctactgtca     | acaatcacta | caacagactc     | atcaatctga | ttcatcgata | cacgaaacat | 5520 |
| tgtttggata     | cctacaatca | gggattggag     | aacctgagag | gtactaacac | tcgccaatgg | 5580 |
| gccaggttca     | atcagttcag | gagagacctt     | acacttactg | tgttagacat | agttgctctc | 5640 |
| tttccgaact     | acgatgttcg | tacctatccg     | attcaaacgt | catcccaact | tacaagggag | 5700 |
| atctacacca     | gttcagtcat | tgaagactct     | ccagtttctg | cgaacatacc | caatggtttc | 5760 |
| aacagggctg     | agtttggagt | cagaccaccc     | catctcatgg | acttcatgaa | ctctttgttt | 5820 |
| gtgactgcag     | agactgttag | atcccaaact     | gtgtggggag | gacacttagt | tagctcacgc | 5880 |
| aacacggctg     | gcaatcgtat | caactttcct     | agttacgggg | tcttcaatcc | cgggggcgcc | 5940 |
| atctggattg     | cagatgaaga | tccacgtcct     | ttctatcgga | ccttgtcaga | tcctgtcttc | 6000 |
| gtccgaggag     | gctttggcaa | tcctcactat     | gtactcggtc | ttaggggagt | ggcctttcaa | 6060 |
| caaactggta     | cgaatcacac | ccgcacattc     | aggaactccg | ggaccattga | ctctctagat | 6120 |
| gagataccac     | ctcaagacaa | cagcggcgca     | ccttggaatg | actactccca | tgtgctgaat | 6180 |
| catgttacct     | ttgtgcgctg | gccaggtgag     | atctcaggtt | ccgactcatg | gagagcacca | 6240 |
| atgttctctt     | ggacgcatcg | tagcgctacc     | cccacaaaca | ccattgatcc | agagagaatc | 6300 |
| actcagattc     | ccttggtgaa | ggcacacaca     | cttcagtcag | gaactacagt | tgtaagaggg | 6360 |
| ccggggttca     | cgggaggaga | cattcttcga     | cgcactagtg | gaggaccatt | cgcgtacacc | 6420 |
| attgtcaaca     | tcaatgggca | acttccccaa     | aggtatcgtg | ccaggatacg | ctatgcctct | 6480 |
| actaccaatc     | taagaatcta | cgttacggtt     | gcaggtgaac | ggatctttgc | tggtcagttc | 6540 |
| aacaagacaa     | tggataccgg | tgatccactt     | acattccaat | ctttctccta | cgccactatc | 6600 |
| aacaccgcgt     | tcacctttcc | aatgagccag     | agcagtttca | cagtaggtgc | tgataccttc | 6660 |
|                | acgaagtgta |                |            |            |            | 6720 |
| _ <del>_</del> | _          | - <del>-</del> |            |            |            |      |

-continued

|            |            |            |            | -COILLI    | iuea       |      |
|------------|------------|------------|------------|------------|------------|------|
| gagtaaggat | ccgtcgacct | gcagccaagc | tttcgcgagc | tcgagatccc | cgacatatgc | 6780 |
| cccggtttcg | ttgcgactaa | catgagttct | tggacaaatt | tgattggacc | tgatgagatg | 6840 |
| atccaacccg | aggatatagc | aaagctcgtt | cgtgcagcaa | tggaacggcc | aaaccgtgct | 6900 |
| tttgtcccca | agaatgaggt | gctatgcatg | aaggaatcta | cccgttgatg | tccaacagtc | 6960 |
| tcagggttaa | tgtctatgta | tcttaaataa | tgttgtcggt | attttgtaat | ctcatataga | 7020 |
| ttttcactgt | gcgacgcaaa | aatattaaat | aaatattatt | attatctacg | ttttgattga | 7080 |
| gatatcatca | atattataat | aaaaatatcc | attaaacacg | atttgataca | aatgacagtc | 7140 |
| aataatctga | tttgaatatt | tattaattgt | aacgaattac | ataaagatcg | aatagaaaat | 7200 |
| actgcactgc | aaatgaaaat | taacacatac | taataaatgc | gtcaaatatc | tttgccaaga | 7260 |
| tcaagcggag | tgagggcctc | atatccggtc | tcagttacaa | gcacggtatc | cccgaagcgc | 7320 |
| gctccaccaa | tgccctcgac | atagatgccg | ggctcgacgc | tgaggacatt | gcctaccttg | 7380 |
| agcatggtct | cagcgccggc | tttaagctca | atcccatccc | aatctgaata | tcctatcccg | 7440 |
| cgcccagtcc | ggtgtaagaa | cgggtctgtc | catccacctc | tgttgggaat | tccggtccgg | 7500 |
| gtcacctttg | tccaccaaga | tggaactgcg | gccgcggacc | gaattcccat | ggagtcaaag | 7560 |
| attcaaatag | aggacctaac | agaactcgcc | gtaaagactg | gcgaacagtt | catacagagt | 7620 |
| ctcttacgac | tcaatgacaa | gaagaaaatc | ttcgtcaaca | tggtggagca | cgacacgctt | 7680 |
| gtctactcca | aaaatatcaa | agatacagtc | tcagaagacc | aaagggcaat | tgagactttt | 7740 |
| caacaaaggg | taatatccgg | aaacctcctc | ggattccatt | gcccagctat | ctgtcacttt | 7800 |
| attgtgaaga | tagtggaaaa | ggaaggtggc | tcctacaaat | gccatcattg | cgataaagga | 7860 |
| aaggccatcg | ttgaagatgc | ctctgccgac | agtggtccca | aagatggacc | cccacccacg | 7920 |
| aggagcatcg | tggaaaaaga | agacgttcca | accacgtctt | caaagcaagt | ggattgatgt | 7980 |
| gatatctcca | ctgacgtaag | ggatgacgca | caatcccact | atccttcgca | agacccttcc | 8040 |
| tctatataag | gaagttcatt | tcatttggag | aggacagggt | acccggggat | ccaccatgtc | 8100 |
| tccggagagg | agaccagttg | agattaggcc | agctacagca | gctgatatgg | ccgcggtttg | 8160 |
| tgatatcgtt | aaccattaca | ttgagacgtc | tacagtgaac | tttaggacag | agccacaaac | 8220 |
| accacaagag | tggattgatg | atctagagag | gttgcaagat | agataccctt | ggttggttgc | 8280 |
| tgaggttgag | ggtgttgtgg | ctggtattgc | ttacgctggg | ccctggaagg | ctaggaacgc | 8340 |
| ttacgattgg | acagttgaga | gtactgttta | cgtgtcacat | aggcatcaaa | ggttgggcct | 8400 |
| aggatccaca | ttgtacacac | atttgcttaa | gtctatggag | gcgcaaggtt | ttaagtctgt | 8460 |
| ggttgctgtt | ataggccttc | caaacgatcc | atctgttagg | ttgcatgagg | ctttgggata | 8520 |
| cacagcccgg | ggtacattgc | gcgcagctgg | atacaagcat | ggtggatggc | atgatgttgg | 8580 |
| tttttggcaa | agggattttg | agttgccagc | tcctccaagg | ccagttaggc | cagttaccca | 8640 |
| gatctgagtc | gacctgcagg | catgcccgct | gaaatcacca | gtctctctct | acaaatctat | 8700 |
| ctctctctat | aataatgtgt | gagtagttcc | cagataaggg | aattagggtt | cttatagggt | 8760 |
| ttcgctcatg | tgttgagcat | ataagaaacc | cttagtatgt | atttgtattt | gtaaaatact | 8820 |
| tctatcaata | aaatttctaa | ttcctaaaac | caaaatccag | tggcgagctc | gaattcgagc | 8880 |
| tcgagcccgg | gtggatcctc | tagagtcgac | ctgcagaagc | ttcggtccgg | cgcgcctcta | 8940 |
| gttgaagaca | cgttcatgtc | ttcatcgtaa | gaagacactc | agtagtcttc | ggccagaatg | 9000 |
| gcctaactca | aggccctcac | tccgcttgat | cttggcaaag | atatttgacg | catttattag | 9060 |
| tatgtgttaa | ttttcatttg | cagtgcagta | ttttctattc | gatctttatg | taattcgtta | 9120 |
|            |            |            |            |            |            |      |

**58** 

| caattaataa | atattcaaat | cagattattg | actgtcattt | gtatcaaatc | gtgtttaatg | 9180  |
|------------|------------|------------|------------|------------|------------|-------|
| gatatttta  | ttataatatt | gatgatatct | caatcaaaac | gtagataata | ataatatta  | 9240  |
| tttaatattt | ttgcgtcgca | cagtgaaaat | ctatatgaga | ttacaaaata | ccgacaacat | 9300  |
| tatttaagaa | acatagacat | taaccctgag | actgttggac | atcaacgggt | agattccttc | 9360  |
| atgcatagca | cctcattctt | ggggacaaaa | gcacggtttg | gccgttccat | tgctgcacga | 9420  |
| acgagctttg | ctatatcctc | gggttggatc | atctcatcag | gtccaatcaa | atttgtccaa | 9480  |
| gaactcatgt | tagtcgcaac | gaaaccgggg | catatgtcgg | gtatctcgag | ctcgcgaaag | 9540  |
| cttggctgca | ggtcgacgga | tccttcaaca | aaagggtacc | tgtacccgaa | accgacacag | 9600  |
| gtgggtaggt | agagaatacc | taggggcgcg | agacaactct | ctctaaggaa | ctcggcaaaa | 9660  |
| tagccccgta | acttcgggag | aaggggtgcc | ccccgctaac | aataaacgaa | tacggtttat | 9720  |
| gtatggattc | cggtaaaata | ccggtactcg | atttcataag | agtcgaatag | gaagttaaga | 9780  |
| tgagggtggt | atcatcataa | aaatggagta | gtatcctaaa | ttatactaat | ccacgtatga | 9840  |
| tatgtatgcc | tttccttatc | aaccggaagt | agtgcaaaaa | aaattctata | ctgcactgct | 9900  |
| ctctttttac | tgagaaatgc | aaaaaataa  | aagtgaagta | agggtgcccc | atagatattt | 9960  |
| gatcttgcct | cctgtccccc | ccccctttt  | ttcatcaaaa | atttccatga | aaaaagaaaa | 10020 |
| gatgaatttg | tccattcatt | gaaccctagt | tcgggactga | cggggctcga | acccgcagct | 10080 |
| tccgcctgtt | cctagccttc | cagggcccag | cgtaagcaat | accagccaca | gcaccctcaa | 10140 |
| cctcagcaac | caaccaaggg | tatctatctt | gcaacctctc | tagatcatca | atccactctt | 10200 |
| gtggtgtttg | tggctctgtc | ctaaagttca | ctgtagacgt | ctcaatgtaa | tggttaacga | 10260 |
| tatcacaaac | cgcggaacac | aagaacgaaa | gcaccttttc | attctttcat | atactagggg | 10320 |
| tttttacttg | gaaaagacaa | tgttccatac | taaaggatag | ctgcagaagc | cgccaccgtc | 10380 |
| ttgaggacct | tccggggagc | cagaccggtc | gaaccgtgcc | tccacttgct | aaggagaaag | 10440 |
| ggaaaatcag | ggccaggaca | tacgaaggag | gagccagaac | gaagatatcc | taagatactt | 10500 |
| actcgctccg | ggccatgatc | aatcatgcct | gtggggaggt | ctctcgcacc | tcgatccatg | 10560 |
| aaggtaccac | cgaggtctgc | cccgccgccg | gcttcggtac | cgtcctcgcc | ttgggcgccc | 10620 |
| gaggcacccg | ggggatggac | tgcccaggcg | cagccacgac | gacccaagga | tcaccctcct | 10680 |
| gcgcagtcgg | cacgagcaat | agttctcggg | gaacaggcag | cttggcctga | ctccccgggg | 10740 |
| tcacctcaac | tacctcggcc | gaggggtcaa | gtaccccctc | agtccgcccc | cgctcttcgg | 10800 |
| accgggaccc | cgacgtcccg | gccccggata | ccgacggcac | cagcccgctc | gggggctggc | 10860 |
| ttgacgaccc | ctggcccagc | ctcagatctg | ggctgaggcc | gaggcaggcg | gccatgtcgt | 10920 |
| cgtcttcatc | atcgtcttca | tcatcgtcgt | cgtcatcagg | cgtctccggc | gacggctccc | 10980 |
| ttgggagccc | ctccctctcc | tgccgacgac | gaagcctttc | caaggcatcc | cgagcccacg | 11040 |
| tccgctcgtg | ggcccgagcc | ttctttgcgt | ccttcttctc | cttcctctc  | tccgcggtga | 11100 |
| ccctccgcgc | agctcggtcc | accgcatcct | ccgggactgg | tggcagggaa | ggcttgtgat | 11160 |
| gccctacctc | ctggagacag | acgaaaagtc | tcagctatga | gaaccgaggg | caatctgacg | 11220 |
| caagaaggaa | gaaggagcgg | atactcacca | gagacacgca | cccgcgatcg | ggacgcatta | 11280 |
| agggctggga | aaaagtgccg | gcctctaatt | tcgctaccgt | gccgtccacc | cacctgtgga | 11340 |
| ggtcatcgat | gggaagggga | a          |            |            |            | 11361 |
|            |            |            |            |            |            |       |

-continued

<211> LENGTH: 6186

<212> TYPE: DNA

<213 > ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: The sequence represents the DNA molecule used to transform maize line TC1507. This sequence represents insert PHI8999A.

<400> SEQUENCE: 25

| caactg | gaag  | agcggttacc | cggaccgaag | cttcggccgg | ggcccatcga | tatccgcggg | 60   |
|--------|-------|------------|------------|------------|------------|------------|------|
| catgcc | tgca  | gtgcagcgtg | acccggtcgt | gcccctctct | agagataatg | agcattgcat | 120  |
| gtctaa | gtta  | taaaaaatta | ccacatattt | tttttgtcac | acttgtttga | agtgcagttt | 180  |
| atctat | cttt  | atacatatat | ttaaacttta | ctctacgaat | aatataatct | atagtactac | 240  |
| aataat | atca  | gtgttttaga | gaatcatata | aatgaacagt | tagacatggt | ctaaaggaca | 300  |
| attgag | tatt  | ttgacaacag | gactctacag | ttttatcttt | ttagtgtgca | tgtgttctcc | 360  |
| ttttt  | tttg  | caaatagctt | cacctatata | atacttcatc | cattttatta | gtacatccat | 420  |
| ttaggg | ttta  | gggttaatgg | tttttataga | ctaattttt  | tagtacatct | attttattct | 480  |
| atttta | .gcct | ctaaattaag | aaaactaaaa | ctctatttta | gtttttttat | ttaataattt | 540  |
| agatat | aaaa  | tagaataaaa | taaagtgact | aaaaattaaa | caaataccct | ttaagaaatt | 600  |
| aaaaaa | acta. | aggaaacatt | tttcttgttt | cgagtagata | atgccagcct | gttaaacgcc | 660  |
| gtcgac | gagt  | ctaacggaca | ccaaccagcg | aaccagcagc | gtcgcgtcgg | gccaagcgaa | 720  |
| gcagac | ggca  | cggcatctct | gtcgctgcct | ctggacccct | ctcgagagtt | ccgctccacc | 780  |
| gttgga | cttg  | ctccgctgtc | ggcatccaga | aattgcgtgg | cggagcggca | gacgtgagcc | 840  |
| ggcacg | gcag  | gcggcctcct | cctcctcta  | cggcacggca | gctacggggg | attcctttcc | 900  |
| caccgc | tcct  | tcgctttccc | ttcctcgccc | gccgtaataa | atagacaccc | cctccacacc | 960  |
| ctcttt | cccc  | aacctcgtgt | tgttcggagc | gcacacacac | acaaccagat | ctccccaaa  | 1020 |
| tccacc | cgtc  | ggcacctccg | cttcaaggta | cgccgctcgt | cctcccccc  | ccccctctc  | 1080 |
| tacctt | ctct  | agatcggcgt | tccggtccat | ggttagggcc | cggtagttct | acttctgttc | 1140 |
| atgttt | gtgt  | tagatccgtg | tttgtgttag | atccgtgctg | ctagcgttcg | tacacggatg | 1200 |
| cgacct | gtac  | gtcagacacg | ttctgattgc | taacttgcca | gtgtttctct | ttggggaatc | 1260 |
| ctggga | tggc  | tctagccgtt | ccgcagacgg | gatcgatttc | atgattttt  | ttgtttcgtt | 1320 |
| gcatag | ggtt  | tggtttgccc | ttttccttta | tttcaatata | tgccgtgcac | ttgtttgtcg | 1380 |
| ggtcat | cttt  | tcatgctttt | ttttgtcttg | gttgtgatga | tgtggtctgg | ttgggcggtc | 1440 |
| gttcta | gatc  | ggagtagaat | tctgtttcaa | actacctggt | ggatttatta | attttggatc | 1500 |
| tgtatg | tgtg  | tgccatacat | attcatagtt | acgaattgaa | gatgatggat | ggaaatatcg | 1560 |
| atctag | gata  | ggtatacatg | ttgatgcggg | ttttactgat | gcatatacag | agatgctttt | 1620 |
| tgttcg | cttg  | gttgtgatga | tgtggtgtgg | ttgggcggtc | gttcattcgt | tctagatcgg | 1680 |
| agtaga | atac  | tgtttcaaac | tacctggtgt | atttattaat | tttggaactg | tatgtgtgtg | 1740 |
| tcatac | atct  | tcatagttac | gagtttaaga | tggatggaaa | tatcgatcta | ggataggtat | 1800 |
| acatgt | tgat  | gtgggtttta | ctgatgcata | tacatgatgg | catatgcagc | atctattcat | 1860 |
| atgctc | taac  | cttgagtacc | tatctattat | aataaacaag | tatgttttat | aattattttg | 1920 |
| atcttg | atat  | acttggatga | tggcatatgc | agcagctata | tgtggatttt | tttagccctg | 1980 |
| ccttca | tacg  | ctatttattt | gcttggtact | gtttcttttg | tcgatgctca | ccctgttgtt | 2040 |
| tggtgt | tact  | tctgcaggtc | gactctagag | gatccaacaa | tggagaacaa | catacagaat | 2100 |
|        |       |            |            |            |            |            |      |

|                       | 03         |            |            |            |      | 04 |
|-----------------------|------------|------------|------------|------------|------|----|
|                       |            |            | -conti     | nued       |      |    |
| cagtgcgtcc cctacaactg | cctcaacaat | cctgaagtag | agattctcaa | cgaagagagg | 2160 |    |
|                       |            |            |            |            | 0000 |    |

togactggca gattgccgtt agacatotoc otgtocotta caogtttoot gttgtotgag 2220 2280 tttgttccag gtgtgggagt tgcgtttggc ctcttcgacc tcatctgggg cttcatcact 2340 ccatctgatt ggagcctctt tcttctccag attgaacagt tgattgaaca aaggattgag 2400 accttggaaa ggaatcgggc catcactacc cttcgtggct tagcagacag ctatgagatc 2460 tacattgaag cactaagaga gtgggaagcc aatcctaaca atgcccaact gagagaagat 2520 gtgcgtatac gctttgctaa cacagatgat gctttgatca cagccatcaa caacttcacc 2580 cttaccagct tcgagatccc tcttctctcg gtctatgttc aagctgctaa cctgcacttg 2640 tcactactgc gcgacgctgt gtcgtttggg caaggttggg gactggacat agctactgtc 2700 aacaatcact acaacagact catcaatctg attcatcgat acacgaaaca ttgtttggat 2760 acctacaatc agggattgga gaacctgaga ggtactaaca ctcgccaatg ggccaggttc 2820 aatcagttca ggagagacct tacacttact gtgttagaca tagttgctct ctttccgaac 2880 tacgatgttc gtacctatcc gattcaaacg tcatcccaac ttacaaggga gatctacacc agttcagtca ttgaagactc tccagtttct gcgaacatac ccaatggttt caacagggct 2940 3000 gagtttggag tcagaccacc ccatctcatg gacttcatga actctttgtt tgtgactgca gagactgtta gatcccaaac tgtgtgggga ggacacttag ttagctcacg caacacggct 3060 3120 ggcaatcgta tcaactttcc tagttacggg gtcttcaatc ccggggggcgc catctggatt 3180 gcagatgaag atccacgtcc tttctatcgg accttgtcag atcctgtctt cgtccgagga 3240 ggctttggca atcctcacta tgtactcggt cttaggggag tggcctttca acaaactggt 3300 acgaatcaca cccgcacatt caggaactcc gggaccattg actctctaga tgagatacca cctcaagaca acagcggcgc accttggaat gactactccc atgtgctgaa tcatgttacc 3420 tttgtgcgct ggccaggtga gatctcaggt tccgactcat ggagagcacc aatgttctct 3480 tggacgcatc gtagcgctac ccccacaaac accattgatc cagagagaat cactcagatt 3540 cccttggtga aggcacacac acttcagtca ggaactacag ttgtaagagg gccggggttc 3600 acgggaggag acattetteg acgeaetagt ggaggaeeat tegegtaeae eattgteaae 3660 atcaatgggc aacttcccca aaggtatcgt gccaggatac gctatgcctc tactaccaat 3720 ctaagaatct acgttacggt tgcaggtgaa cggatctttg ctggtcagtt caacaagaca 3780 atggataccg gtgatccact tacattccaa tctttctcct acgccactat caacaccgcg 3840 ttcacctttc caatgagcca gagcagtttc acagtaggtg ctgatacctt cagttcaggc 3900 aacgaagtgt acattgacag gtttgagttg attccagtta ctgccacact cgagtaagga 3960 teegtegaee tgeageeaag ettttegega getegagate eeegaeatat geeeeggttt 4020 cgttgcgact aacatgagtt cttggacaaa tttgattgga cctgatgaga tgatccaacc 4080 cgaggatata gcaaagctcg ttcgtgcagc aatggaacgg ccaaaccgtg cttttgtccc caagaatgag gtgctatgca tgaaggaatc tacccgttga tgtccaacag tctcagggtt aatgtctatg tatcttaaat aatgttgtcg gtattttgta atctcatata gattttcact 4200 gtgcgacgca aaaatattaa ataaatatta ttattatcta cgttttgatt gagatatcat 4260 caatattata ataaaaatat ccattaaaca cgatttgata caaatgacag tcaataatct 4320 gatttgaata tttattaatt gtaacgaatt acataaagat cgaatagaaa atactgcact 4380 gcaaatgaaa attaacacat actaataaat gcgtcaaata tctttgccaa gatcaagcgg 4440

-continued

| agtgagggcc  | tcatatccgg | tctcagttac | aagcacggta | tccccgaagc | gcgctccacc | 4500 |
|-------------|------------|------------|------------|------------|------------|------|
| aatgccctcg  | acatagatgc | cgggctcgac | gctgaggaca | ttgcctacct | tgagcatggt | 4560 |
| ctcagcgccg  | gctttaagct | caatcccatc | ccaatctgaa | tatcctatcc | cgcgcccagt | 4620 |
| ccggtgtaag  | aacgggtctg | tccatccacc | tctgttggga | attccggtcc | gggtcacctt | 4680 |
| tgtccaccaa  | gatggaactg | cggccgcgga | ccgaattccc | atggagtcaa | agattcaaat | 4740 |
| agaggaccta  | acagaactcg | ccgtaaagac | tggcgaacag | ttcatacaga | gtctcttacg | 4800 |
| actcaatgac  | aagaagaaaa | tcttcgtcaa | catggtggag | cacgacacgc | ttgtctactc | 4860 |
| caaaaatatc  | aaagatacag | tctcagaaga | ccaaagggca | attgagactt | ttcaacaaag | 4920 |
| ggtaatatcc  | ggaaacctcc | tcggattcca | ttgcccagct | atctgtcact | ttattgtgaa | 4980 |
| gatagtggaa  | aaggaaggtg | gctcctacaa | atgccatcat | tgcgataaag | gaaaggccat | 5040 |
| cgttgaagat  | gcctctgccg | acagtggtcc | caaagatgga | ccccaccca  | cgaggagcat | 5100 |
| cgtggaaaaa  | gaagacgttc | caaccacgtc | ttcaaagcaa | gtggattgat | gtgatatctc | 5160 |
| cactgacgta  | agggatgacg | cacaatccca | ctatccttcg | caagaccctt | cctctatata | 5220 |
| aggaagttca  | tttcatttgg | agaggacagg | gtacccgggg | atccaccatg | tctccggaga | 5280 |
| ggagaccagt  | tgagattagg | ccagctacag | cagctgatat | ggccgcggtt | tgtgatatcg | 5340 |
| ttaaccatta  | cattgagacg | tctacagtga | actttaggac | agagccacaa | acaccacaag | 5400 |
| agtggattga  | tgatctagag | aggttgcaag | atagataccc | ttggttggtt | gctgaggttg | 5460 |
| agggtgttgt  | ggctggtatt | gcttacgctg | ggccctggaa | ggctaggaac | gcttacgatt | 5520 |
| ggacagttga  | gagtactgtt | tacgtgtcac | ataggcatca | aaggttgggc | ctaggatcca | 5580 |
| cattgtacac  | acatttgctt | aagtctatgg | aggcgcaagg | ttttaagtct | gtggttgctg | 5640 |
| ttataggcct  | tccaaacgat | ccatctgtta | ggttgcatga | ggctttggga | tacacagccc | 5700 |
| ggggtacatt  | gcgcgcagct | ggatacaagc | atggtggatg | gcatgatgtt | ggtttttggc | 5760 |
| aaagggattt  | tgagttgcca | gctcctccaa | ggccagttag | gccagttacc | cagatctgag | 5820 |
| tcgacctgca  | ggcatgccgc | tgaaatcacc | agtctctctc | tacaaatcta | tctctcta   | 5880 |
| taataatgtg  | tgagtagttc | ccagataagg | gaattagggt | tcttataggg | tttcgctcat | 5940 |
| gtgttgagca  | tataagaaac | ccttagtatg | tatttgtatt | tgtaaaatac | ttctatcaat | 6000 |
| aaaatttcta  | attcctaaaa | ccaaaatcca | gtggcgagct | cgaattcgag | ctcgagcccg | 6060 |
| ggtggatcct  | ctagagtcga | cctgcagaag | cttcggtccg | gcgcgcctct | agttgaagac | 6120 |
| acgttcatgt  | cttcatcgta | agaagacact | cagtagtctt | cggccagaat | ggcctaactc | 6180 |
| aaggcc      |            |            |            |            |            | 6186 |
|             |            |            |            |            |            |      |
| <210> SEQ ] |            |            |            |            |            |      |

<400> SEQUENCE: 26

| 60  | ggtgccattc | gtttatggaa | ccccaccgac | cgtgccccta | tagcccgcgt | actagtttcc |
|-----|------------|------------|------------|------------|------------|------------|
| 120 | gtaaatttcc | gtaaaatccg | taaatggtcg | cctaaggatg | cgtggccgcc | cacggttctt |
| 180 | cagaaaacga | tcgggatata | ttcggattta | tccagccgtt | accagatttt | ggtaccgttt |
| 240 | acaaacttac | aaacggtgag | acggtacggt | tttttcgaaa | gaataggttt | gacggaaacg |

<sup>&</sup>lt;211> LENGTH: 3830

<sup>&</sup>lt;212> TYPE: DNA

<sup>&</sup>lt;213> ORGANISM: Artificial Sequence

<sup>&</sup>lt;220> FEATURE:

<sup>&</sup>lt;223> OTHER INFORMATION: Sequence that represents part of the PHI8999A insert as well as flanking sequence 5' to the insert.

|            |            | 67         | •          |            |             |      | 6 |
|------------|------------|------------|------------|------------|-------------|------|---|
|            |            |            |            | -contir    | iuea<br>——— |      |   |
| cgtccgtttt | cgtatttctc | gggaaactct | ggtatattcc | cgtatttgtc | ccgtattttc  | 300  |   |
| ccgacccacg | gacctgccaa | tcaaccatca | gccagtcagc | ccatccccac | agctatggcc  | 360  |   |
| catggggcca | tgttggccac | atgcccacgc | aacgcaaggc | agtaaggctg | gcagcctggc  | 420  |   |
| acgcattgac | gcatgtggac | acacacagcc | gccgcctgtt | cgtgtttctg | tgccgttgtg  | 480  |   |
| cgagactgtg | actgcgagtg | gcggagtcgg | cgaacggcga | ggcgtctccg | gagtctggac  | 540  |   |
| tgcggctgtg | gacagcgacg | ctgtgacggc | gactcggcga | agccccaagc | taccaagccc  | 600  |   |
| ccaagtcccc | atccatctct | gcttctctgg | tcatctcctt | cccctggtcg | atctgcaggc  | 660  |   |
| gccagaccgg | ccgaagcatc | acgaaacgca | ctaagacctc | gaaggagtca | aaccactcct  | 720  |   |
| ccgaggcctc | gggggctaca | cccggcgggt | gcgctcgcgc | gcacccaccg | gaacaaaatg  | 780  |   |
| taaccgagaa | aggtcggtcc | ccttgcaaaa | aaagtgcgac | aaaagcctcc | aagcgagtat  | 840  |   |
| taacactcac | tttgaggctc | gggggctact | gtcggggacc | ataattaggg | gtacccccaa  | 900  |   |
| gactcctaat | ctcagctggt | aacccccatc | agcacaaagc | tgcaaaggcc | tgatgggtgc  | 960  |   |
| gattaagtca | aggctcggtc | cactcaaggg | acacgatctc | gcctcgcccg | agcccagcct  | 1020 |   |
| cgggcaaggg | cggccgaccc | cgaggattca | cgtctcgccc | gagggcccc  | tcaagcgacg  | 1080 |   |
| ggcacacctt | cggctcgccc | gaggcccatt | cttcgccgag | aagcaacctt | ggccagatcg  | 1140 |   |
| ccacaccgac | cgaccgtatc | gcaggagcat | ttaatgcgag | gatcgcctga | caccttatcc  | 1200 |   |
| tgacgcgcgc | tcttcagtcg | acagagccga | agtgaccgca | atcacttcgc | cgctccactg  | 1260 |   |
| accgacctga | caagaagaca | gcgccgcctg | cgtcgctccg | actgctgtgc | cactcgacag  | 1320 |   |
| agtgaggctg | acagcagcca | agtccggcct | cgggcgccat | aggaagctcc | gcctcgcccg  | 1380 |   |
| accctagggc | tcggactcgg | cctcggctcc | ggaagacgac | gaactacgct | tcgcccgacc  | 1440 |   |
| ccagggcttg | gactcagcct | cggctccgga | agacgacgaa | ttccgcctcg | cccgacccca  | 1500 |   |
| gggctcggac | tcggcctcgg | ctccagaaga | cgacgaactc | cgcctcgccc | gaccccaggg  | 1560 |   |
| ctcggactca | gcctcggctc | cggaagacga | cgaactccgc | ctcgcccgac | cccagggctc  | 1620 |   |
| ggactcagcc | tcggcctcag | acgatggtct | ccgcctcgcc | cgacccgggg | ctcggactcg  | 1680 |   |
| acctttctat | cggaccttgt | cagatcctgt | cttcgtccga | ggaggctttg | gcaatcctca  | 1740 |   |
| ctatgtactc | ggtcttaggg | gagtggcctt | tcaacaaact | ggtacgaatc | acacccgcac  | 1800 |   |
| attcaggaac | tccgggacca | ttgactctct | agatgagata | ccacctcaag | acaacagcgg  | 1860 |   |
| cgcaccttgg | aatgactact | cccatgtgct | gaatcatgtt | acctttgtgc | gctggccagg  | 1920 |   |
| tgagatctca | ggttccgact | catggagagc | accaatgttc | tcttggacgc | atcgtagcgc  | 1980 |   |
| tacccccaca | aacaccattg | atccagagag | aatcactcat | tcttcaagaa | ctgcatatct  | 2040 |   |
| tgccgagatc | ctcatcccta | aaggtacttg | acaatagtat | tattggagtc | gatacacaac  | 2100 |   |
| tcacaaaaaa | tacaagaagt | cgactaggtg | gattggtccg | agtgaagaga | aaaaaagcc   | 2160 |   |
| atacagaact | caaaatcttt | tccggagata | ttcattttcc | tgaagaggcg | gataagatat  | 2220 |   |
| taggtggcag | tttgatacca | ccagaaagag | aaaaaaaga  | ttctaaggaa | tcaaaaaaaa  | 2280 |   |
| ggaaaaattg | ggtttatgtt | caacggaaaa | aatttctcaa | aagcaaggaa | aagtattgtg  | 2340 |   |
|            |            |            | ccgcggtttg |            |             | 2400 |   |
|            |            |            |            |            |             | 2460 |   |
|            |            |            | agccacaaac |            |             |      |   |
|            |            |            | ggttggttgc |            |             | 2520 |   |
| ctggtattgc | ttacgctggg | ccctggaagg | ctaggaaccc | tcaacctcag | caaccaacca  | 2580 |   |

atggtatcta tcttgcaacc tctctagatc atcaatccac tcttgtggtg tttgtggctc

2640

**68** 

-continued

| tgtcctaaag | ttcactgtag | acgtctcaat | gtaatggtta | acgatatcac | aaaccgagag | 2700 |
|------------|------------|------------|------------|------------|------------|------|
| aagagggatc | tcgaagcttc | ggccggggcc | catcgatatc | cgcgggcatg | cctgcagtgc | 2760 |
| agcgtgaccc | ggtcgtgccc | ctctctagag | ataatgagca | ttgcatgtct | aagttataaa | 2820 |
| aaattaccac | aactggaaga | gcggttaccc | ggaccgaagc | ttcggccggg | gcccatcgat | 2880 |
| atccgcgggc | atgcctgcag | tgcagcgtga | cccggtcgtg | cccctctcta | gagataatga | 2940 |
| gcattgcatg | tctaagttat | aaaaattac  | cacatatttt | ttttgtcaca | cttgtttgaa | 3000 |
| gtgcagttta | tctatcttta | tacatatatt | taaactttac | tctacgaata | atataatcta | 3060 |
| tagtactaca | ataatatcag | tgttttagag | aatcatataa | atgaacagtt | agacatggtc | 3120 |
| taaaggacaa | ttgagtattt | tgacaacagg | actctacagt | tttatcttt  | tagtgtgcat | 3180 |
| gtgttctcct | tttttttgc  | aaatagcttc | acctatataa | tacttcatcc | attttattag | 3240 |
| tacatccatt | tagggtttag | ggttaatggt | ttttatagac | taatttttt  | agtacatcta | 3300 |
| ttttattcta | ttttagcctc | taaattaaga | aaactaaaac | tctattttag | ttttttatt  | 3360 |
| taataattta | gatataaaat | agaataaaat | aaagtgacta | aaaattaaac | aaataccctt | 3420 |
| taagaaatta | aaaaactaa  | ggaaacattt | ttcttgtttc | gagtagataa | tgccagcctg | 3480 |
| ttaaacgccg | tcgacgagtc | taacggacac | caaccagcga | accagcagcg | tcgcgtcggg | 3540 |
| ccaagcgaag | cagacggcac | ggcatctctg | tcgctgcctc | tggacccctc | tcgagagttc | 3600 |
| cgctccaccg | ttggacttgc | tccgctgtcg | gcatccagaa | attgcgtggc | ggagcggcag | 3660 |
| acgtgagccg | gcacggcagg | cggcctcctc | ctcctctcac | ggcacggcag | ctacggggga | 3720 |
| ttcctttccc | accgctcctt | cgctttccct | tcctcgcccg | ccgtaataaa | tagacacccc | 3780 |
| ctccacaccc | tctttcccca | acctcgtgtt | gttcggagcg | cacacacaca |            | 3830 |

<sup>&</sup>lt;210> SEQ ID NO 27

## <400> SEQUENCE: 27

| cccactatcc | ttcgcaagac | ccttcctcta | tataaggaag | ttcatttcat | ttggagagga | 60  |
|------------|------------|------------|------------|------------|------------|-----|
| cagggtaccc | ggggatccac | catgtctccg | gagaggagac | cagttgagat | taggccagct | 120 |
| acagcagctg | atatggccgc | ggtttgtgat | atcgttaacc | attacattga | gacgtctaca | 180 |
| gtgaacttta | ggacagagcc | acaaacacca | caagagtgga | ttgatgatct | agagaggttg | 240 |
| caagatagat | acccttggtt | ggttgctgag | gttgagggtg | ttgtggctgg | tattgcttac | 300 |
| gctgggccct | ggaaggctag | gaacgcttac | gattggacag | ttgagagtac | tgtttacgtg | 360 |
| tcacataggc | atcaaaggtt | gggcctagga | tccacattgt | acacacattt | gcttaagtct | 420 |
| atggaggcgc | aaggttttaa | gtctgtggtt | gctgttatag | gccttccaaa | cgatccatct | 480 |
| gttaggttgc | atgaggcttt | gggatacaca | gcccggggta | cattgcgcgc | agctggatac | 540 |
| aagcatggtg | gatggcatga | tgttggtttt | tggcaaaggg | attttgagtt | gccagctcct | 600 |
| ccaaggccag | ttaggccagt | tacccagatc | tgagtcgacc | tgcaggcatg | cccgctgaaa | 660 |
| tcaccagtct | ctctctacaa | atctatctct | ctctataata | atgtgtgagt | agttcccaga | 720 |
| taagggaatt | agggttctta | tagggtttcg | ctcatgtgtt | gagcatataa | gaaaccctta | 780 |
| gtatgtattt | gtatttgtaa | aatacttcta | tcaataaaat | ttctaattcc | taaaaccaaa | 840 |

<sup>&</sup>lt;211> LENGTH: 3347

<sup>&</sup>lt;212> TYPE: DNA

<sup>&</sup>lt;213> ORGANISM: Artificial Sequence

<sup>&</sup>lt;220> FEATURE:

<sup>&</sup>lt;223> OTHER INFORMATION: Sequence that represents part of the PHI8999A
insert as well as flanking sequence 3' to the insert.

-continued

**71** 

| atccagtggc | gagctcgaat | tcgagctcga | gcccgggtgg | atcctctaga | gtcgacctgc | 900  |
|------------|------------|------------|------------|------------|------------|------|
| agaagcttcg | gtccggcgcg | cctctagttg | aagacacgtt | catgtcttca | tcgtaagaag | 960  |
| acactcagta | gtcttcggcc | agaatggcct | aactcaaggc | cctcactccg | cttgatcttg | 1020 |
| gcaaagatat | ttgacgcatt | tattagtatg | tgttaatttt | catttgcagt | gcagtatttt | 1080 |
| ctattcgatc | tttatgtaat | tcgttacaat | taataaatat | tcaaatcaga | ttattgactg | 1140 |
| tcatttgtat | caaatcgtgt | ttaatggata | tttttattat | aatattgatg | atatctcaat | 1200 |
| caaaacgtag | ataataataa | tatttattta | atatttttgc | gtcgcacagt | gaaaatctat | 1260 |
| atgagattac | aaaataccga | caacattatt | taagaaacat | agacattaac | cctgagactg | 1320 |
| ttggacatca | acgggtagat | tccttcatgc | atagcacctc | attcttgggg | acaaaagcac | 1380 |
| ggtttggccg | ttccattgct | gcacgaacga | gctttgctat | atcctcgggt | tggatcatct | 1440 |
| catcaggtcc | aatcaaattt | gtccaagaac | tcatgttagt | cgcaacgaaa | ccggggcata | 1500 |
| tgtcgggtat | ctcgagctcg | cgaaagcttg | gctgcaggtc | gacggatcct | tcaacaaaag | 1560 |
| ggtacctgta | cccgaaaccg | acacaggtgg | gtaggtagag | aatacctagg | ggcgcgagac | 1620 |
| aactctctct | aaggaactcg | gcaaaatagc | cccgtaactt | cgggagaagg | ggtgccccc  | 1680 |
| gctaacaata | aacgaatacg | gtttatgtat | ggattccggt | aaaataccgg | tactcgattt | 1740 |
| cataagagtc | gaataggaag | ttaagatgag | ggtggtatca | tcataaaaat | ggagtagtat | 1800 |
| cctaaattat | actaatccac | gtatgatatg | tatgcctttc | cttatcaacc | ggaagtagtg | 1860 |
| caaaaaaaat | tctatactgc | actgctctct | ttttactgag | aaatgcaaaa | aaataaaagt | 1920 |
| gaagtaaggg | tgccccatag | atatttgatc | ttgcctcctg | tccccccc   | ccttttttca | 1980 |
| tcaaaaattt | ccatgaaaaa | agaaaagatg | aatttgtcca | ttcattgaac | cctagttcgg | 2040 |
| gactgacggg | gctcgaaccc | gcagcttccg | cctgttccta | gccttccagg | gcccagcgta | 2100 |
| agcaatacca | gccacagcac | cctcaacctc | agcaaccaac | caagggtatc | tatcttgcaa | 2160 |
| cctctctaga | tcatcaatcc | actcttgtgg | tgtttgtggc | tctgtcctaa | agttcactgt | 2220 |
| agacgtctca | atgtaatggt | taacgatatc | acaaaccgcg | gaacacaaga | acgaaagcac | 2280 |
| cttttcattc | tttcatatac | taggggtttt | tacttggaaa | agacaatgtt | ccatactaaa | 2340 |
| ggatagctgc | agaagccgcc | accgtcttga | ggaccttccg | gggagccaga | ccggtcgaac | 2400 |
| cgtgcctcca | cttgctaagg | agaaagggaa | aatcagggcc | aggacatacg | aaggaggagc | 2460 |
| cagaacgaag | atatcctaag | atacttactc | gctccgggcc | atgatcaatc | atgcctgtgg | 2520 |
| ggaggtctct | cgcacctcga | tccatgaagg | taccaccgag | gtctgccccg | ccgccggctt | 2580 |
| cggtaccgtc | ctcgccttgg | gcgcccgagg | cacccggggg | atggactgcc | caggcgcagc | 2640 |
| cacgacgacc | caaggatcac | cctcctgcgc | agtcggcacg | agcaatagtt | ctcggggaac | 2700 |
| aggcagcttg | gcctgactcc | ccggggtcac | ctcaactacc | tcggccgagg | ggtcaagtac | 2760 |
| cccctcagtc | cgcccccgct | cttcggaccg | ggaccccgac | gtcccggccc | cggataccga | 2820 |
| cggcaccagc | ccgctcgggg | gctggcttga | cgacccctgg | cccagcctca | gatctgggct | 2880 |
| gaggccgagg | caggcggcca | tgtcgtcgtc | ttcatcatcg | tcttcatcat | cgtcgtcgtc | 2940 |
| atcaggcgtc | tccggcgacg | gctcccttgg | gagcccctcc | ctctcctgcc | gacgacgaag | 3000 |
| cctttccaag | gcatcccgag | cccacgtccg | ctcgtgggcc | cgagccttct | ttgcgtcctt | 3060 |
| cttctccttc | ctcttctccg | cggtgaccct | ccgcgcagct | cggtccaccg | catcctccgg | 3120 |
| gactggtggc | agggaaggct | tgtgatgccc | tacctcctgg | agacagacga | aaagtctcag | 3180 |
|            |            |            |            |            |            |      |

**72** 

|  | 70              |              |              |             |      |  |
|--|-----------------|--------------|--------------|-------------|------|--|
|  |                 |              | -contir      | nued        |      |  |
| ctatgagaac cgagggcaat  | ctgacgcaag      | aaggaagaag   | gagcggatac   | tcaccagaga  | 3240 |  |
| acgcacccg cgatcgggad   | gcattaaggg      | ctgggaaaaa   | gtgccggcct   | ctaatttcgc  | 3300 |  |
| accgtgccg tccacccacc   | tgtggaggtc      | atcgatggga   | aggggaa      |             | 3347 |  |
| <210> SEQ ID NO 28<br><211> LENGTH: 669<br><212> TYPE: DNA<br><213> ORGANISM: Zea m  | nays            |              |              |             |      |  |
| <400> SEQUENCE: 28   |                 |              |              |             |      |  |
| actagtttcc tagcccgcgt  | cgtgccccta      | ccccaccgac   | gtttatggaa   | ggtgccattc  | 60   |  |
| cacggttctt cgtggccgcc  | cctaaggatg      | taaatggtcg   | gtaaaatccg   | gtaaatttcc  | 120  |  |
| ggtaccgttt accagatttt  | tccagccgtt      | ttcggattta   | tcgggatata   | cagaaaacga  | 180  |  |
| gacggaaacg gaataggttt  | tttttcgaaa      | acggtacggt   | aaacggtgag   | acaaacttac  | 240  |  |
| cgtccgtttt cgtatttctc  | gggaaactct      | ggtatattcc   | cgtatttgtc   | ccgtattttc  | 300  |  |
| ccgacccacg gacctgccaa  | tcaaccatca      | gccagtcagc   | ccatccccac   | agctatggcc  | 360  |  |
| catggggcca tgttggccac  | atgcccacgc      | aacgcaaggc   | agtaaggctg   | gcagcctggc  | 420  |  |
| acgcattgac gcatgtggac  | acacacagcc      | gccgcctgtt   | cgtgtttctg   | tgccgttgtg  | 480  |  |
| cgagactgtg actgcgagtg  | gcggagtcgg      | cgaacggcga   | ggcgtctccg   | gagtctggac  | 540  |  |
| tgcggctgtg gacagcgacg  | ctgtgacggc      | gactcggcga   | agccccaagc   | taccaagccc  | 600  |  |
| ccaagtcccc atccatctct  | gcttctctgg      | tcatctcctt   | cccctggtcg   | atctgcaggc  | 660  |  |
| ccagaccg   |                 |              |              |             | 669  |  |
| <210> SEQ ID NO 29<br><211> LENGTH: 200<br><212> TYPE: DNA<br><213> ORGANISM: Zea m  | nays            |              |              |             |      |  |
| <400> SEQUENCE: 29   |                 |              |              |             |      |  |
| gccgaagcat cacgaaacgo  | actaagacct      | cgaaggagtc   | aaaccactcc   | tccgaggcct  | 60   |  |
| cgggggctac acccggcggg  | tgcgctcgcg      | cgcacccacc   | ggaacaaaat   | gtaaccgaga  | 120  |  |
| aggtcggtc cccttgcaaa   | aaaagtgcga      | caaaagcctc   | caagcgagta   | ttaacactca  | 180  |  |
| etttgagget egggggetae  | !               |              |              |             | 200  |  |
| <pre>&lt;210&gt; SEQ ID NO 30 &lt;211&gt; LENGTH: 812 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Zea m &lt;220&gt; FEATURE: &lt;221&gt; NAME/KEY: misc_ &lt;222&gt; LOCATION: (0)</pre> | feature<br>.(0) | t of maize 1 | Huck-1 retro | otransposon |      |  |
| <400> SEQUENCE: 30   |                 |              |              |             |      |  |
| tgtcgggac cataattagg   | ggtaccccca      | agactcctaa   | tctcagctgg   | taacccccat  | 60   |  |
| cagcacaaag ctgcaaaggo  | ctgatgggtg      | cgattaagtc   | aaggctcggt   | ccactcaagg  | 120  |  |
| gacacgatct cgcctcgccc  | gagcccagcc      | tcgggcaagg   | gcggccgacc   | ccgaggattc  | 180  |  |
| acgtctcgcc cgagggcccc  | ctcaagcgac      | gggcacacct   | tcggctcgcc   | cgaggcccat  | 240  |  |
| tcttcgccga gaagcaacct  | tggccagatc      | gccacaccga   | ccgaccgtat   | cgcaggagca  | 300  |  |
| ttaatgcga ggatcgcctg   | , acaccttatc    | ctgacgcgcg   | ctcttcagtc   | gacagagccg  | 360  |  |
| aadgaccgc aatcacttco   |                 | anaanaata    | 242242424    | agggggggt   | 420  |  |

aagtgaccgc aatcacttcg ccgctccact gaccgacctg acaagaagac agcgccgcct

420

-continued

```
gcgtcgctcc gactgctgtg ccactcgaca gagtgaggct gacagcagcc aagtccggcc
                                                                     480
                                                                     540
tegggegeea taggaagete egeetegeee gaeeetaggg eteggaeteg geeteggete
                                                                     600
cggaagacga cgaactacgc ttcgcccgac cccagggctt ggactcagcc tcggctccgg
                                                                     660
aagacgacga attccgcctc gcccgacccc agggctcgga ctcggcctcg gctccagaag
                                                                     720
acgacgaact ccgcctcgcc cgaccccagg gctcggactc agcctcggct ccggaagacg
                                                                     780
acgaactccg cctcgcccga ccccagggct cggactcagc ctcggcctca gacgatggtc
                                                                     812
teegeetege eegaeeeggg geteggaete ga
<210> SEQ ID NO 31
<211> LENGTH: 335
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Represents part of PHI8999A insert sequence -
      fragment of Cry1F gene
<400> SEQUENCE: 31
                                                                      60
cettetate ggacettgte agateetgte ttegteegag gaggetttgg caateeteae
                                                                     120
tatgtactcg gtcttagggg agtggccttt caacaaactg gtacgaatca cacccgcaca
                                                                     180
ttcaggaact ccgggaccat tgactctcta gatgagatac cacctcaaga caacagcggc
                                                                     240
gcaccttgga atgactactc ccatgtgctg aatcatgtta cctttgtgcg ctggccaggt
                                                                     300
gagateteag gtteegaete atggagagea eeaatgttet ettggaegea tegtageget
                                                                     335
acccccacaa acaccattga tccagagaga atcac
<210> SEQ ID NO 32
<211> LENGTH: 321
<212> TYPE: DNA
<213 > ORGANISM: Zea mays
<220> FEATURE:
<221> NAME/KEY: misc_feature
<222> LOCATION: (0) ...(0)
<223> OTHER INFORMATION: Fragment of maize chloroplast rpoC2 gene
<400> SEQUENCE: 32
                                                                      60
tcattcttca agaactgcat atcttgccga gatcctcatc cctaaaggta cttgacaata
                                                                     120
gtattattgg agtcgataca caactcacaa aaaatacaag aagtcgacta ggtggattgg
tccgagtgaa gagaaaaaa agccatacag aactcaaaat cttttccgga gatattcatt
                                                                     180
                                                                     240
ttcctgaaga ggcggataag atattaggtg gcagtttgat accaccagaa agagaaaaaa
                                                                     300
aagattctaa ggaatcaaaa aaaaggaaaa attgggttta tgttcaacgg aaaaaatttc
                                                                     321
tcaaaagcaa ggaaaagtat t
<210> SEQ ID NO 33
<211> LENGTH: 17
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Represents part of PHI8999A insert
      sequence-fragment of ubiZM1(2) promoter; also a
      fragment of the maize chloroplast trn1 gene
<400> SEQUENCE: 33
                                                                      17
gtggctattt atctatc
```

<210> SEQ ID NO 34 <211> LENGTH: 201 -continued

77

```
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Represents part of PHI8999A insert
      sequence-fragment of pat gene
<400> SEQUENCE: 34
gcagctgata tggccgcggt ttgtgatatc gttaaccatt acattgagac gtctacagtg
                                                                     120
aactttagga cagagccaca aacaccacaa gagtggattg atgatctaga gaggttgcaa
                                                                     180
gatagatacc cttggttggt tgctgaggtt gagggtgttg tggctggtat tgcttacgct
                                                                      201
gggccctgga aggctaggaa c
<210> SEQ ID NO 35
<211> LENGTH: 138
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Represents part of PHI8999A insert
      sequence-fragment of pat gene (complement)
<400> SEQUENCE: 35
cctcaacctc agcaaccaac caatggtatc tatcttgcaa cctctctaga tcatcaatcc
                                                                      60
actcttgtgg tgtttgtggc tctgtcctaa agttcactgt agacgtctca atgtaatggt
                                                                     120
                                                                     138
taacgatatc acaaaccg
<210> SEQ ID NO 36
<211> LENGTH: 15
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Represents part of PHI8999A insert
      sequence-fragment of cry1F gene (complement)
<400> SEQUENCE: 36
                                                                      15
agagaagagg gatct
<210> SEQ ID NO 37
<211> LENGTH: 118
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Represents part of PHI8999A insert
      sequence-fragment of polylinker
<400> SEQUENCE: 37
                                                                      60
cgaagetteg geegggeee ategatatee gegggeatge etgeagtgea gegtgaeeeg
gtcgtgcccc tctctagaga taatgagcat tgcatgtcta agttataaaa aattacca
                                                                     118
<210> SEQ ID NO 38
<211> LENGTH: 550
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Represents part of PHI8999A insert
      sequence-fragment of ORF25 terminator (complement)
<400> SEQUENCE: 38
                                                                      60
ctcactccgc ttgatcttgg caaagatatt tgacgcattt attagtatgt gttaattttc
atttgcagtg cagtattttc tattcgatct ttatgtaatt cgttacaatt aataaatatt
                                                                     120
                                                                     180
caaatcagat tattgactgt catttgtatc aaatcgtgtt taatggatat ttttattata
                                                                     240
atattgatga tatctcaatc aaaacgtaga taataataat atttatttaa tatttttgcg
```

```
300
tcgcacagtg aaaatctata tgagattaca aaataccgac aacattattt aagaaacata
                                                                   360
gacattaacc ctgagactgt tggacatcaa cgggtagatt ccttcatgca tagcacctca
                                                                   420
ttcttgggga caaaagcacg gtttggccgt tccattgctg cacgaacgag ctttgctata
                                                                   480
teetegggtt ggateatete ateaggteea ateaaatttg teeaagaaet eatgttagte
                                                                   540
gcaacgaaac cggggcatat gtcgggtatc tcgagctcgc gaaagcttgg ctgcaggtcg
                                                                   550
acggatcctt
<210> SEQ ID NO 39
<211> LENGTH: 128
<212> TYPE: DNA
<213 > ORGANISM: Zea mays
<220> FEATURE:
<221> NAME/KEY: misc_feature
<222> LOCATION: (0) ...(0)
<223> OTHER INFORMATION: Fragment of maize chloroplast rps12 rRNA gene
      (complement)
<400> SEQUENCE: 39
                                                                    60
caacaaaagg gtacctgtac ccgaaaccga cacaggtggg taggtagaga atacctaggg
                                                                   120
gcgcgagaca actctctcta aggaactcgg caaaatagcc ccgtaacttc gggagaaggg
                                                                   128
gtgccccc
<210> SEQ ID NO 40
<211> LENGTH: 392
<212> TYPE: DNA
<213 > ORGANISM: Zea mays
<220> FEATURE:
<221> NAME/KEY: misc_feature
<222> LOCATION: (0)...(0)
<223> OTHER INFORMATION: Fragment of maize chloroplast genome
<400> SEQUENCE: 40
ctaacaataa acgaatacgg tttatgtatg gattccggta aaataccggt actcgatttc
                                                                    60
ataagagtcg aataggaagt taagatgagg gtggtatcat cataaaaatg gagtagtatc
                                                                   120
                                                                   180
ctaaattata ctaatccacg tatgatatgt atgcctttcc ttatcaaccg gaagtagtgc
                                                                   240
aaaaaaaatt ctatactgca ctgctctctt tttactgaga aatgcaaaaa aataaaagtg
300
                                                                   360
caaaaatttc catgaaaaaa gaaaagatga atttgtccat tcattgaacc ctagttcggg
                                                                   392
actgacgggg ctcgaacccg cagcttccgc ct
<210> SEQ ID NO 41
<211> LENGTH: 188
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Represents part of PHI8999A insert
     sequence-fragment of pat gene (complement)
<400> SEQUENCE: 41
                                                                    60
gttcctagcc ttccagggcc cagcgtaagc aataccagcc acagcaccct caacctcagc
                                                                   120
aaccaaccaa gggtatctat cttgcaacct ctctagatca tcaatccact cttgtggtgt
                                                                   180
ttgtggctct gtcctaaagt tcactgtaga cgtctcaatg taatggttaa cgatatcaca
                                                                   188
aaccgcgg
```

-continued

**81 82** 

```
<211> LENGTH: 81
<212> TYPE: DNA
<213> ORGANISM: Zea mays
<220> FEATURE:
<221> NAME/KEY: misc_feature
<222> LOCATION: (0) ...(0)
<223> OTHER INFORMATION: Fragment of maize chloroplast ORF241
      (complement)
<400> SEQUENCE: 42
cacaagaacg aaagcacctt ttcattcttt catatactag gggtttttac ttggaaaaga
                                                                      60
                                                                      81
caatgttcca tactaaagga t
<210> SEQ ID NO 43
<211> LENGTH: 254
<212> TYPE: DNA
<213 > ORGANISM: Zea mays
<400> SEQUENCE: 43
                                                                      60
agctgcagaa gccgccaccg tcttgaggac cttccgggga gccagaccgg tcgaaccgtg
                                                                     120
cctccacttg ctaaggagaa agggaaaatc agggccagga catacgaagg aggagccaga
                                                                     180
acgaagatat cctaagatac ttactcgctc cgggccatga tcaatcatgc ctgtggggag
                                                                     240
gtctctcgca cctcgatcca tgaaggtacc accgaggtct gccccgccgc cggcttcggt
                                                                     254
accgtcctcg cctt
<210> SEQ ID NO 44
<211> LENGTH: 749
<212> TYPE: DNA
<213 > ORGANISM: Zea mays
<400> SEQUENCE: 44
gggcgcccga ggcacccggg ggatggactg cccaggcgca gccacgacga cccaaggatc
                                                                     120
accetectge geagteggea egageaatag ttetegggga acaggeaget tggeetgaet
                                                                     180
ccccggggtc acctcaacta cctcggccga ggggtcaagt acccctcag tccgccccg
                                                                     240
ctcttcggac cgggaccccg acgtcccggc cccggatacc gacggcacca gcccgctcgg
                                                                     300
gggctggctt gacgacccct ggcccagcct cagatctggg ctgaggccga ggcaggcggc
                                                                     360
catgtcgtcg tcttcatcat cgtcttcatc atcgtcgtcg tcatcaggcg tctccggcga
                                                                     420
eggeteeett gggageeet eeeteteetg eegaegaega ageettteea aggeateeeg
                                                                     480
agcccacgtc cgctcgtggg cccgagcctt ctttgcgtcc ttcttctcct tcctcttctc
                                                                     540
cgcggtgacc ctccgcgcag ctcggtccac cgcatcctcc gggactggtg gcagggaagg
                                                                     600
cttgtgatgc cctacctcct ggagacagac gaaaagtctc agctatgaga accgagggca
                                                                     660
atctgacgca agaaggaaga aggagcggat actcaccaga gacacgcacc cgcgatcggg
                                                                     720
acgcattaag ggctgggaaa aagtgccggc ctctaatttc gctaccgtgc cgtccaccca
                                                                     749
cctgtggagg tcatcgatgg gaaggggaa
<210> SEQ ID NO 45
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 5' flanking
      sequence; junction between regions 3 and 4
```

<400> SEQUENCE: 45

-continued

```
<210> SEQ ID NO 46
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 5' flanking
      sequence; junction between regions 4 and 5
<400> SEQUENCE: 46
agagaatcac tcattcttca
                                                                       20
<210> SEQ ID NO 47
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 5' flanking
      sequence; junction between regions 5 and 6
<400> SEQUENCE: 47
                                                                      20
gaaaagtatt gtggctattt
<210> SEQ ID NO 48
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 5' flanking
      sequence; junction between regions 6 and 7a
<400> SEQUENCE: 48
                                                                      20
tctcaaggcc gcagctgata
<210> SEQ ID NO 49
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 5' flanking
      sequence; junction between regions 7a and 7b
<400> SEQUENCE: 49
ggctaggaac cctcaacctc
                                                                      20
<210> SEQ ID NO 50
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 5' flanking
      sequence; junction between regions 7b and 7c
<400> SEQUENCE: 50
                                                                       20
tcacaaaccg agagaagagg
<210> SEQ ID NO 51
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 5' flanking
      sequence; junction between regions 7c and 8
<400> SEQUENCE: 51
```

agagggatct cgaagcttcg

```
<210> SEQ ID NO 52
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 5' flanking
      sequence; junction between regions 8 and 9
<400> SEQUENCE: 52
                                                                      20
aaaattacca caactggaag
<210> SEQ ID NO 53
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 3' flanking
      sequence; junction between regions 9 and 10
<400> SEQUENCE: 53
agctatgttt ctcactccgc
                                                                      20
<210> SEQ ID NO 54
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 3' flanking
      sequence; junction between regions 10 and 11
<400> SEQUENCE: 54
                                                                      20
acggatcctt caacaaagg
<210> SEQ ID NO 55
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 3' flanking
      sequence; junction between regions 11 and 12
<400> SEQUENCE: 55
                                                                      20
gtgcccccg ctaacaataa
<210> SEQ ID NO 56
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 3' flanking
      sequence; junction between regions 12 and 13
<400> SEQUENCE: 56
                                                                      20
gcttccgcct gttcctagcc
<210> SEQ ID NO 57
<211> LENGTH: 20
<212> TYPE: DNA
<213 > ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Biolistic transformation event 3' flanking
      sequence; junction between regions 13 and 14
<400> SEQUENCE: 57
                                                                      20
aaaccgcgga acacaagaac
```

What is claimed is:

- 1. A method for screening a seed lot for the presence of corn event TC1507 event, the method comprising: performing nucleic acid amplification of a sample of said seed lot using a specific primer pair, wherein one primer hybridizes to corn genomic DNA flanking the insertion site of PHI8999A of corn event TC1507, and the other primer hybridizes either to a sequence of PHI8999A or hybridizes to corn genomic DNA on the opposite side of the insertion site of PHI8999A of corn event TC1507, wherein said 10 specific primer pair is capable of producing an amplicon diagnostic of corn event TC1507; and detecting said diagnostic amplicon.
- 2. The method of claim 1, wherein one of the primers of the specific primer pair recognizes a sequence within SEQ 15 ID NO: 21.
- 3. The method of claim 2, wherein one of the primers of the specific primer pair comprises the sequence of SEQ ID NO: 1, SEQ ID NO: 2, or SEQ ID NO: 23.
- 4. The method of claim 1, wherein one of the primers of 20 the specific primer pair recognizes a sequence within SEQ ID NO: 22.
- 5. The method of claim 4, wherein one of the primers of the specific primer pair comprises the sequence of SEQ ID NO: 3, SEQ ID NO: 4, or SEQ ID NO: 5.
- 6. The method of claim 1, wherein the amplicon is generated from the specific primer pair: SEQ ID NO: 1 and SEQ ID NO: 2.

\* \* \* \* \*