



US006372430B1

(12) **United States Patent**  
**Morrison et al.**(10) **Patent No.:** US 6,372,430 B1  
(45) **Date of Patent:** Apr. 16, 2002(54) **NUCLEIC ACIDS FOR DETECTING ASPERGILLUS SPECIES AND OTHER FILAMENTOUS FUNGI**(75) Inventors: **Christine J. Morrison**, Decatur; **Errol Reiss**, Chamblee, both of GA (US); **Liliana Aidorevich**, Maracay Edo Aragun (VE); **Jong Soo Choi**, Taegu (KR)(73) Assignee: **The United States of America as represented by the Department of Health and Human Services**, Washington, DC (US)

(\*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 0 days.

(21) Appl. No.: **09/423,233**(22) PCT Filed: **May 1, 1998**(86) PCT No.: **PCT/US98/08926**§ 371 Date: **Jun. 27, 2000**§ 102(e) Date: **Jun. 27, 2000**(87) PCT Pub. No.: **WO98/50584**PCT Pub. Date: **Nov. 12, 1998****Related U.S. Application Data**

(60) Provisional application No. 60/045,400, filed on May 2, 1997.

(51) **Int. Cl.<sup>7</sup>** ..... **C12Q 1/68**; C07H 21/04; C12P 19/34(52) **U.S. Cl.** ..... **435/6**; 435/91.1; 536/23.1; 536/24.32; 536/24.3; 536/23.7(58) **Field of Search** ..... 435/6, 91.1; 536/23.1, 536/24.32, 24.3, 23.7(56) **References Cited**

## U.S. PATENT DOCUMENTS

5,426,027 A	6/1995	Lott et al.	.....	435/6
5,585,238 A	* 12/1996	Ligon et al.	.....	435/6
5,958,693 A	* 9/1999	Sanhu et al.	.....	435/6

## FOREIGN PATENT DOCUMENTS

WO 96/21741 7/1996

## OTHER PUBLICATIONS

Geiser et al. Genbank Accession No. L76747, Nov. 1996.\*  
LoBuglio et al. "Phylogeny and PCR identification of the human pathogenic fungus Penicillium marneffei" J. of Clinical Microbiology, vol. 33, No. 1, p. 85-89, Jan. 1995.\*

Geiser et al. Genbank Accession No. L76748, Nov. 1996.\*

Geiser et al. Genbank Accession No. L76774, Nov. 1996.\*

Borsuk et al. Genbank Accession No. U03521, Feb. 1995.\*

Borsuk et al. Genbank Accession No. U03523, Feb. 1995.\*

Borsuk et al. Genbank Accession No. U03519, Feb. 1995.\*

Kumeda et al. "Single-Strand Conformation polymorphism analysis of PCR-amplified ribosomal DNA ITS to differentiate species of Aspergillus section flavi" Applied Environ. Microbiology, vol. 62, No. 8, p. 2947-2952, Aug. 1996.\*

White et al. "Amplification and Direct sequencing of fungal ribosomal RNA genes for Phylogenetics" PCR Protocols: A guide to Methods and Applications, p. 315-322, Dec. 1989.\*

Lu, J-J, et al., "Typing of *Pneumocystis carinii* Strains with Type-Specific Oligonucleotide Probes Derived from Nucleotide Sequences of Internal Transcribed Spacers of rRNA Genes," *J. Clinical Microbiol.*, vol. 33, No. 11, p. 2973-2977, (Nov. 1995).Gaskell, G.J. et al., "Analysis of the internal transcribed spacer regions of ribosomal DNA in common airborne allergenic fungi," *Electrophoresis*, vol. 18, pp. 1567-1569 (1997).

\* cited by examiner

*Primary Examiner*—Lisa B. Arthur*Assistant Examiner*—Jeanine Goldberg(74) *Attorney, Agent, or Firm*—Klarquist Sparkman, LLC(57) **ABSTRACT**

Nucleic acids for detecting Aspergillus species and other filamentous fungi are provided. Unique internal transcribed space 2 coding regions permit the development of nucleic acid probes specific for five different species of Aspergillus, three species of Fusarium, four species of Mucor, two species of Penecillium, five species of Rhizopus, one species of Rhizomucor, as well as probes for *Absidia corymbifer*, *Cunninghamella elagans*, *Pseudallescheria boydii*, and *Sporothrix schenckii*. Methods are disclosed for the species-specific detection and diagnosis of infection by Aspergillus, Fusarium, Mucor, Penecillium, Rhizomucor, *absidia*, Cunninghamella, Pseudallescheria or Sporthrix in a subject. Furthermore, genus-specific probes are also provided for Aspergillus, Fusarium and Mucor, in addition to an all-fungus nucleic acid probe.

**29 Claims, No Drawings**

**1**

**NUCLEIC ACIDS FOR DETECTING  
ASPERGILLUS SPECIES AND OTHER  
FILAMENTOUS FUNGI**

**PRIORITY CLAIM**

This application claims priority to PCT/US98/08926, filed May 1, 1998, which claims the benefit of U.S. Provisional Application No. 60/045,400, filed May 2, 1997.

This invention was made in the Centers for Disease Control Mycotic Diseases Laboratories, an agency of the United States Government.

**TECHNICAL FIELD**

This application relates in general to the field of diagnostic microbiology. In particular, the invention relates to the species-specific detection of Aspergillus, Fusarium, Mucor, Penicillium, Rhizopus, Rhizomucor, Absidia, Cunninghamella, *Pseudallescheria boydii* (*Scedosporium apiospermum*), and Sporothrix species.

**BACKGROUND OF THE INVENTION**

In recent years, chemotherapy for hematological malignancies, and high-dose corticosteroid treatment for organ transplant recipients, along with the spread of AIDS, have greatly increased the number of immunocompromised patients (1, 12, 14, 43). Saprophytic filamentous fungi, such as Aspergillus, Rhizopus, and Mucor species, found in the environment and considered to be of low virulence, are now responsible for an increasing number of infections in the immunocompromised host (17, 20, 43). In addition, these infections are often fulminant and rapidly fatal in immunocompromised patients (7, 11, 12, 20, 44). Morbidity and mortality is extremely high; for example, aspergillosis has a mortality rate of approximately 90% (8, 11).

To complicate matters, diagnosis is difficult and symptoms are often non-specific (18, 27, 29, 42, 44). Antibody-based tests can be unreliable due to the depressed or variable immune responses of immunocompromised patients (2, 9, 18, 46). Antigen detection tests developed to date have fallen short of the desired sensitivity (2, 9, 38). Radiographic evidence can be non-specific and inconclusive (5, 29, 36), although some progress in diagnosis has been made with the advent of computerized tomography (40). However, definitive diagnosis still requires either a positive blood or tissue culture or histopathological confirmation (3, 21). An added complication is that the invasive procedures necessary to obtain biopsy materials are often not recommended in thrombocytopenic patient populations (37, 41).

Even when cultures of blood, lung or rhinocerebral tissues are positive, morphological and biochemical identification of filamentous fungi can require several days for adequate growth and sporulation to occur, delaying targeted drug therapy. Some atypical isolates may never sporulate, making identification even more difficult (23). When histopathology is performed on tissue biopsy sections, the morphological similarities of the various filamentous fungi in tissue make differentiation difficult (16). Fluorescent antibody staining of histopathological tissue sections is not specific unless cross-reactive epitopes are absorbed out which can make the resultant antibody reactions weak (14, 19). Therapeutic choices vary (7, 41, 44) making a test to rapidly and specifically identify filamentous fungi urgently needed for the implementation of appropriately targeted therapy. Early and accurate diagnosis and treatment can decrease morbidity and increase the chances for patient survival (6, 27, 39).

**2**

Furthermore, identification of filamentous fungi to at least the species level would be epidemiologically useful (24, 31, 43, 47).

PCR-based methods of detection, which show promise as rapid, sensitive means to diagnose infections, have been used in the identification of DNA from Candida species (13, 15, 30) and some other fungi, particularly Aspergillus species (31, 33, 45). However, most of these tests are only genus-specific (28, 38) or are directed to detect only single-copy genes (4, 35). Others have designed probes to detect multi-copy genes so as to increase test sensitivity (31, 33) but in doing so have lost test specificity because they have used highly conserved genes, which detect one or a few species but which are also plagued with cross-reactivities to human, fungal or even viral DNA (25, 31, 33).

Therefore, it is an object of the invention to provide improved materials and methods for detecting and differentiating Aspergillus and other filamentous fungal species in the clinical and laboratory settings.

**SUMMARY OF THE INVENTION**

The present invention relates to nucleic acids for detecting Aspergillus, Fusarium, Mucor, Penicillium, Rhizopus, Rhizomucor, Absidia, Cunninghamella, *Pseudallescheria* (*Scedosporium*), and Sporothrix species. Unique internal transcribed spacer 2 coding regions permit the development of probes specific for five different Aspergillus species, *A. flavus*, *A. fumigatus*, *A. niger*, *A. terreus*, and *A. nidulans*. The invention thereby provides methods for the species-specific detection and diagnosis of Aspergillus infection in a subject. In addition, species probes have been developed for three Fusarium, four Mucor, two Penicillium, five Rhizopus and one Rhizomucor species, as well as probes for *Absidia corymbifera*, *Cunninghamella elegans*, *Pseudallescheria boydii* (*Scedosporium apiospermum*), and *Sporothrix schenckii*. Generic probes for Aspergillus, Fusarium, and Mucor species have also been developed.

These and other objects, features and advantages of the present invention will become apparent after a review of the following detailed description of the disclosed embodiments and the appended claims.

**DETAILED DESCRIPTION OF THE INVENTION**

This invention provides a simple, rapid, and useful method for differentiating filamentous fungal species from each other and from other medically important fungi. This invention enables a rapid, simple and useful method to isolate fungal DNA from host samples, and to apply the species- and genus-specific probes for the diagnosis of a disease. Ultimately, these probes can be used for in situ hybridization or in situ PCR diagnostics so that the morphology of host tissue, and microorganisms, remain intact.

The invention provides nucleic acids containing regions of specificity for five Aspergillus, three Fusarium, four Mucor, two Penicillium, five Rhizopus and one Rhizomucor species as well as probes for *Absidia corymbifera*, *Cunninghamella elegans*, *Pseudallescheria boydii* (*Scedosporium apiospermum*), and *Sporothrix schenckii*. These nucleic acids are from the internal transcribed spacer 2 ("ITS2") region of ribosomal deoxyribonucleic acid (rDNA) of the genome of the aforementioned filamentous fungi. The ITS2 region is located between the 5.8S rDNA region and the 28S rDNA region.

In particular, the invention provides nucleic acids from *Aspergillus flavus* (SEQ ID NO:1), *Aspergillus fumigatus*

(SEQ ID NO:2), *Aspergillus niger* (SEQ ID NO:3), *Aspergillus terreus* (SEQ ID NO:4), *Aspergillus nidulans* (SEQ ID NO:5), *Fusarium solani* (SEQ ID NO:6), *Fusarium moniliforme* (SEQ ID NO:7), *Mucor rouxii* (SEQ ID NO:8), *Mucor racemosus* (SEQ ID NO:9), *Mucor plumbeus* (SEQ ID NO:10), *Mucor indicus* (SEQ ID NO:11), *Mucor circinelloides f. circinelloides* (SEQ ID NO:12), *Rhizopus oryzae* (SEQ ID NO:13 and NO:14), *Rhizopus microsporitis* (SEQ ID NO:15 and 16), *Rhizopus circinans* (SEQ ID NO:17 and 18). *Rhizopus stolonifer* (SEQ ID NO:19), *Rhizomucor pusillus* (SEQ ID NO:20), *Absidia corymbifera* (SEQ ID NO:21 and 22), *Cunninghamella elegans* (SEQ ID NO:23), *Pseudallescheria boydii* (teleomorph of *Scedosporium apiospermum*) (SEQ ID NO:24, 25, 26, and 27), *Penicillium notatum* (SEQ ID NO:28), and *Sporothrix schenckii* (SEQ ID NO:29). These sequences can be used to identify and distinguish the respective species of *Aspergillus*, *Fusarium*, *Mucor*, *Rhizopus*, and *Penicillium*, and identify and distinguish these species from each other and from *Absidia corymbifera*, *Cunninghamella elegans*, *Pseudallescheria boydii* (*Scedosporium apiospermum*), and *Sporothrix schenckii*.

Furthermore, the invention provides isolated nucleic acid probes derived from GenBank nucleic acid sequences (for *Penicillium marneffei* and *Fusarium oxysporum* only) or from the above nucleic acid sequences which may be used as species-specific identifiers of *Aspergillus flavus* (SEQ ID NO:30 and 31), *Aspergillus fumigatus* (SEQ ID NO:32), *Aspergillus niger* (SEQ ID NO:33), *Aspergillus terreus* (SEQ ID NO:34), *Aspergillus nidulans* (SEQ ID NO:35), *Mucor rouxii* (SEQ ID NO:36), *Mucor plumbeus* (SEQ ID NO:37), *Mucor indicus* (SEQ ID NO:38), *Mucor circinelloides f. circinelloides* (SEQ ID NO:39), *Mucor racemosus* (SEQ ID NO:40), *Rhizopus oryzae* (SEQ ID NO:41), *Rhizopus circinans* (SEQ ID NO:42), *Rhizomucor pusillus* (SEQ ID NO:43), *Rhizopus stolonifer* (SEQ ID NO:44), *Pseudallescheria boydii* (*Scedosporium apiospermum*) (SEQ ID NO:45), *Penicillium notatum* (SEQ ID NO:46), *Penicillium marneffei* (SEQ ID NO:47 and 48), *Fusarium moniliforme* (SEQ ID NO:49), *Fusarium oxysporum* (SEQ ID NO:50), *Fusarium solani* (SEQ ID NO:51), *Cunninghamella elegans* (SEQ ID NO:52, 53, and 54), *Absidia corymbifera* (SEQ ID NO:55), *Sporothrix schenckii* (SEQ ID NO:56), and *Rhizopus microsporus* (SEQ ID NO:57). Such probes can be used to selectively hybridize with samples containing nucleic acids from species of *Aspergillus*, *Fusarium*, *Mucor*, *Rhizopus* (or *Rhizomucor*), *Penicillium*, or from *Absidia corymbifera*, *Cunninghamella elegans*, *Pseudallescheria boydii* (*Scedosporium apiospermum*), and *Sporothrix schenckii*. These fungi can be detected after polymerase chain reaction or ligase chain reaction amplification of fungal DNA and specific probing of amplified DNA with DNA probes labeled with digoxigenin, reacted with anti-digoxigenin antibodies labeled with horseradish peroxidase and a colorimetric substrate, for example. Additional probes can routinely be derived from the sequences given in SEQ ID NOS: 1–29, which are specific for the respective species. Therefore, the probes shown in SEQ ID NOS:30–57 are only provided as examples of the species-specific probes that can be derived from SEQ ID NOS: 1–29.

Generic probes for *Aspergillus* (SEQ ID NO:58), *Fusarium*, (SEQ ID NO:59) and *Mucor* (SEQ ID NO:60) species have also been developed to identify all members of their respective species which are listed above as well as an all-fungus biotinylated probe (SEQ ID NO:61) to capture all species-specific and generic probes listed above for their detection.

By “isolated” is meant nucleic acid free from at least some of the components with which it naturally occurs. By “selective” or “selectively” is meant a sequence which does not hybridize with other nucleic acids to prevent adequate determination of an *Aspergillus*, *Fusarium*, *Mucor*, *Penicillium*, *Rhizopus* or *Rhizomucor* genus or species or of *Absidia corymbifera*, *Cunninghamella elegans*, *Pseudallescheria boydii* (*Scedosporium apiospermum*), or *Sporothrix schenckii* species.

The hybridizing nucleic acid should have at least 70% complementarity with the segment of the nucleic acid to which it hybridizes. As used herein to describe nucleic acids, the term “selectively hybridizes” excludes the occasional randomly hybridizing nucleic acids and thus has the same meaning as “specifically hybridizing”. The selectively hybridizing nucleic acids of the invention can have at least 70%, 80%, 85%, 90%, 95%, 97%, 98%, and 99% complementarity with the segment of the sequence to which it hybridizes.

The invention contemplates sequences, probes and primers which selectively hybridize to the complementary, or opposite, strand of DNA as those specifically provided herein. Specific hybridization with nucleic acid can occur with minor modifications or substitutions in the nucleic acid, so long as functional species-specific or genus-specific hybridization capability is maintained. By “probe” is meant nucleic acid sequences that can be used as probes or primers for selective hybridization with complementary nucleic acid sequences for their detection or amplification, which probes can vary in length from about 5 to 100 nucleotides, or preferably from about 10 to 50 nucleotides, or most preferably about 18 nucleotides. The invention provides isolated nucleic acids that selectively hybridize with the species-specific nucleic acids under stringent conditions and should have at least 5 nucleotides complementary to the sequence of interest. See generally, Maniatis (26).

If used as primers, the invention provides compositions including at least two nucleic acids which hybridize with different regions so as to amplify a desired region. Depending on the length of the probe or primer, target region can range between 70% complementary bases and full complementarity and still hybridize under stringent conditions. For example, for the purpose of diagnosing the presence of the *Aspergillus*, the degree of complementarity between the hybridizing nucleic acid (probe or primer) and the sequence to which it hybridizes (e.g., *Aspergillus* DNA from a sample) is at least enough to distinguish hybridization with a nucleic acid from other yeasts and filamentous fungi. The invention provides examples of nucleic acids unique to each filamentous fungus in the listed sequences so that the degree of complementarity required to distinguish selectively hybridizing from nonselectively hybridizing nucleic acids under stringent conditions can be clearly determined for each nucleic acid.

Alternatively, the nucleic acid probes can be designed to have homology with nucleotide sequences present in more than one species of the fungi listed above. Such a nucleic acid probe can be used to selectively identify a group of species such as the generic probes listed for *Aspergillus* (SEQ ID NO:58), *Fusarium* (SEQ ID NO:59), and *Mucor* (SEQ ID NO:60) as well as all fungi listed (SEQ ID NO:61). Additionally, the invention provides that the nucleic acids can be used to differentiate the filamentous fungi listed in general from other filamentous fungi and yeasts, such as *Candida* species. Such a determination is clinically significant, since therapies for these infections differ.

The invention further provides methods of using the nucleic acids to detect and identify the presence of the

filamentous fungi listed, or particular species thereof. The method involves the steps of obtaining a sample suspected of containing filamentous fungi. The sample may be taken from an individual, such as blood, saliva, lung lavage fluids, vaginal mucosa, tissues, etc., or taken from the environment. The filamentous fungal cells can then be lysed, and the DNA extracted and precipitated. The DNA is preferably amplified using universal primers derived from the internal transcribed spacer regions, 18S, 5.8S and 28S regions of the filamentous fungal rDNA. Examples of such universal primers are shown below as ITS1 (SEQ ID NO:62), ITS3 (SEQ ID NO:63), ITS4 (SEQ ID NO:64). Detection of filamentous fungal DNA is achieved by hybridizing the amplified DNA with a species-specific probe that selectively hybridizes with the DNA. Detection of hybridization is indicative of the presence of the particular genus (for generic probes) or species (for species probes) of filamentous fungus.

Preferably, detection of nucleic acid (e.g. probes or primers) hybridization can be facilitated by the use of detectable moieties. For example, the species-specific or generic probes can be labeled with digoxigenin, and an all-fungus probe, such as described in SEQ ID NO:61, can be labeled with biotin and used in a streptavidin-coated microtiter plate assay. Other detectable moieties include radioactive labeling, enzyme labeling, and fluorescent labeling, for example.

The invention further contemplates a kit containing one or more species-specific probes, which can be used for the detection of particular filamentous fungal species and genera in a sample. Such a kit can also contain the appropriate reagents for hybridizing the probe to the sample and detecting bound probe. The invention may be further demonstrated by the following non-limiting examples.

#### EXAMPLES

In this example, PCR assay employing universal, fungus-specific primers and a simple, rapid EIA-based format for amplicon detection were used.

##### Extraction of Filamentous Fungal DNA

A mechanical disruption method was used to obtain DNA from filamentous fungal species and an enzymatic disruption method described previously (13) was used to obtain DNA from yeasts. Filamentous fungi were grown for 4 to 5 days on Sabouraud dextrose agar slants (BBL, division of Becton Dickinson, Cockeysville, Md.) at 35° C. Two slants were then washed by vigorously pipeting 5 mls of 0.01 M potassium phosphate buffered saline (PBS) onto the surface of each slant and the washes were transferred to 500 ml Erlenmeyer flasks containing 250 ml of Sabouraud dextrose broth (BBL). Flasks were then incubated for 4 to 5 days on a rotary shaker (140 rpm) at ambient temperature. Growth was then harvested by vacuum filtration through a sterile Whatman #1 filter paper which had been placed into a sterile Buchner funnel attached to a 2 L side-arm flask. The resultant cellular mat was washed on the filtration apparatus three times with sterile distilled water, removed from the filter paper by gentle scraping with a rubber policeman, and placed into a sterile Petri plate which was then sealed with parafilm and frozen at -20° C. until used.

Just prior to use, a portion of the frozen cellular mat, equal in size to a quarter, was removed and placed into a cold mortar (6" diameter). Liquid nitrogen was added to cover the mat which was then ground into a powder with a pestle. Additional liquid nitrogen was added as needed to keep the mat frozen during grinding.

DNA was then purified using proteinase K and RNase treatment, multiple phenol extractions, and ethanol precipitation by conventional means (26).

##### PCR amplification

The fungus-specific, universal primer pair ITS3 (5'-GCA TCG ATG AAG AAC GCA GC-3') (SEQ ID NO:63) and ITS4 (5'-TCC TCC GCT TAT TGA TAT GC-3') (SEQ ID NO:64) was used to amplify a portion of the 5.8S rDNA region, the entire ITS2 region, and a portion of the 28S rDNA region for each species as previously described (13, 34). DNA sequencing used this primer pair and also the fungus-specific, universal primer pair ITS1 (5'-TCC GTA GGT GAA CCT GCG G-3') (SEQ ID NO: 62) and ITS4 to amplify a portion of the 18S rDNA region, the entire 5.8S region, the entire ITS1 and ITS2 regions, and a portion of the 28S rDNA region.

A DNA reagent kit (TaKaRa Biomedicals, Shiga, Japan) was used for PCR amplification of genomic DNA. PCR was performed using 2 µl of test sample in a total PCR reaction volume of 100 µl consisting of 10 µl of 10×Ex Tag buffer, 2.5 mM each of dATP, dGTP, dCTP, and dTTP, in 8 µl 0.2 µM of each primer, and 0.5 U of TaKaRa Ex Tag DNA polymerase. Thirty cycles of amplification were performed in a Perkin-Elmer 9600 thermal cycler (Emeryville, Calif.) after initial denaturation of DNA at 95° C. for 5 minutes. Each cycle consisted of a denaturation step at 95° C. for 30 seconds, an annealing step at 58° C. for 30 seconds, and an extension step at 72° C. for 1 minute. A final extension at 72° C. for 5 minutes followed the last cycle. After amplification, samples were stored at -20° C. until used.

TABLE 1

Synthetic Universal Oligonucleotides Used in PCR and Hybridization Analyses		
Primers or Probes	Nucleotide Sequence (5' to 3')	Chemistry and Location
ITS3	GCA TCG ATG AAG AAC GCA GC (SEQ ID NO:63)	5.8S rDNA universal 5' primer
ITS4	TCC TCC GCT TAT TGA TAT GC (SEQ ID NO:64)	28S rDNA universal 3' primer
ITS1	TCC GTA GGT GAA CCT GCG G (SEQ ID NO:62)	18S rDNA universal 5' primer

##### DNA Sequencing

Primary DNA amplifications were conducted as described above. The aqueous phase of the primary PCR reaction was purified using QIAquick Spin Columns (Quiagen, Chatsworth, Calif.). DNA was eluted from each column with 50 µl of heat-sterilized Tris-EDTA buffer (10 mM Tris, 1 mM EDTA, pH 8.0).

Purified DNA was labeled using a dye terminator cycle sequencing kit (ABI PRISM, Perkin Elmer, Foster City, Calif.). One mix was made for each of the primers so that sequencing could be performed in both the forward and reverse directions. The reaction volume (20 µl) contained 9.5 µl Terminator Premix, 2 µl (1 ng) DNA template, 1 µl primer (3.2 pmol) and 7.5 µl heat-sterilized distilled H<sub>2</sub>O. The mixture was then placed into a pre-heated (96° C.) Perkin Elmer 9600 thermal cycler for 25 cycles of 96° C. for 10 seconds, 50° C. for 5 seconds, 60° C. for 4 minutes. The PCR product was then purified before sequencing using CentriSep spin columns (Princeton Separations, Adelphia, N.J.). DNA was then vacuum dried, resuspended in 6 µl of formamide-EDTA (5 µl deionized formamide plus 1 µl 50 mM EDTA, pH 8.0), and denatured for 2 min at 90° C. prior to sequencing using an automated capillary DNA sequencer (ABI Systems, Model 373, Bethesda, Md.).

The sequencing results were as follows:

*Aspergillus flavus* 5.8S ribosomal RNA gene, partial sequence, internal transcribed spacer 2, complete sequence, and 28S ribosomal RNA gene, partial sequence.

(SEQ ID NO:1)  
 GCTGCCCATC AAGCACGGC TTGTGTGTTG GGTCGTCGTC  
 CCCTCTCCGG GGGGGACGGG CCCCAAAGGC AGCGGCAGCA  
 CCGCGTCCGA TCCTCGAGCG TATGGGGCTT TGTCACCCGC  
 TCTGTAGGCC CGGCCGGCGC TTGCCGAACG CAAATCAATC  
 TTTTCCAGG TTGACCTCGG ATCAGGTAGG GATAACCGCT  
 GAACTTCAA

*Aspergillus fumigatus* 5.8S ribosomal RNA gene, partial sequence, internal transcribed spacer 2, complete sequence, and 28S ribosomal RNA gene, partial sequence.

(SEQ ID NO:2)  
 AAACTTCAA CAATGGATCT CTTGGTTCCG GCATCGATGA  
 AGAACGCAGC GAAATGCGAT AACTAATGTG AATTGCAGAA  
 TTCAGTGAAT CATCGAGTCT TTGAACGCAC ATTGCGCCCC  
 CTGGTATTCC GGGGGCATG CCTGTCCGAG CGTCATTGCT  
 GCCCATCAAAG CACGGCTTGT GTGTTGGGCC CCCGTCCCCC  
 TCTCCCAGGG GACGGGCCCG AAAGGCAGCG GCGGCACCGC  
 GTCCGGTCCT CGAGCGTATG GGGCTTGCA CCTGCTCTGT  
 AGGCCCGGCC GGCGCCAGCC GACACCAAC TTTATTTTC  
 TAAGGTTGAC CTGGATCAG GTAGGGATAC CCGCTGAAC TAAA

*Aspergillus niger* 5.8S ribosomal RNA gene, partial sequence, internal transcribed spacer 2, complete sequence, and 28S ribosomal RNA gene, partial sequence.

(SEQ ID NO:3)  
 AAACTTCAA CAATGGATCT CTTGGTTCCG GCATCGATGA  
 AGAACGCAGC GAAATGCGAT AACTAATGTG AATTGCAGAA  
 TTCAGTGAAT CATCGAGTCT TTGAACGCAC ATTGCGCCCC  
 CTGGTATTCC GGGGGCATG CCTGTCCGAG CGTCATTGCT  
 GCCCTCAAGC ACGGCTTGTG TGTTGGGTG CGTCATCGCT  
 CTCCCGGGGG ACGGGCCGA AAGGCAGCGG CGGCACCGCG  
 TCCGATCCTC GAGCGTATGG GGCTTGCA CCTGCTCTGT  
 AGGCCCGGCC GGCGCCTGCC GACGTTATCC AACCAATTTT  
 TTCCAGGTTG ACCTCGGATC AGGTAGGGAT ACCCGCTGAA CTTAA

*Aspergillus terreus* 5.8S ribosomal RNA gene, partial sequence, internal transcribed spacer 2, complete sequence, and 28S ribosomal RNA gene, partial sequence.

(SEQ ID NO:4)  
 AAACTTCAA CAATGGATCT CTTGGTTCCG GCATCGATGA  
 AGAACGCAGC GAAATGCGAT AACTAATGTG AATTGCAGAA  
 TTCAGTGAAT CATCGAGTCT TTGAACGCAC ATTGCGCCCC

-continued

CTGGTATTCC GGGGGGGCAT GCCTGTCCGA GCGTCATTGC  
 5 TGCCCTCAAG CCCGGCTTGT GTGTTGGGCC CTCGTCCCCC  
 GGCTCCCGGG GGACGGGCC GAAAGGCAGC GGCGGCACCG  
 CGTCCGGTCC TCGAGCGTAT GGGGCTTCGT CTTCCGCTCC  
 10 GTAGGCCCGG CGGGCGCCCG CGAACGCAT TTATTTGCAA  
 CTTGTTTTTT TTTCCAGGTT GACCTCGGAT CAGGT

*Aspergillus nidulans* 5.8S ribosomal RNA gene, partial sequence, internal transcribed spacer 2, complete sequence, and 28S ribosomal RNA gene, partial sequence.

(SEQ ID NO:5)  
 20 AAACTTCAA CAATGGATCT CTTGGTTCCG GCATCGATGA  
 AGAACGCAGC GAACTGCGAT AAGTAATGTG AATTGCAGAA  
 TTCAGTGAAT CATCGAGTCT TTGAACGCAC ATTGCGCCCC  
 25 CTGGCATTCC GGGGGCATG CCTGTCCGAG CGTCATTGCT  
 GCCCTCAAGC CGGGCTTGTG TGTTGGGTG CGTCCCGCC  
 CCCCGGGGG ACGGGCCGAA AGGCAGCGC GGCACCGGTC  
 CGGTCCCTGA GCGTATGGGG CTTGGTCACC CGCTCGATTA  
 30 GGGCCGGCCG GGCGCCAGCC GGCGTCTCCA ACCTTATCTT  
 TCTCAGGTTG ACCTCGGATC AGGTAGGGAT ACCCGCTGAA CTTAA

35 *Fusarium solani* (strain ATCC62877) internal transcribed spacer 2 and adjacent regions.

(SEQ ID NO:6)  
 40 GAAAATGCGA TAAGTAATGT GAATTGCAGA ATTCAAGTGA  
 TCATCGAATC TTTGAACGCA CATTGCGCCC GCCAGTATT  
 TGGCGGGCAT GCCTGTTCGA GCGTCATTAC AACCTCAGG  
 45 CCCCCGGGCC TGGCGTTGGG GATCGCGGA AGCCCCCTGC  
 GGGCACAAACG CGTCCCCCA AATACAGTGG CGGTCCCGCC  
 GCAGCTTCCA TTGCGTAGTA GCTAACACCT CGCAACTGGA  
 GAGCGGCGCG GCCACGCCGT AAAACACCCA ACTTCTGAAT  
 50 GTTGACCTCG AATCAGGTAG GAATACCCGC TGAACCTAA

55 *Fusarium moniliforme* (strain ATCC38519) internal transcribed spacer 2 and adjacent regions.

(SEQ ID NO:7)  
 60 AAATGCGATA AGTAATGTGA ATTGCAAAAT TCAGTGAATC  
 ATCGAATCTT TGAACGCACA TTGCGCCCGC CAGTATTCTG  
 GCGGGCATGC CTGTTGAGC GTCATTCAA CCCTCAAGCC  
 CCCGGGTTTG GTGTTGGGG A CGGGCAAGCC CTTGCGGCAA  
 65 GCCGGCCCG AAATCTAGTG GCGGTCTCGC TGCAGCTTCC  
 ATTGCGTAGT AGTAAAACCC TCGCAACTGG TACGCGGCGC

**-continued**

GGCCAAGCCG TTAAACCCCC AACTCTGAA TGTTGACCTC  
GGATCAGGTA GGAATACCCG CTGAACTTAA 5 (SEQ ID NO:11)  
AAAGTGCAT AACTAGTGTG AATTGCATAT TCAGTGAATC  
ATCGAGTCTT TGAACGCAAC TTGCGCTCAT TGGTATTCCA  
ATGAGCACGC CTGTTTCAGT ATCAAAACAA ACCCTCTATC  
CAGCATTTCAGT TTGAATAGGA ATACTGAGAG TCTCTTGATC  
TATTCTGATC TCGAACCTCT TGAAATGTAC AAAGGCCTGA  
TCTTGTAA ATGCCTGAAC TTTTTTTAA TATAAAGAGA 10 (SEQ ID NO:8)  
AGCTCTTGC GTAAACTGTG CTGGGGCCTC CCAAATAATA  
CTCTTTAA ATTTGATCTG AAATCAGGCG GGATTACCCG  
CTGAACTTAA 15 (SEQ ID NO:9)  
AAAGTGCAT AACTAGTGTG AATTGCATAT TCAGTGAATC  
ATCGAGTCTT TGAACGCAAC TTGCGCTCAT TGGTATTCCA  
ATGAGCACGC CTGTTTCAGT ATCAAAACAA ACCCTCTATC  
CAACTTTGTTGT TGTATAGGAT TATTGGGGC CTCTCGATCT  
GTATAGATCT TGAAATCCCT GAAATTACT AAGGCCTGAA  
CTTGTAA TGCCTGAAC TTTTTTAAT ATAAAGGAAA 20  
GCTCTTGTAA TTGACTTTGA TGGGGCCTCC CAAATAAATC  
TCTTTAAAT TTGATCTGAA ATCAGGCGGG ATTACCCGCT  
GAACCTAA 25  
*Mucor rouxii* (strain ATCC24905) internal transcribed spacer 2 and adjacent regions.  
*Mucor circinelloides f. circinelloides* (strain ATCC1209B) internal transcribed spacer 2 and adjacent regions. 20  
AAAGTGCAT AACTAGTGTG AATTGCATAT TCAGTGAATC  
ATCGAGTCTT TGAACGCAAC TTGCGCTCAT TGGTATTCCA  
ATGAGCACGC CTGTTTCAGT ATCAAAACAA ACCCTCTATC  
CAACATTTT GTTGAATAGG ATGACTGAGA GTCTCTTGAT  
CTATTCTGAT CTCGAAGCTC TTGAAATGTA CAAAGGCCTG  
ATCTTGTAA AATGCCTGAA CTTTTTTAA ATATAAAGAG 30  
AAGCTCTTGC GGTAAACTGT GCTGGGGCCT CCCAAATAAC  
ACATCTTAA ATTTGATCTG AAATCAGGT GGGACTACCC  
GCTGAACTT AA (SEQ ID NO:12)  
*Rhizopus oryzae* (strain ATCC34965) internal transcribed spacer 2 and adjacent regions. 35  
AAAGTGCAT AACTAGTGTG AATTGCATAT TCAGTGAATC  
ATCGAGTCTT TGAACGCAAC TTGCGCTCAT TGGTATTCCA  
ATGAGCACGC CTGTTTCAGT ATCAAAACAA ACCCTCTATC  
CAACTTTGTTGT TGTATAGGAT TATTGGGGC CTCTCGATCT  
GTATAGATCT TGAAACCCTT GAAATTACT AAGGCCTGAA  
CTTGTAA GCCTGAACCTT TTTTTAATA TAAAGGAAAG 40  
CTCTTGTAA TGACTTTGAT GGGGCCTCCC AAATAAATCT  
TTTTAAATT TGATCTGAAA TCAGGTGGGA TTACCCGCTG  
GAACCTAA 45  
*Mucor plumbeus* (strain ATCC4740) internal transcribed spacer 2 and adjacent regions.  
*Rhizopus oryzae* (strain ATCC11886) internal transcribed spacer 2 and adjacent regions. 50  
AAAGTGCAT AACTAGTGTG AATTGCATAT TCAGTGAATC  
ATCGAGTCTT TGAACGCAAC TTGCGCTCAT TGGTATTCCA  
ATGAGCACGC CTGTTTCAGT ATCAAAACAA ACCCTCTATC  
CAACTTTGTTGT TGTATAGGAT TATTGGGGC CTCTCGATCT  
GTATAGATCT TGAAACCCTT GAAATTACT AAGGCCTGAA  
CTTGTAA GCCTGAACCTT TTTTTAATA TAAAGGAAAG 55  
CTCTTGTAA TGACTTTGAT GGGGCCTCCC AAATAAATCT  
TTTTAAATT TGATCTGAAA TCAGGTGGGA TTACCCGCTG  
GAACCTAA  
AAAGTGCAT AACTAGTGTG AATTGCATAT TCAGTGAATC  
ATCGAGTCTT TGAACGCAAC TTGCGCTCAT TGGTATTCCA  
ATGAGCACGC CTGTTTCAGT ATCAAAACAA ACCCTCTATC  
CAACTTTGTTGT TGTATAGGAT TATTGGGGC CTCTCGATCT  
GTATAGATCT TGAAACCCTT GAAATTACT AAGGCCTGAA  
CTTGTAA GCCTGAACCTT TTTTTAATA TAAAGGAAAG 60  
CTCTTGTAA TGACTTTGAT GGGGCCTCCC AAATAAATCT  
TTTTAAATT TGATCTGAAA TCAGGTGGGA TTACCCGCTG  
GAACCTAA (SEQ ID NO:13)  
AAAGTGCAT AACTAGTGTG AATTGCATAT TCAGTGAATC  
ATCGAGTCTT TGAACGCAAC TTGCGCTCAT TGGTATTCCA  
ATGAGCACGC CTGTTTCAGT ATCAAAACAA ACCCTCTATC  
CAACATTTT GTTGAATAGG ATGACTGAGA GTCTCTTGAT  
CTATTCTGAT CTCGAAGCTC TTGAAATGTA CAAAGGCCTG  
ATCTTGTAA AATGCCTGAA CTTTTTTAA ATATAAAGAG 65  
AAGCTCTTGC GGTAAACTGT GCTGGGGCCT CCCAAATAAC  
ACATCTTAA ATTTGATCTG AAATCAGGT GGGACTACCC  
GCTGAACTT AA (SEQ ID NO:14)  
*Rhizopus oryzae* (strain ATCC11886) internal transcribed spacer 2 and adjacent regions.  
AAAGTGCAT AACTAGTGTG AATTGCATAT TCAGTGAATC  
ATCGAGTCTT TGAACGCAAC TTGCGCTCAT TGGTATTCCA  
ATGAGCACGC CTGTTTCAGT ATCAAAACAA ACCCTCTATC  
CAACATTTT GTTGAATAGG ATGACTGAGA GTCTCTTGAT  
CTATTCTGAT CTCGAAGCTC TTGAAATGTA CAAAGGCCTG  
ATCTTGTAA AATGCCTGAA CTTTTTTAA ATATAAAGAG  
AAGCTCTTGC GGTAAACTGT GCTGGGGCCT CCCAAATAAC  
ACATCTTAA ATTTGATCTG AAATCAGGT GGGACTACCC  
GCTGAACTT AA (SEQ ID NO:15)  
*Mucor indicus* (strain ATCC4857) internal transcribed spacer 2 and adjacent regions.

**-continued**

- TTGCTAGGCA GGAATATTAC GCTGGTCTCA GGATCTTTT  
 CTTGGTTCG CCCAGGAAGT AAAGTACAAG AGTATAATCC  
 AGCAACTTTC AAACTATGAT CTGAAGTCAG GTGGGATTAC  
 CCGCTGAAC TAA (SEQ ID NO:14)
- Rhizopus microsporus* (strain ATCC14056) internal transcribed spacer 2 and adjacent regions.
- AAAGTGCAT AACTAGTGTG AATTGCATAT TCGTGAATCA  
 TCGAGTCTTT GAACGCAGCT TGCACTCTAT GGATCTTCTA  
 TAGAGTACGC TTGCTTCAGT ATCATAACCA ACCCACACAT  
 AAAATTATT TTATGTGGTG ATGGACAAGC TCGGTTAAAT  
 TTAATTATTA TACCGATTGT CTAAAATACA GCCTCTTGT  
 AATTTCATT AAATTACGAA CTACCTAGCC ATCGTGCTTT  
 TTTGGTCAA CCAAAAAACA TATAATCTAG GGGTTCTGCT  
 AGCCAGCAGA TATTTAACG ATCTTTAATC ATGATCTGAA  
 GTCAAGTGGG ACTACCCGCT GAACTTAA (SEQ ID NO:15)
- Rhizopus microsporus* (strain ATCC12276) internal transcribed spacer 2 and adjacent regions.
- AAAGTGCAT AACTAGTGTG AATTGCATAT TCGTGAATCA  
 TCGAGTCTTT GAACGCAGCT TGCACTCTAT GGATCTTCTA  
 TAGAGTACGC TTGCTTCAGT ATCATAACCA ACCCACACAT  
 AAAATTATT TTATGTGGTG ATGGACAAGC TCGGTTAAAT  
 TTAATTATTA TACCGATTGT CTAAAATACA GCCTCTTGT  
 AATTTCATT AAATTACGAA CTACCTAGCC ATCGTGCTTT  
 TTTGGTCAA CCAAAAAACA TATAATCTAG GGGTTCTGCT  
 AGCCAGCAA TATTTAACG ATCTTTAACG TATGATCTGA  
 AGTCAAGTGG GACTACCCGCT TGAACCTAA (SEQ ID NO:16)
- Rhizopus circinans* (strain ATCC34106) internal transcribed spacer 2 and adjacent regions.
- AAATTGCGAT AACTAGTGTG AATTGCATTT TCAGTGAATC  
 ATCGAGTCTT TGAACGCAT CTTGCGCTCT TGGGATTCTT  
 CCCTAGAGCA CACTGCTTC AGTATCATAA CAAAACCCTC  
 ACCTAATATT TTTTTTTTTT AAAAAAAA TATTAGAGTG  
 GTATTGGGT CTCTTGGTA ATTCTTGTA ATTATAAAAG  
 TACCCCTAAA TGTCATAAAC AGGTTAGCTT TAGCTTGCCT  
 TTAAAGATCT TCTTAGGGTA TCATTACTTT TCGTAAATCT  
 TTAATAGGCC TGTCACATAA TTCTACCCCTT AAATTCTTA  
 AACCTTGATC TGAAGTCAAG TGGGAGTACC CGCTGAACCTT AA  
 (SEQ ID NO:17)
- Rhizopus circinans* (strain ATCC34101) internal transcribed spacer 2 and adjacent regions.
- AAATTGCGAT AACTAGTGTG AATTGCATTT TCAGTGAATC  
 ATCGAGTCTT TGAACGCATC TTGCGCTCTT GGGATTCTTC  
 CCTAGAGCAC ACTTGCTTCA GTATCATAAC AAAACCTCA  
 CCTAATATT TTGTTTAAAAA AAAAAAAATA TTAGAGTGGT  
 ATTGGGTCT CTTGGTAAT TCTTGTAAT TATAAAAGTA  
 CCCTTAAATG TCATAAACAG GTTAGCTTA GCTTGCCTTT  
 AAAGATCTTC TTAGGGTATC ATTACTTTTC GTAAATCTTT  
 AATAGGCCTG TCACATAATT CTACCCTTAA ATTTCTTAA  
 CCTTGATCTG AAGTCAAGTG GGAGTACCCG CTGAACCTAA  
 (SEQ ID NO:18)
- 20 *Rhizous stolonifer* (strains ATCC14037 and 6227A) internal transcribed spacer 2 and adjacent regions.
- AAAGTGCAT AACTAGTGTG AATTGCATAT TCGTGAATC  
 ATCGAGTCTT TGAACGCAAC TTGCACTCTA TGGTTTCCG  
 TAAAGTACGC TTGCTTCAGT ATCATAAAAGA CCCCACCTG  
 ATTATTATT TTGTTTAAATAATT TTGGAGATAA  
 TAAAAATGAG GCTCTTCTT TTCTTTTTTT TTTTTTAA  
 AAAAAGGGGG GGAAAGGGTC TTTTAAATG GGCAAATTCT  
 GGGTTTTTA CTAAACCTGA ACTCCCCCA AAAATTCAA  
 AAAAAAAA TGGTTTTAC CAAATTTTT TTTTTTTCT  
 CCTTTTGTG TAGTTAATAC TCTATTAAAT TTATTTACTT  
 GGTATTATAA CGATTATGCA AGAAGGGAGA GAACAAAGAA  
 40 TAATGAAAGA GAGTTTTAA ATAAATTCTT TTTTCATTT  
 TTCAATCAAT GATCTGAAGT CAAGTGGGAT TACCCGCTGA  
 ACTTAA (SEQ ID NO:19)
- 45 *Rhizomucor pusillus* (strain ATCC36606) internal transcribed spacer 2 and adjacent regions.
- 50 AAATTGCGAA AAGTAATGCG ATCTGCAGCC TTTGCGAATC  
 ATCGAATTCT CGAACGCACC TTGCACCCCTT TGGTCATCC  
 ATTGGGTACG TCTAGTTCAAG TATCTTATT AACCCCTAAA  
 55 GGTTTATTCTT TTGATAAAATC TTTGGATTG CGGTGCTGAT  
 GGATTTTCAT CCGTTCAAGC TACCCGAACA ATTTGTATGT  
 TGTTGACCT TGATATTCC TTGAGGGCTT GCATTGGTAT  
 60 CTAATTTTT ACCAGTGTGC TTCGAGATGA TCAAGTATAA  
 AGGTCAATCA ACCACAAATA AATTCAACT ATGGATCTGA  
 ACTTAGATGG GATTACCCGC TGAACCTAA (SEQ ID NO:20)
- 65 *Absidia corymbifera* (strain ATCC46774) internal transcribed spacer 2 and adjacent regions.

AAAGTGCAT AATTATTGCG ACTTGCATTC ATAGCGAATC  
 ATCGAGTTCT CGAACGCATC TTGCGCCTAG TAGTCAATCT  
 5 ACTAGGCACA GTTGTTCAG TATCTGCAAC TACCAATCAG  
 TTCAACTTGG TTCTTTGAAC CTAAGCGAGC TGGAAATGGG  
 CTTGTGTTGA TGGCATTCAAG TTGCTGTCAT GGCCTTAAAT  
 10 ACATTTAGTC CTAGGCAATT GGCTTAGTC ATTTGCCGGA  
 TGTAGACTCT AGAGTGCCTG AGGAGCAACG ACTTGTTAG  
 TGAGTTCATA ATTCCAAGTC AATCAGTCTC TTCTTGAACT  
 15 AGGTCTTAAT CTTTATGGAC TAGTGAGAGG ATCTAACTTG  
 GGTCTTCTCT TAAAACAAAC TCACATCTAG ATCTGAAATC  
 AACTGAGATC ACCCGCTGAA CTTAA (SEQ ID NO:21)

*Absidia corymbifera* (strain ATCC46773) internal transcribed spacer 2 and adjacent regions.

AAAGTGCAT AATTATTGCG ACTTGCATTC ATAGTGAATC  
 ATCGAGTTCT TGAACGCATC TTGCGCCTAG TAGTCAATCT  
 ACTAGGCACA GTTGTTCAG TATCTGCAATC CACCAATCAA  
 25 CTTAACCTTT TGTGTTGAGT TGGAACTGGG CTTCTAGTTG  
 ATGGCATTAA GTTGCTGTCA TGGCCTTAAA TCAATGTCCT  
 AGGTGTTAGA ACATCTAACCA CCGGATGGAA ACTTTAGAGC  
 GCTTTAAGAG CAGCTTGGTT AGTGAGTTCA ATAATTCAA  
 GCATTAAGTC TTTTAATGAA CTAGCTTTTC TATCTATGGG  
 30 ACACTACTTG GAGAAATCCA AGTAACCTTT AAAACTCCCAT  
 TTAGATCTGA AATCAACTGA GACCACCCGC TGAACCTAA  
 (SEQ ID NO:22)

*Cunninghamella elegans* (strain ATCC42113) internal transcribed spacer 2 and adjacent regions.

AAATCGCGAT ATGTAATGTG ACTGCCTATA GTGAATCATC  
 AAATCTTGA AACGCATCTT GCACCTTATG GTATTCCATA  
 AGGTACGTCT GTTCAGTAC CACTAATAAA TCTCTCTCTA  
 45 TCCTTGATGA TAGAAAAAAA AAAAATAATT TTTACTGGGC  
 CCGGGGAATC CTTTTTTTT TTTAATAAAA AGGACCAATT  
 TTGGCCAAA AAAAAGGGTT GAACTTTTT TACAGATCT  
 TGCATCTAGT AAAAACCTAG TCGGCTTAA TAGATTTTA  
 50 TTTCTATTA AGTTATAGC CATTCTTATA TTTTTTAA  
 TCTTGGCCTG AAATCAGATG GGATACCCGC TGAACCTAA  
 (SEQ ID NO:23)

*Pseudallescheria boydii* (strain ATCC44328) internal transcribed spacer 2 and adjacent regions (teleomorph of *Scedosporium apiospermum*). 65

AAATGCGATA AGTAATGTAA ATTGCAAAT TCAGTGAATC  
 ATCGAATCTT TGAAACGCAC ATTGCCCGGC CAGTAATCT  
 5 GCCGGGCATG CCTGTCCGAG CGTCATTCA ACCCTCGAAC  
 CTCCGTTTC CTTAGGGAAG CCTAGGGTCG GTGTTGGGC  
 GCTACGGCAA GTCTCGCAA CCCCCGTAGG CCCTGAAATA  
 CAGTGGCGGT CCCGCCGCG TTGCCTCTG CGTAGTAAGT  
 10 CTCTTTGCA AGCTCGCATT GGGTCCGGC GGAGGCCTGC  
 CGTCAAACCA CCTAACAACT CCAGATGGTT TGACCTCGGA  
 15 TCAGGTAGGG TTACCCGCTG AACTAA (SEQ ID NO:24)

*Pseudallescheria boydii* (strain ATCC36282) internal transcribed spacer 2 and adjacent regions (teleomorph of *Scedosporium apiospermum*).

GAAATGCGAT AAGTAATGTG AATTGCAGAA TTCAGTGAAT  
 CATCGAATCTT TGAAACGCAC CATTGCCGCC GGCAGTAATC  
 25 TGCCGGGCAT GCCTGTCCGA GCGTCATTTC ACCCTCGAA  
 CCTCCGTTTC CTCAGGGAAG CTCAGGGTCG GTGTTGGGC  
 GCTACGGCAA GTCTCGCAA CCCTCCGTAG GCCCTGAAAT  
 30 ACAGTGGCGG TCCCGCCGCG GTTGCCTCT GCGTAGAAGT  
 CTCTTTGCA AGCTCGCATT GGGTCCGGC GGAGGCCTGC  
 CGTCAAACCA CCTATAACTC CAAATGGTTT GACCTCGGAT  
 35 CAGGTAGGGT TACCCGCTGA ACTTAA (SEQ ID NO:25)

*Scedosporium apiospermum* (strain ATCC64215) internal transcribed spacer 2 and adjacent regions.

GAAATGCGAT AAGTAATGTG AATTGCAGAA TTCAGTGAATC  
 ATCGAATCTT TGAAACGCAC TTGCGCCGG CAGTAATCTG  
 40 CCGGGCATGC CTGTCCGAGC GTCATTCAA CCCTCGAAC  
 TCCGTTTCCT CAGGGAAGCT CAGGGTCGGT GTTGGGGCGC  
 TACGGCGAGT CTTCGCGACC CTCCGTAGGC CCTGAAATAC  
 AGTGGCGGTC CCGCCGCGGT TGCCTCTGC GTAGTAAGTC  
 45 TCTTTGCAA GCTCGCATTG GGTCCCGCG GAGGCCTGCC  
 GTCAAACAC CTATAACTCC AGATGGTTT ACCTCGGATC  
 AGGTAGGTAC CCGCTGAAC TAA (SEQ ID NO:26)

*Scedosporium apiospermum* (strain ATCC46173) internal transcribed spacer 2 and adjacent regions.

AAATGCGATA AGTAATGTGA ATTGCAGAAAT TCAGTGAATC  
 ATCGAATCTT TGAAACGCAC TTGCGCCGG CAGTAATCTG  
 60 CCGGGCATGC CTGTCCGAGC GTCATTCAA CCCTCGAAC  
 TCCGTTTCCT CAGGGAAGCT CAGGGTCGGT GTTGGGGCGC  
 TACGGCGAGT CTTCGCGACC CTCCGTAGGC CCTGAAATAC

-continued

```

AGTGGCGGTC CCGCCGCGGT TGCCTTCTGC GTAGTAAGTC
TCTTTGCAA GCTCGCATTG GGTCCCGGCG GAGGCCTGCC
GTCAAACAC CTATAACTCC AGATGGTTG ACCTCGGATC
AGGTAGGTAC CCGCTGAAC TAA (SEQ ID NO:27)

```

*Penicillium notatum* (strain ATCC10108) internal transcribed spacer 2 and adjacent regions.

```

AAATGCGATA CGTAATGTGA ATTGCAAATT CAGTGAATCA
TCGAGTCTT TGAACGCACA TTGCGCCCCC TGGTATTCCG
GGGGGCATGC CTGTCCGAGC GTCATTGCTG CCCTCAAGCA
CGGCTTGTGT GTTGGGCCCG GTCCTCCGAT CCCGGGGGAC
GGGCCCAGAA GGCAGCGCG GCACCGCGTC CGGTCCCTCGA
GCGTATGGGG CTTTGTCAAC CGCTCTGTAG GCCCGGCCGG
CGCTTGCCGA TCAACCCAAA TTTTATCCA GGTTGACCTC
GGATCAGGTA GGGATACCCG CTGAACCTAA (SEQ ID NO:28)

```

*Sporothrix schenckii* (strain ATCC14284) internal transcribed spacer 2 and adjacent regions.

```

GAAATGCGAT ACTAATGTGA ATTGCAGAAT TCAGCGAAC
ATCGAATCTT TGAACGCACA TTGCGCCCGC CAGCATTCTG
GCGGGCATGC CTGTCCGAGC GTCATTCCC CCCTCACCGC
CCCCGTTGCG CGCTGGTGT GGGGCGCCCT CGGCCTGGCG
GGGGGCCCGC GAAAGCGAGT GGCGGGCCCT GTGGAAGGCT
CCGAGCGCAG TACCGAACGC ATGTTCTCCC CTCGCTCCGG
AGGCCCCCAGA GGCGCCCTGC CGGTGAAAAC GCGCATGACG
CGCAGCTCTT TTTACAAGGT TGACCTCGGA TCAGGTGAGG 2
ATACCCGCTG ACTTAA (SEQ ID NO:29)

```

#### Contamination Precautions

Precautions were taken to avoid possible contamination of PCR samples by following the guidelines of Fujita and Kwok (13, 22). All buffers and distilled water used for PCR assays were autoclaved and fresh PCR reagents were aliquoted prior to use. Physical separation of laboratory areas used to prepare PCR assays and to analyze PCR products, and the use of aerosol-resistant pipette tips, reduced possible cross-contamination of samples by aerosols. Appropriate negative controls were included in each test run, including controls omitting either the primer or the DNA template during PCR assays.

#### Agarose gel Electrophoresis

Gel electrophoresis was conducted in TBE buffer (0.1 M Tris, 0.09 M boric acid, 1 mM EDTA, pH 8.4) at 80 V for 1 to 2 hours using gels composed of 1% (w/vol) agarose (International Technologies, New Haven, Conn.) and 1%

(w/vol) NuSieve agar (FMC Bioproducts, Rockland, Me.). Gels were stained with 0.5 µg of ethidium bromide (EtBr) per ml of distilled H<sub>2</sub>O for 10 minutes followed by three serial washes for 10 minutes each with distilled H<sub>2</sub>O.

#### Microtitration Plate Enzyme Immunoassay for the Detection of PCR Products

Amplicons were detected using species-specific and genus probes labeled with digoxigenin and an all-filamentous fungal probe labeled with biotin in a streptavidin-coated microtiter plate format (13, 34). Ten µl of PCR product was added to each 1.5 ml Eppendorf tube. Single-stranded DNA was then prepared by heating the tubes at 95° C. for 5 minutes and cooling immediately on ice. Two-tenths of a ml of hybridization solution [4×SSC (saline sodium citrate buffer, 0.6 M NaCl, 0.06 M trisodium citrate, pH 7.0) containing 20 mM Hepes, 2 mM EDTA, and 0.15% (vol/vol) Tween 20] supplemented with 50 ng/ml each of the all-Aspergillus biotinylated probe and a species-specific digoxigenin-labeled probe was added to each tube containing denatured PCR product. Tubes were mixed by inversion and placed in a water bath at 37° C. to allow probes to anneal to PCR product DNA. After 1 hour, 100 µl of each sample was added to duplicate wells of a commercially prepared streptavidin-coated microtitration plate (Boehringer Mannheim, Indianapolis, Ind.). The plate was incubated at ambient temperature for 1 hour with shaking, using a microtitration plate shaker (manufactured for Dynatech by CLTI, Middletown, N.Y.). Plates were washed 6 times with 0.01 M potassium phosphate buffered saline, pH 7.2, containing 0.05% Tween 20 (PBST). Each well then received 100 µl of horseradish peroxidase-conjugated, anti-digoxigenin Fab fragment (Boehringer Mannheim) diluted 1:1000 in hybridization buffer. After incubation at ambient temperature for 30 minutes with shaking, the plate was washed 6 times with PBST. One hundred µl of a mixture of one volume of 3, 3', 5, 5'-tetramethyl benzidine peroxidase substrate (Kirkegaard and Perry Laboratories, Inc., Gaithersberg, Md.) and one volume of peroxidase solution (Kirkegaard and Perry Laboratories) was added to each well and the plate was placed at ambient temperature for 10 minutes for color development. The A<sub>650</sub> nm of each well was determined with a microtitration plate reader (UV Max, Molecular Devices, Inc., Menlo Park, Calif.). The absorbance value for the reagent blank, where DNA was absent but replaced with distilled H<sub>2</sub>O, was subtracted from each test sample.

#### Statistical Analysis

The Student's t test was used to determine differences between sample means. Means are expressed as the mean plus or minus the standard error from the mean. Differences were considered significant when P<0.05.

The following probes were used to detect and distinguish each species.

TABLE 2

<u>Probe Sequences</u>		
PROBES	5' to 3' OLIGONUCLEOTIDE SEQUENCE	
<u>Generic Biotin Probe</u>	5' end-labeled biotinylated	
B-58	probe 5.8S region of rDNA GAA TCA TCG A(AG)T CTT TGA ACG	SEQ ID NO 61
Digoxigenin-probe	5' end-labeled digoxigenin probe ITS2 region of rDNA	
<u>Aspergillus species</u>		
<i>A. flavus</i> 22	GCA AAT CAA TCT TTT TCC	SEQ ID NO 30
<i>A. flavus</i> 23	GAA CGC AAA TCA ATC TTT	SEQ ID NO 31
<i>A. fumigatus</i>	CCG ACA CCC ATC TTT ATT	SEQ ID NO 32
<i>A. niger</i>	GAC GTT ATC CAA CCA TTT	SEQ ID NO 33
<i>A. nidulans</i>	GGC GTC TCC AAC CTT ATC	SEQ ID NO 35
<i>A. terreus</i>	GCA TTT ATT TGC AAC TTG	SEQ ID NO 34
<u>Fusarium species</u>		
<i>F. moniliforme</i>	TCT AGT GAC GGT CTC GCT	SEQ ID NO 49
<i>F. oxysporum</i>	CGT TAA TTC GCG TTC CTC	SEQ ID NO 50
<i>F. solani</i>	CTA ACA CCT CGC AAC TGG AGA	SEQ ID NO 51
<u>Mucor species</u>		
<i>M. circinelloides</i>	AAC ATT TTT GTG AAT AGG ATG	SEQ ID NO 39
<i>M. indicus</i>	CGT GGA TTG AGT GCC GAT	SEQ ID NO 38
<i>M. plumbeus</i>	GAA ACC CTT GAA ATT	SEQ ID NO 37
<i>M. rouxii</i>	GAA TAG GAA TAC TGA GAG	SEQ ID NO 36
<i>M. racemosus</i>	GAA ATC CCT GAA ATT	SEQ ID NO 40
<u>Penicillium species</u>		
<i>Penicillium marneffei</i> 1	GGG TTG GTC ACC ACC ATA	SEQ ID NO 47
<i>Penicillium marneffei</i> 2	TGG TCA CCA CCA TAT TTA	SEQ ID NO 48
<i>Penicillium notatum</i>	GAT CAA CCC AAA TTT TTA	SEQ ID NO 46
<u>Rhizopus species</u>		
<i>R. circinans</i>	CTT AGG GTA TCA TTA CTT	SEQ ID NO 42
<i>R. microsporus</i>	CAT ATA ATC TAG GGG TTC	SEQ ID NO 57
<i>R. oryzae</i>	GAG TAT AAT CCA G(CT)A ACT	SEQ ID NO 41
<i>R. stolonifer</i>	CTT GGT ATT ATA ACG ATT	SEQ ID NO 44
<i>Rhizomucor pusillus</i>	TCC TTG AGG GCT TGC ATT	SEQ ID NO 43
<u>Other Genera</u>		
<i>Absidia corymbifera</i>	GTT GCT GTC ATG GCC TTA	SEQ ID NO 55
<i>Cunninghamella elegans</i> 4	TAG TCG GCT TTA ATA GAT	SEQ ID NO 52
<i>Cunninghamella elegans</i> 5	TAT TAA GTT TAT AGC CAT	SEQ ID NO 53
<i>Cunninghamella elegans</i> 6	TAA GTt TAT AGC CAT TCT	SEQ ID NO 54
<i>Pseudallescheria boydii</i>	AAG TCT CTT TTG CAA GCT	SEQ ID NO 45
<i>Sporothrix schoenckii</i>	GAC GCG CAG CTC TTT TTA	SEQ ID NO 56
<u>Genus Probes</u>		
G-ASPERGILLUS	CCT CGA GCG TAT GGG GCT	SEQ ID NO 58
G-FUSARIUM	CCC AAC TTC TGA ATG TTG	SEQ ID NO 59
G-MUCOR	(AC)TG GGG CCT CCC AAA TAA	SEQ ID NO 60

Species-specific probes to the ITS2 region of rDNA for *Aspergillus fumigatus* (SEQ ID NO:32), *A. flavus* (SEQ ID NO:31), *A. niger* (SEQ ID NO:33), *A. terreus* (SEQ ID NO:34), and *A. nidulans* (SEQ ID NO:35) correctly identified each of the respective species ( $P<0.001$ ), and gave no

false-positive reactions with *Rhizopus*, *Mucor*, *Fusarium*, *Penicillium*, or *Candida* species. The *A. flavus* probe also recognized *A. oryzae*, which belongs to the *A. flavus* group. Identification time was reduced from a mean of 5 days by conventional methods to 8 hours.

TABLE 3

Fungus	<u>Aspergillus Probes</u>				
	<i>A. fumigatus</i>	<i>A. nidulans</i>	<i>A. niger</i>	<i>A. terreus</i>	<i>A. flavus</i>
<i>A. fumigatus</i> (n = 6)	2.197 ± 0.187	0.002	0.000	0.001	0.001
<i>A. nidulans</i> (n = 3)	0.001	1.315 ± 0.464	0.002	0.000	0.001
<i>A. niger</i> (n = 5)	0.000	0.000	1.242 ± 0.471	0.001	0.003
<i>A. terreus</i> (n = 4)	0.001	0.000	0.001	1.603 ± 0.378	0.001
<i>A. flavus</i> (n = 6)	0.001	0.001	0.000	0.001	2.043 ± 0.390
<i>A. oryzae</i> (n = 2)	0.001	0.002	0.001	0.001	2.445 ± 0.106
<i>A. parasitica</i> (n = 1)	0.001	0.002	0.002	0.002	0.051
<i>A. clavus</i> (n = 1)	0.005	0.005	0.006	0.005	0.003
<i>C. albicans</i> (n = 1)	0.002	0.001	0.002	0.000	0.000
<i>C. parasilosis</i> (n = 1)	0.001	0.002	0.002	0.002	0.001
<i>C. glabrata</i> (n = 1)	0.001	0.003	0.001	0.001	0.005
<i>C. krusei</i> (n = 1)	0.002	0.002	0.002	0.001	0.001
<i>C. tropicalis</i> (n = 1)	0.002	0.002	0.001	0.000	0.001
<i>F. moniliforme</i> (n = 1)	0.003	0.003	0.001	0.001	0.001
<i>F. solani</i> (n = 1)	0.006	0.002	0.001	0.000	0.001
<i>R. oryzae</i> (n = 1)	0.001	0.001	0.001	0.001	0.001
<i>M. racemosus</i> (n = 1)	0.001	0.002	0.005	0.002	0.000
<i>P. notatum</i> (n = 1)	0.001	0.002	0.002	0.002	0.000
Avg ± SD	0.001 ±	0.001 ±	0.000 ±	0.000 ±	0.002 ±
negative controls	0.002	0.001	0.002	0.002	0.010

Species-specific probes to the ITS2 region of rDNA for *Fusarium oxysporum*, *F. solani*, and *F. moniliforme*, correctly identified each of the respective species (P<0.001), and gave no false-positive reactions with Blastomycetes,

<sup>40</sup> Apophysomyces, Candida, Aspergillus, Mucor, Penecillium, Rhizopus, Rhizomucor, Absidia, Cunninghamella, Pseudallescheria, Sporothrix, or Neosartorya. Empty boxes in Table 4 represent zero probe reactivity.

TABLE 4

Fungus	<u>Fusarium Probes</u>			Generic Fusarium
	<i>F. oxysporum</i>	<i>F. solani</i>	<i>F. moniliforme</i>	
<i>F. oxysporum</i> (n = 3)	1.40 ± 0.13			1.76 ± 0.27
<i>F. solani</i> (n = 5)		1.57 ± 0.07		1.35 ± 0.28
<i>F. moniliforme</i> (n = 2)			1.40 ± 0.91	1.34 ±
Negative control				
<i>A. fumigatus</i>				
<i>A. flavus</i>				
<i>A. niger</i>				
<i>A. nidulans</i>				
<i>A. terreus</i>				
<i>A. parasiticus</i>				
<i>A. clavatus</i>				
<i>P. marneffei</i>		0.01	0.01	0.01
<i>P. notatum</i>	0.01	0.01	0.01	0.01
<i>Rhizopus oryzae</i>		0.03	0.01	0.01
<i>Rhizopus microsporus</i>		0.01	0.01	0.01

TABLE 4-continued

Fungus	<u>Fusarium</u> Probes			Generic Fusarium
	<i>F. oxysporum</i>	<i>F. solani</i>	<i>F. moniliforme</i>	
<i>Rhizopus circinans</i>		0.01	0.01	
<i>Rhizopus stolonifer</i>	0.01	0.01		
<i>Rhizomucor pusillus</i>	0.03	0.02		
<i>M. racemosus</i>				
<i>M. circinelloides</i>				
<i>M. rouxii</i>				
<i>M. plumbeus</i>				
<i>M. indicus</i>				
<i>Absidia corymbifera</i>		0.01	0.01	
<i>Cunninghamella elegans</i>		0.01	0.02	
<i>P. boydii</i>	0.02			
<i>Sporothrix schenckii</i>		0.01	0.01	
<i>C. albicans</i>				
<i>C. tropicalis</i>				
<i>C. krusei</i>				
<i>C. parasilosis</i>				
<i>C. glabrata</i>				
<i>Neosartorya fischeri</i>		0.01		
<i>Blastomyces dermatitidis</i>				
<i>Apophysomyces elegans</i>				
Average of negative controls	0.001 ± 0.002	0.005 ± 0.01	0.004 ± 0.006	

Species-specific probes to various other zygomycetes are presented in Table 5, showing correct identification of each species and no false positives. The exceptions are that the *M. circinelloides* probe hybridized with the *M. rouxii* DNA and the *M. plumbeus* probe hybridized with the *M. racemosus*

DNA. However, the *M. rouxii* probe did not hybridize with *M. circinelloides* DNA, nor did the *M. racemosus* probe hybridize with *M. plumbeus* DNA. Therefore, by a process of elimination, each species can be correctly identified. Empty boxes in Table 5 represent zero probe reactivity.

TABLE 5

FUNGUS	<u>Zygomycetes</u> Probes											
	D-probes RORY	RMIC	RCIR	RSTOL	RPUS	MRACE	MCIR	MRX	MPLUM	MIND	ABS	CUN
<i>R. oryzae</i> (n = 5)	1.50 ± 0.48					0.01						
<i>R. microsporus</i> (n = 5)		0.96 ± 0.61										
<i>R. circinans</i> (n = 3)			1.56 ± 0.19									
<i>R. stolonifer</i> (n = 5)				2.53 ± 0.07				0.01				
<i>Rhizomucor pusillus</i> (n = 2)					1.10 ± 0.68							
<i>M. racemosus</i> (n = 6)				0.01		2.02 ± 0.34			0.29 ± 0.52			
<i>M. circinelloides</i> (n = 3)							1.63 ± 0.37	0.01	0.02			
<i>M. rouxii</i> (n = 1)							1.77	0.76				
<i>M. plumbeus</i> (n = 2)									2.14 ± 0.25			
<i>M. indicus</i> (n = 1)	0.01									1.70 ± 0.04		
<i>Absidia corymbifera</i> (n = 2)					0.01				0.01		1.61 ± 0.08	
<i>Cunninghamella elegans</i> (n = 2)	0.01											2.26 ± 0.03
Negative control												
<i>A. fumigatus</i>									0.01	0.02		
<i>A. flavus</i>					0.01					0.05		
<i>A. niger</i>								0.01				
<i>A. nidulans</i>									0.01	0.01		
<i>A. terreus</i>	0.01											
<i>A. parasiticus</i>					0.01					0.03		

TABLE 5-continued

FUNGUS	<u>Zygomycetes Probes</u>											
	D-probes RORY	RMIC	RCIR	RSTOL	RPUS	MRACE	MCIR	MRX	MPLUM	MIND	ABS	CUN
<i>A. clavatus</i>										0.02		
<i>P. marneffei</i>				0.01								
<i>P. notatum</i>										0.03		
<i>F. oxysporum</i>								0.01				
<i>F. solani</i>												
<i>F. moniliforme</i>	0.01				0.01				0.01		0.01	
<i>P. boydii</i>	0.02											
<i>Sporothrix schenckii</i>												
<i>C. albicans</i>												
<i>C. tropicalis</i>												
<i>C. krusei</i>												
<i>C. parasilosis</i>												
<i>C. glabrata</i>												
<i>Neosartorya fischeri</i>			0.01									
<i>Blastomyces dermatitidis</i>												
<i>Apophysomyces elegans</i>												
Average	0.001 ± .004	0.001 ± 0.02	0.000 ± 0.002	0.000 ± 0.003	0.001 ± 0.003	0.001 ± 0.002	0.001 ± 0.002	0.001 ± 0.003	0.003 ± 0.005	0.005 ± 0.01	0.001 ± 0.001	

25

Species-specific probes to various other fungi are presented in Table 6, showing correct identification of each

species and no false positives. Empty boxes in Table 6 represent zero probe reactivity.

TABLE 6

<u>Pseudallescheria and Sporothrix Probes</u>				
Fungus	<i>P. boydii</i>	<i>P. marneffei</i>	<i>P. notatum</i>	<i>Sporothrix schenckii</i>
<i>P. boydii</i> (n = 4)	1.65 ± 0.48			
<i>P. marneffei</i> (n = 3)	0.01	1.24 ± 0.12		
<i>P. notatum</i> (n = 3)			1.93 ± 0.25	
<i>Sporothrix schenckii</i> (n = 3)	0.01			1.94 ± 0.25
Negative control				
<i>A. fumigatus</i>	0.01			
<i>A. flavus</i>				
<i>A. niger</i>				
<i>A. nidulans</i>				
<i>A. terreus</i>				
<i>A. parasiticus</i>				
<i>A. clavatus</i>			0.11	
<i>F. oxysporum</i>		0.10		
<i>F. solani</i>		0.14		
<i>F. moniliforme</i>		0.08		
<i>R. oryzae</i>	0.01			
<i>R. microsporus</i>	0.01			
<i>R. circinans</i>	0.01			
<i>R. stolonifer</i>	0.01			
<i>Rhizomucor pusillus</i>				
<i>M. racemosus</i>		0.04		
<i>M. circinelloides</i>	0.01	0.09		
<i>M. rouxii</i>	0.01			
<i>M. plumbeus</i>		0.05		
<i>M. indicus</i>				
<i>Absidia corymbifera</i>	0.01			
<i>Cunninghamella bertholletiae</i>	0.01			
<i>C. albicans</i>				
<i>C. tropicalis</i>		0.02		
<i>C. krusei</i>				
<i>C. parasilosis</i>				
<i>C. glabrata</i>				
<i>Neosartorya pseudofischeri</i>			0.03	
<i>Blastomyces dermatitidis</i>	0.01			

TABLE 6-continued

<u>Pseudallescheria and Sporothrix Probes</u>				
Fungus	<i>P. boydii</i>	<i>P. marneffei</i>	<i>P. notatum</i>	<i>Sporothrix schenckii</i>
<i>Apophysomyces elegans</i>	0.01			
Average Negative Controls	0.004 ± 0.002	0.013 ± 0.03	0.002 ± 0.019	0.001 ± 0.002

All of the references mentioned in this Specification are hereby incorporated by reference in their entirety.

## References

1. Ampel, N. M. 1996. Emerging disease issues and fungal pathogens associated with HIV infection. *Emerg. Infect. Dis.* 2:109–116.
2. Andriole, V. T. 1996. Aspergillus infections: problems in diagnosis and treatment. *Infect. Agents and Dis.* 5:47–54.
3. Andriole, V. T. 1993. Infections with Aspergillus species. *Clin. Infect. Dis.* 17 Suppl 2:S481–S486.
4. Bir, N., A. Paliwal, K. Muralidhar, P. Reddy, and P. U. Sarma. 1995. A rapid method for the isolation of genomic DNA from *Aspergillus fumigatus*. *Prep. Biochem.* 25:171–181.
5. Blum, U., M. Windfuhr, C. Buitrago-Tellez, G. Sigmund, E. W. Herbst, and M. Langer. 1994. Invasive pulmonary aspergillosis. MRI, CT, and plain radiographic findings and their contribution for early diagnosis. *Chest* 106:1156–1161.
6. Caillot, D., O. Casasnovas, A. Bernard, J. F. Couaillier, C. Durand, B. Cuisenier, E. Solary, F. Piard, T. Petrella, A. Bonnin, G. Couaillault, M. Dumas, and H. Guy, 1997. Improved management of invasive pulmonary aspergillosis in neutropenic patients using early thoracic computed tomographic scan and surgery. *J. Clin. Oncol.* 15:139–147.
7. Denning, D. W. Therapeutic outcome in invasive aspergillosis. *Clin. Infect. Dis.* 23:608–615.
8. Denning, D. W. Diagnosis and management of invasive aspergillosis. *Curr. Clin. Topics Inf. Dis.* 16:277–299.
9. de Repentigny, L., L. Kaufman, G. T. Cole, D. Kruse, J. P. Latge, and R. C. Matthews. 1994. Immunodiagnosis of invasive fungal infections. *J. Med. Vet. Mycol.* 32 Suppl 1239–252.
10. Dupont, B., D. W. Denning, D. Marriott, A. Sugar, M. A. Viviani, and T. Sirisanthana. 1994. Mycoses in AIDS patients. *J. Med. Vet. Mycol.* 32 Suppl 1:221–239.
11. Fisher, B. D., D. Armstrong, B. Yu, and J. W. M. Gold. 1981. Invasive aspergillosis: progress in early diagnosis and treatment. *Am. J. Med.* 71:571–577.
12. Fridkin, S. K. and W. R. Jarvis. 1996. Epidemiology of nosocomial fungal infections. *Clin. Microbiol. Rev.* 9:499–511.
13. Fujita, S.-I., B. A. Lasker, T. J. Lott, E. Reiss, and C. J. Morrison. 1995. Micro titration plate enzyme immunoassay to detect PCR-amplified DNA from *Candida* species in blood. *J. Clin. Microbiol.* 33:962–967.
14. Gordon, M. A., E. W. Lapa, and J. Kane. 1977. Modified indirect fluorescent antibody test for aspergillosis. *J. Clin. Microbiol.* 6:161–165.
15. Holmes, A. R., R. D. Cannon, M. G. Shepard, and H. F. Jenkinson. 1994. Detection of *Candida albicans* and other yeast in blood by PCR. *J. Clin. Microbiol.* 32:228–231.
16. Hung, C. C., S. C. Chang, P. C. Yang, W. C. Hseigh. 1994. Invasive pulmonary pseudallescheriasis with direct invasion of the thoracic spine in an immunocompromised patient. *Eur. J. Clin. Microbiol. Infect. Dis.* 13:749–751.
17. Kappe, R., and H. P. Seeliger. 1993. Serodiagnosis of deep-seated fungal infections. *Curr. Topics Med. Mycol.* 5:247–280.
18. Kappe, R., A. Schulze-Berge, H. G. Sonntag. 1996. Evaluation of eight antibody tests and one antigen test for the diagnosis of invasive aspergillosis. *Mycoses* 39:13–23.
19. Kaufman and Reiss, Manual of Clinical Microbiology.
20. Kremery, V., Jr., E. Kunova, Z. Jesenska, J. Trupi, S. Spanik, J. Mardiak, M. Studena, and E. Kukuckova. 1996. Invasive mold infections in cancer patients: 5 years' experience with *Aspergillus*, *Mucor*, *Fusarium* and *Acremonium* infections. *Supportive Care in Cancer* 4:39–45.
21. Khoo, S. H., and D. W. Denning. 1994. Invasive aspergillosis in patients with AIDS. *Clin. Infect. Dis.* 19 Suppl 1: S41–S48.
22. Kwok, S., and R. Higuichi. 1989. Avoiding false positives with PCR. *Nature (London)* 339:237–238.
23. Larone, D. H. *Medically Important Fungi: A Guide to Identification.* 3rd ed. ASM Press, Washington, D. C. 1995.
24. Leenders, A., A. van Belkum, S. Janssen, S. de Marie, J. Kluytmans, J. Wielenga, B. Lowenberg, and H. Verbrugh. 1996. Molecular epidemiology of apparent outbreak of invasive aspergillosis in a hematology ward. *J. Clin. Microbiol.* 34:345–351.
25. Makimura, K., S. Y. Murayama, H. Yamaguchi. 1994. Specific detection of *Aspergillus* and *Penicillium* species from respiratory specimens by polymerase chain reaction (PCR). *Jap. J. Med. Sci. Biol.* 47:141–156.
26. Maniatis, T., E. F. Fritsch, and J. Sambrook. 1982. *Molecular cloning: a laboratory manual.* Cold Spring Harbor Laboratory, Cold Spring Harbor, N.Y.
27. Martino, P., and C. Girmenia. 1993. Diagnosis and treatment of invasive fungal infections in cancer patients. *Supportive Care in Cancer.* 1:240–244.
28. Melchers, W. J., P. E. Verweij, P. van den Hurk, A. van Belkum, B. E. De Pauw, J. A. Hoogkamp-Korstanje, and J. F. Meis. 1994. General primer-mediated PCR for detection of *Aspergillus* species. *J. Clin. Microbiol.* 32:1710–1717.
29. Miller, W. T. J., G. J. Sals, I. Frank, W. B. Gefter, M. Aronchick, W. T. Miller. 1994. Pulmonary aspergillosis patients with AIDS. Clinical and radiographic correlations. *Chest* 105:37–44.
30. Miyakawa, Y., T. Mabuchi, and Y. Fukazawa. 1993. New method for detection of *Candida albicans* in human blood by polymerase chain reaction. *J. Clin. Microbiol.* 31:3344–3347.
31. Montone, K. T., and L. A. Litzky. 1995. Rapid method for detection of *Aspergillus* 5S ribosomal RNA using a genus-specific oligonucleotide probe. *J. Clin. Microbiol.* 103:48–51.
32. Rogers, T. R., K. A. Haynes, and R. A. Barnes. 1990. Value of antigen detection in predicting invasive aspergillosis. *Lancet* 336:1210–1213.

33. Sandhu, G. S., B. C. Kline, L. Stockman, and G. D. Roberts. 1995. Molecular probes for diagnosis of fungal infections. *J. Clin. Microbiol.* 33:2913–2919.
34. Shin, J. H., F. S. Nolte, and C. J. Morrison. 1997. Rapid identification of *Candida* species in blood cultures using a clinically useful PCR method. *J. Clin. Microbiol.* in press.
35. Tang, C. M., D. W. Holden, A. Aufauvre-Brown, and J. Cohen. The detection of *Aspergillus* spp. by the polymerase chain reaction and its evaluation in bronchoalveolar lavage fluid. *Amer. Rev. Respir. Dis.* 148:1313–1317.
36. Thompson, B. H., W. Stanford, J. R. Galvin, and Y. Kurlhara. 1995. Varied radiologic appearances of pulmonary aspergillosis. *Radiographics* 15:1273–1284.
37. Tierney, Jr. L. M. Aspergillosis. In *Current Medical Diagnosis and Treatment*. 33rd ed. Norwalk, Conn.: Appleton and Lange, 1994.
38. Verweij, P. E., J. P. Latge, A. J. Rijs, W. J. Melchers, B. E. De Pauw, J. A. Hoogkamp-Korstanje, and J. F. Mels. 1995. Comparison of antigen detection and PCR assay using bronchoalveolar lavage fluid for diagnosing invasive pulmonary aspergillosis in patients receiving treatment for hematological malignancies. *J. Clin. Microbiol.* 33:3150–3153.
39. von Eiff, M., N. Roos, R. Schulten, M. Hesse, M. Zuhisdorf, and J. van de Loo. 1995. Pulmonary aspergillosis: early diagnosis improves survival. *Respiration* 62:341–347.

40. von Eiff, M., N. Roos, W. Fegeler, C. von Eiff, R. Schulten, M. Hesse, M. Zuhisdorf, and J. van de Loo. 1996. Hospital acquired *Candida* and *Aspergillus* pneumonia-diagnostic approaches and clinical findings. *J. Hosp. Infect.* 32:17–28.
41. Walsh, T. J. 1993. Management of immunocompromised patients with evidence of an invasive mycosis. *Hematol. Oncol. Clin. N. Amer.* 7:1003–1026.
42. Walsh, T. J., C. Gonzalez, C. A. Lyman, S. J. Chanock, and P. A. Pizzo. 1996. Invasive fungal infections in children: recent advances in diagnosis and treatment. *Adv. Ped. Inf. Dis.* 11:187–290.
43. Walsh, T. J., B. De Pauw, E. Anaissle, and P. Martino. 1994. Recent advances in the epidemiology, prevention, and treatment of invasive fungal infections in neutropenic patients. *J. Med. Vet. Mycol.* 32 Suppl 1:33–51.
44. Warnock, D. W. 1995. Fungal complications of transplantation: diagnosis, treatment, and prevention. *J. Antimicrob. Chemother.* 36 Suppl B:73–90.
45. Yamakami, Y., A. Hashimoto, I. Tokimatsu, and M. Nasu. 1996. PCR detection of DNA specific for *Aspergillus* species in serum of patients with invasive aspergillosis. *J. Clin. Microbiol.* 34:2464–2468.
46. Young, R. C., and J. E. Bennett. 1971. Invasive aspergillosis: absence of detectable antibody response. *Amer. Rev. Respir. Dis.* 104:710–716.
47. Zervos, M. J. and J. A. Vasquez. 1996. DNA analysis in the study of fungal infections in the immunocompromised host. *Clin. Lab. Med.* 16:73–88.

## SEQUENCE LISTING

&lt;160&gt; NUMBER OF SEQ ID NOS: 61

<210> SEQ ID NO 1  
<211> LENGTH: 208  
<212> TYPE: DNA  
<213> ORGANISM: *Aspergillus flavus*

&lt;400&gt; SEQUENCE: 1

gctgcccatac aagcacggct tgggtttgg gtcgtcgatcc cctctccggg ggggacgggc	60
ccccaaaggca gcggcggcac cgcgatccat cctcgagcgt atggggcttt gtcacccgct	120
ctgttaggccc ggccggcgct tgccgaacgc aaatcaatct ttttccaggt tgacctcgga	180
tcaaggtaggg ataccgcgtg aacttcaa	208

<210> SEQ ID NO 2  
<211> LENGTH: 364  
<212> TYPE: DNA  
<213> ORGANISM: *Aspergillus fumigatus*

&lt;400&gt; SEQUENCE: 2

aaactttcaa caatggatct ctgggttccg gcatcgatga agaacgcagc gaaatgcgt	60
aactaatgtg aattgcagaa ttcaagtgaat catcgagtct ttgaacgcac attgcgcggc	120
ctggatttcc ggggggcattg cctgtccgag cgtcattgtt gcccattcaag cacggcttgt	180
gtgttgggcc cccgtccccc tctccgggg gacgggccc aaaggcagcg gcggcaccgc	240
gtccggtcct cgagcgtatg gggcttgta cctgctctgt aggccggcc ggcggccagcc	300
gacacccaac ttatattttc taagggttgc acggatcag gttagggatac ccgctgaact	360
taaa	364

-continued

---

<210> SEQ ID NO 3  
<211> LENGTH: 365  
<212> TYPE: DNA  
<213> ORGANISM: Aspergillus niger

<400> SEQUENCE: 3

```

aaactttcaa caatggatct cttggttccg gcatcgatga agaacgcagc gaaatgcgat      60
aactaatgtg aattgcagaa ttcaagtgaat catcgagtct ttgaacgcac attgcgcccc      120
ctggtattcc ggggggcatg cctgtccgag cgtcattgtc gccctcaagc acggcttgcg      180
tgttgggtcg ccgtccccct ctcccggggg acggggccga aaggcagcgg cggcaccgcg      240
tccgatcctc gagcgtatgg ggctttgtca cctgctctgt aggccggcc ggcgcctgcc      300
gacgttatcc aaccatttt ttccaggttg acctcggatc aggttagggat acccgctgaa      360
cttaa                                              365

```

<210> SEQ ID NO 4  
<211> LENGTH: 355  
<212> TYPE: DNA  
<213> ORGANISM: Aspergillus terreus

<400> SEQUENCE: 4

```

aaactttcaa caatggatct cttggttccg gcatcgatga agaacgcagc gaaatgcgat      60
aactaatgtg aattgcagaa ttcaagtgaat catcgagtct ttgaacgcac attgcgcccc      120
ctggtattcc ggggggcat gctgtccga gcgtcattgc tgccctcaag cccggcttgc      180
gtgttgggcc ctgtccccc ggctccggg ggacggggcc gaaaggcagc ggccggaccgc      240
cgtccggtcc tcgagcgtat ggggcttcgt cttccgctcc gtaggccgg ccggcgcccg      300
ccgaacgcat ttatggcaa ttgttttt tttccaggtt gacctcggat caggt               355

```

<210> SEQ ID NO 5  
<211> LENGTH: 365  
<212> TYPE: DNA  
<213> ORGANISM: Aspergillus nidulans

<400> SEQUENCE: 5

```

aaactttcaa caatggatct cttggttccg gcatcgatga agaacgcagc gaactgcgat      60
aagtaatgtg aattgcagaa ttcaagtgaat catcgagtct ttgaacgcac attgcgcccc      120
ctggcattcc ggggggcatg cctgtccgag cgtcattgtc gccctcaagc ccggcttgcg      180
tgttgggtcg tcgtccccc ccccgaaaaa cggggccgaa aggcagcggc ggcaccggc      240
cggtcctcga gcgtatgggg ttggtcacc cgctcgatta gggccggccg ggcgccagcc      300
ggcgtctcca accttatctt tctcaggttg acctcggatc aggttagggat acccgctgaa      360
cttaa                                              365

```

<210> SEQ ID NO 6  
<211> LENGTH: 319  
<212> TYPE: DNA  
<213> ORGANISM: Fusarium solani

<400> SEQUENCE: 6

```

gaaaatgcga taagtaatgt gaattgcaga attcagtgaa tcatcgatc tttgaacgca      60
cattgcgccc gccagtattc tggcgggcat gcctgttcga gcgtcattac aaccctcagg      120
cccccgccgc tggcggtgg gatcggcgga agccccctgc gggcacaacg ccgtccccc      180
aatacagtgg cggtccgcgc gcagcttcca ttgcgttagta gctaacacct cgcaactgga      240

```

-continued

gagcggcgcg gccacgccgt aaaacaccca acttctgaat gttgacacctg aatcaggtag 300  
qaatacccq tqaacttaa 319

<210> SEQ ID NO 7  
<211> LENGTH: 310  
<212> TYPE: DNA  
<213> ORGANISM: Fusarium moniliforme  
  
<400> SEQUENCE: 7

aaatgcgata agtaatgtga attgcaaaaat tcagtgaatc atcgaatctt tgaacgcaca 60  
ttgcgcggcg cagtattctg gcgggcatgc ctgttcgagc gtcatttcaa ccctcaagcc 120  
cccgggttg gtgttgggaa tcggcaagcc cttgcggcaa gccggccccg aaatcttagtg 180  
gcgggtctcgc tgcaagcttcc attgcgtagt agtaaaaaccc tcgcaactgg tacgcggcgc 240  
ggccaagccg ttaaaacccccc aacttctgaa tgttgacctc ggatcaggta ggaataccccg 300  
ctgaacttaa 310

<210> SEQ\_ID NO 8  
<211> LENGTH: 330  
<212> TYPE: DNA  
<213> ORGANISM: Mucor rouxii  
  
<400> SEQUENCE: 8

aaagtgcgat aactagtgtg aattgcata tcagtgaatc atcgagtctt tgaacgcaac 60  
ttgcgcgtcat tggttattcca atgagcacgc ctgtttcagt atcaaaaacaa accctctatc 120  
cagcattttg ttgaatagga atactgagag tctcttgatc tattctgatc tcgaacctct 180  
tgaatgtac aaaggcctga tcttgttaa atgcctgaac ttttttttaa tataaagaga 240  
agctcttgcg gtaaaactgtg ctggggcctc ccaaataata ctcttttaa attgatctg 300  
aaatcaggcg ggattaccgg ctgaacttaa 330

<210> SEQ ID NO 9  
<211> LENGTH: 328  
<212> TYPE: DNA  
<213> ORGANISM: Mucor racemosus  
  
<400> SEQUENCE: 9

aaagtgcgat aactagtgtg aattgcatat tcagtgaatc atcgagtctt tgaacgcaac 60  
ttgcgcctcat tggttattcca atgagcacgc ctgtttcagt atcaaaaacaa accctctatc 120  
caactttgt tgtataggat tattgggggc ctctcgatct gtatacatct tgaatccct 180  
gaaatttact aaggcctgaa cttgtttaaa tgcctgaact ttttttaat ataaaggaaa 240  
gctcttgtaa ttgactttga tggggcctcc caaataaatc tcttttaat ttgatctgaa 300  
atcaggcgaa attaccggct gaacttaa 328

<210> SEQ ID NO 10  
<211> LENGTH: 327  
<212> TYPE: DNA  
<213> ORGANISM: Mucor plumbeus  
  
<400> SEQUENCE: 10  
  
aaagtgcgat aactagtgtg aattgcata tcagtgaatc atcgagtctt tgaacgcaac 60  
ttgcgcgtcat tggttattcca atgagcacgc ctgtttcagt atcaaaaaca accctctatc 120  
caactttgt tgtataggat tattgggggc ctctcgatct gtatacatct taaaaaccctt 180  
qaaatttact aaqqcctqaa cttqtttaat qcctqaactt ttttttaata taaaqqaaaq 240

---

-continued

---

ctcttgataat tgactttgat ggggcctccc aaataaatct ttttaaatt tgatctgaaa	300
tcaggtggaa ttacccgctg aacttaa	327

<210> SEQ ID NO 11	
<211> LENGTH: 322	
<212> TYPE: DNA	
<213> ORGANISM: Mucor indicus	
<400> SEQUENCE: 11	
aaagtgcgat aactagtgtg aattgcataat tcagtgaatc atcgaggctt tgaacgcac	60
ttgcactcaa tggattcca ttgagtcgc ctgtttcagt atcaaaaaca acccttattc	120
aaaattcttt ttttgaatag atatgagtgt agcaaccta caagttgaga cattttaaat	180
aaagtcaggc catatcgtgg attgagtgcc gatactttt taatttgaa aaggtaaagc	240
atgttcatgt ccgccttttgc ggcctccaa ataactttt aaacttgatc tggaaatcagg	300
tgggattacc cgctgaactt aa	322

<210> SEQ ID NO 12	
<211> LENGTH: 330	
<212> TYPE: DNA	
<213> ORGANISM: Mucor circinelloides f.	
<400> SEQUENCE: 12	
aaagtgcgat aactagtgtg aattgcataat tcagtgaatc atcgaggctt tgaacgcac	60
ttgcgctcat tggattcca atgagcacgc ctgtttcagt atcaaaaacaa accctctatc	120
caacattttt gttgaatagg atgactgaga gtctcttgcat ctattctgat ctcgaagctc	180
ttgaaatgta caaaggcctg atcttgcatttgc aatgcctgaa cttttttta atataaagag	240
aagctcttgc ggtaaaactgt gctggggcct cccaaataac acatctttaa atttgatctg	300
aaatcaggtg ggactacccg ctgaacttaa	330

<210> SEQ ID NO 13	
<211> LENGTH: 333	
<212> TYPE: DNA	
<213> ORGANISM: Rhizopus oryzae	
<400> SEQUENCE: 13	
agtgcgataa cttagtgtgaa ttgcataattc agtgaatcat cgagtctttg aacgcagctt	60
gcactctatg gttttctat agagtacgcc tgcttcagta tcatacataaa cccacacata	120
acatttgcattt atgtggatgt gggtcgcac tcgttttat tacagtgagc acctaaaatg	180
tgtgtgattt tctgtctggc ttgctaggca ggaatattac gctggctca ggatctttt	240
ttttgggtcg cccaggaagt aaagtacaag agtataatcc agtaacttcc aaactatgat	300
ctgaagtca gttggattac ccgcgtgaact taa	333

<210> SEQ ID NO 14	
<211> LENGTH: 333	
<212> TYPE: DNA	
<213> ORGANISM: Rhizopus oryzae	
<400> SEQUENCE: 14	
agtgcgataa cttagtgtgaa ttgcataattc agtgaatcat cgagtctttg aacgcagctt	60
gcactctatg gttttctat agagtacgcc tgcttcagta tcatacataaa cccacacata	120
acatttgcattt atgtggatgt gggtcgcac tcgttttat tacagtgagc acctaaaatg	180

**-continued**


---

tgtgtgattt tctgtctggc ttgcttaggca ggaatattac gctggctca ggatctttt	240
ctttggttcg cccaggaagt aaagtacaag agtataatcc agcaacttgc aaactatgtat	300
ctgaagtcag gtgggattac ccgctgaact taa	333

<210> SEQ ID NO 15  
<211> LENGTH: 348  
<212> TYPE: DNA  
<213> ORGANISM: Rhizopus microsporus

<400> SEQUENCE: 15

aaagtgcgat aactagtgtg aattgcataat tcgtgaatca tcgagtcattt gaacgcagct	60
tgcactctat ggatcttcta tagagtacgc ttgcttcagt atcataacca acccacacat	120
aaaatttatt ttatgtggtg atggacaaggc tcggtaaat ttaattatta taccgattgt	180
ctaaaataca gcctctttgt aattttcatt aaattacgaa ctacctagcc atcgtgcattt	240
tttggtccaa ccaaaaaaca tataatctag gggttctgct agccagcaga tattttatg	300
atctttaact atgatctgaa gtcaagtggg actacccgct gaacttaa	348

<210> SEQ ID NO 16  
<211> LENGTH: 349  
<212> TYPE: DNA  
<213> ORGANISM: Rhizopus microsporus

<400> SEQUENCE: 16

aaagtgcgat aactagtgtg aattgcataat tcgtgaatca tcgagtcattt gaacgcagct	60
tgcactctat ggatcttcta tagagtacgc ttgcttcagt atcataacca acccacacat	120
aaaatttatt ttatgtggtg atggacaaggc tcggtaaat ttaattatta taccgattgt	180
ctaaaataca gcctctttgt aattttcatt aaattacgaa ctacctagcc atcgtgcattt	240
tttggtccaa ccaaaaaaca tataatctag gggttctgct agccagcaga tattttatg	300
atctttaacc tatgatctga agtcaagtgg gactacccgc tgaacttaa	349

<210> SEQ ID NO 17  
<211> LENGTH: 361  
<212> TYPE: DNA  
<213> ORGANISM: Rhizopus circinans

<400> SEQUENCE: 17

aaattgcgat aactagtgtg aattgcattt tcagtgaatc atcgagtcattt tgaacgcattc	60
ttgcgcctttt gggattcttc cctagagcac acttgctca gtatcataac aaaaccctca	120
cctaataattt tttttttta aaaaaaaaaat attagagtgg tattgggttc tctttggtaa	180
ttctttgtaa ttataaaagt acccttaaat gtcataaaca ggttagctt agcttgcattt	240
taaagatctt ctttaggatcat cattactttt cgtaaatctt taataggcct gtcacataat	300
tctaccctta aatttcttaa accttgatct gaagtcaagt gggagtagccc gctgaactta	360
a	361

<210> SEQ ID NO 18  
<211> LENGTH: 360  
<212> TYPE: DNA  
<213> ORGANISM: Rhizopus circinans

<400> SEQUENCE: 18

aaattgcgat aactagtgtg aattgcattt tcagtgaatc atcgagtcattt tgaacgcattc	60
ttgcgcctttt gggattcttc cctagagcac acttgctca gtatcataac aaaaccctca	120

---

-continued

---

cctaataattt ttttttaaaa aaaaaaaata tttagagtggt attgggtct ctttgtaat	180
tctttgtaat tataaaaagta cccttaaatg tcataaacag gttagctta gcttcctt	240
aaagatctc ttagggatc attactttc gtaaatctt aataggcctg tcacataatt	300
ctacccttaa atttcttaaa ctttgatctg aagtcaagtg ggagtacccg ctgaacttaa	360

<210> SEQ ID NO 19  
<211> LENGTH: 486  
<212> TYPE: DNA  
<213> ORGANISM: Rhizopus stolonifer

<400> SEQUENCE: 19	
aaagtgcgat aactagtgtg aattgcataat tcagtgaatc atcgagtctt tgaacgcaac	60
ttgcactcta tggtttccg taaagtacgc ttgcttcagt atcataaaga ccccatcctg	120
attattatTT ttttattaaa ataattaatt ttggagataa taaaaatgag gctctttctt	180
ttcttttttt tttttttaaa aaaaaggggg ggaaagggtc ttttaaaatg ggcaaattct	240
gggttttta ctaaacctga actccccca aaaattcaaa aaaaaaaaaa tggttttac	300
caaatttttt tttttttct ctttttgtg tagttaatac tctattaaat ttatTTactt	360
ggtattataa cgattatgca agaaggaga gaacaaagaa taatgaaaga gagttttaa	420
ataaattctt tttcatttt ttcaatcaat gatctgaagt caagtggat taccgctga	480
acttaa	486

<210> SEQ ID NO 20  
<211> LENGTH: 349  
<212> TYPE: DNA  
<213> ORGANISM: Rhizomucor pusillus

<400> SEQUENCE: 20	
aaattgcgaa aagtaatgcg atctgcagcc tttgcgaatc atcgaattct cgaacgcacc	60
ttgcaccctt tggttcatcc attgggtacg tctagttcag tatctttatt aaccctaaa	120
ggtttatttt ttgataaatac tttggatttg cgggtctgat ggattttcat ccgttcaagc	180
tacccgaaca atttgtatgt tggatccct tgatatttcc ttgagggctt gcattggat	240
ctaatttttt accagtgtgc ttcgagatga tcaagtataa aggtcaatca accacaaata	300
aatttcaact atggatctga acttagatgg gattacccgc tgaacttaa	349

<210> SEQ ID NO 21  
<211> LENGTH: 425  
<212> TYPE: DNA  
<213> ORGANISM: Absidia corymbifera

<400> SEQUENCE: 21	
aaagtgcgat aattattgcg acttgcattc atagcgaatc atcgaggctt cgaacgcattc	60
ttgcgcctag tagtaatct actaggcaca gttgtttcag tatctgcac taccaatcag	120
ttcaacttgg ttcttgaac ctaagcgagc tggaaatggg cttgtgttga tggcattcag	180
ttgctgtcat ggccttaat acatTTatgc ctaggcaatt ggcttttagtc atttgcggaa	240
tgttagactct agagtgcctg aggagcaacg acttgggttag tgagttcata attccaagtc	300
aatcagtctc ttcttgaact aggtcttaat ctttatggac tagtgagagg atctaacttg	360
ggtcttcctt taaaacaaac tcacatctag atctgaaatc aactgagatc acccgctgaa	420
cttaa	425

-continued

<210> SEQ ID NO 22  
<211> LENGTH: 399  
<212> TYPE: DNA  
<213> ORGANISM: Absidia corymbifera

<400> SEQUENCE: 22

aaagtgcgat aattattgcg acttgcattc atagtgaatc atcgagttct tgaacgcattc 60  
ttgcgcctag tagtcaatct actaggcaca gttgttttag tatctgcattc caccaatcaa 120  
cttaaccttt tgtgttgagt tggaactggg cttctagttg atggcattta gttgctgtca 180  
tggccttaaa tcaatgtcct aggtgttaga acatctaaca ccggatggaa acttttagagc 240  
gctttaagag cagcttggtt agtgagttca ataattccaa gcattaagtc ttttaatgaa 300  
ctagcttttc tatctatggg acactacttg gagaaatcca agtaaccttt aaactcccat 360  
ttagatctga aatcaactga gaccaccgc tgaacttaa 399

<210> SEQ ID NO 23  
<211> LENGTH: 359  
<212> TYPE: DNA  
<213> ORGANISM: Cunninghamella elegans

<400> SEQUENCE: 23

aaatcgcgat atgtaatgtg actgcctata gtgaatcatc aaatcttga aacgcatttt 60  
gcacccattatg gtattccata aggtacgtct gtttcagttac cactaataaaa tctctctcta 120  
tccttcatatgtga tagaaaaaaa aaaaataatt tttactgggc ccggggaaatc cttttttttt 180  
tttaataaaaa aggaccaatt ttggcccaaa aaaaagggtt gaactttttt taccagatct 240  
tgcatctagt aaaaacctag tcggctttaa tagattttta ttttcttattt agtttatagc 300  
cattcttata tttttttttt tcttggcctg aaatcagatg ggataccgc tgaacttaa 359

<210> SEQ ID NO 24  
<211> LENGTH: 346  
<212> TYPE: DNA  
<213> ORGANISM: Pseudallescheria boydii

<400> SEQUENCE: 24

aaatgcgata agtaatgtaa attgcaaaat tcagtgaatc atcgaatctt tgaaacgcac 60  
attgcgcccc gcagtaatct gccgggcatg cctgtccgag cgtcatttca accctcgaac 120  
ctccgtttcc ttagggaagc ctagggtcgg tgttggggcg ctacggcaag tcctcgcaac 180  
ccccgttaggc cctgaaatac agtggcggtc ccggccgcggc tgccttctgc gtagtaagtc 240  
tctttgcaa gctcgatttgc ggtccggcg gaggcctgcc gtcaaaccac ctaacaactc 300  
cagatggttt gacctcgat caggtagggt tacccgctga acttaa 346

<210> SEQ ID NO 25  
<211> LENGTH: 346  
<212> TYPE: DNA  
<213> ORGANISM: Pseudallescheria boydii

<400> SEQUENCE: 25

gaaatgcgat aagtaatgtg aattgcagaa ttcaatgtat catcgaatct ttgaaacgcac 60  
cattgcgccc ggcagtaatc tgccggcat gcctgtccga gcgtcatttc aaccctcgaac 120  
cctccgtttcc ctcagggaaat ctcagggtcgg gtgttggggcg gctacggcaag gtctcgcaac 180  
ccctccgttag ggcctgaaat acagtggcg tcccgccgcg gttgccttct gcgtagaagt 240  
ctctttqca aqctcqcat qqqtcccqqc qqaqqcctqc cqtcaaaccac cttataactc 300

-continued

caaatggttt gacctcgat caggtagggt taccgcgtga acttaa 346

<210> SEQ ID NO 26

<211> LENGTH: 344

<212> TYPE: DNA

<213> ORGANISM: Scedosporium apiospermum

<400> SEQUENCE: 26

gaaatgcgt aagtaatgtg aattgcagaat ttcaagtgaat catcgaatct ttgaacgcac	60
attgcgcggc gcagtaatctt gccgggcattt cctgtccgag cgtcatttca accctcgaaac	120
ctccgtttcc tcagggaaagc tcagggtcgg tggggggcg ctacggcgag tcttcgcac	180
cctccgttagg ccctgaaata cagtgccgggt cccggccggg ttgccttctg cgttagtaagt	240
ctcttttgcg agctcgattt gggtcccggc ggaggcctgc cgtcaaacca cctataactc	300
caatggttt gacctcgat caggtaggta cccgctgaac tttaa	344

<210> SEQ ID NO 27

<211> LENGTH: 343

<212> TYPE: DNA

<213> ORGANISM: Scedosporium apiospermum

<400> SEQUENCE: 27

aaatgcgata agtaatgtga attgcagaat tcagtgaaatc atcgaatctt tgaacgcaca	60
ttgcgcggc cagtaatctt ccgggcatttgc ctgtccgagc gtcatttca accctcgaaac	120
tccgtttcc cagggaaagc cagggtcgggt gttggggcgcc tacggcgagt ctgcgcacc	180
cctccgttaggc cctgaaatac agtggcggtt ccggccgggt tgcccttctgc gtagtaagtc	240
tcttttgcgaa gctcgattt ggtcccggcg gaggcctgcc gtcaaacccac ctataactcc	300
agatggtttgc acctcgatc aggttaggtac ccgctgaact taa	343

<210> SEQ ID NO 28

<211> LENGTH: 309

<212> TYPE: DNA

<213> ORGANISM: Penicillium notatum

<400> SEQUENCE: 28

aaatgcgata cgtaatgtga attgcaattt cagtgaaatca tcgagttttt gaacgcacat	60
tgcgcggccctt ggtattccgg ggggcatttgc ttgtccgagcg tcatttgcgc cctcaagcac	120
ggcttgtgtt ttggggccccg tcctccgatc ccggggggacg ggcccgaaag gcagcggcg	180
caccgcgtcc ggtcctcgag cgtatggggc tttgtcaccc gctctgttagg cccggccggc	240
gcttgccat caacccaaat ttttatccag gttgacctcg gatcaggttag ggataacccgc	300
tgaacttaa	309

<210> SEQ ID NO 29

<211> LENGTH: 336

<212> TYPE: DNA

<213> ORGANISM: Sporothrix schenckii

<400> SEQUENCE: 29

gaaatgcgt actaatgtga attgcagaat tcagcgaacc atcgaatctt tgaacgcaca	60
ttgcgcggc cagcattcttgc ggggcatttgc ctgtccgagc gtcatttccc ccctcacgc	120
ccccgttgcg cgctgggtttt ggggcgcctt ccgcctggcg gggggcccccc gaaagcgagt	180
ggcggggccctt gtggaaaggctt ccgagcgcag taccgaacgc atgttctccc ctcgcgtccgg	240

**-continued**


---

aggcccccca ggcccccgc cggtaaaaac gcgcatacg cgcaatctt tttacaaggt	300
tgacctcgga tcaggtgagg ataccgctg acttaa	336
<210> SEQ ID NO 30	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Aspergillus flavus	
<400> SEQUENCE: 30	
gcaaataat ctttttcc	18
<210> SEQ ID NO 31	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Aspergillus fumigatus	
<400> SEQUENCE: 31	
gaacgcaaat caatcttt	18
<210> SEQ ID NO 32	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Aspergillus fumigatus	
<400> SEQUENCE: 32	
ccgacaccca tctttatt	18
<210> SEQ ID NO 33	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Aspergillus niger	
<400> SEQUENCE: 33	
gacgttatcc aaccattt	18
<210> SEQ ID NO 34	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Aspergillus terreus	
<400> SEQUENCE: 34	
gcatttattt gcaacttg	18
<210> SEQ ID NO 35	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Aspergillus nidulans	
<400> SEQUENCE: 35	
ggcgtctcca accttatac	18
<210> SEQ ID NO 36	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Mucor rouxii	
<400> SEQUENCE: 36	
gaataggaat actgagag	18
<210> SEQ ID NO 37	
<211> LENGTH: 15	
<212> TYPE: DNA	
<213> ORGANISM: Mucor indicus	

**-continued**

---

<400> SEQUENCE: 37	
gaaacccttg aaatt	15
<210> SEQ ID NO 38	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Mucor indicus	
<400> SEQUENCE: 38	
cgtggattga gtgccgat	18
<210> SEQ ID NO 39	
<211> LENGTH: 21	
<212> TYPE: DNA	
<213> ORGANISM: Mucor circinelloides f.	
<400> SEQUENCE: 39	
aacatttttg tgaataggat g	21
<210> SEQ ID NO 40	
<211> LENGTH: 15	
<212> TYPE: DNA	
<213> ORGANISM: Mucor racemosus	
<400> SEQUENCE: 40	
gaaatccctg aaatt	15
<210> SEQ ID NO 41	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Rhizopus oryzae	
<400> SEQUENCE: 41	
gagtataatc cagyaact	18
<210> SEQ ID NO 42	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Rhizopus circinans	
<400> SEQUENCE: 42	
cttagggtat cattactt	18
<210> SEQ ID NO 43	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Rhizomucor pusillus	
<400> SEQUENCE: 43	
tccttgaggg ctgcatt	18
<210> SEQ ID NO 44	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Rhizopus stolonifer	
<400> SEQUENCE: 44	
cttggattta taacgatt	18
<210> SEQ ID NO 45	
<211> LENGTH: 18	
<212> TYPE: DNA	

-continued

---

<213> ORGANISM: Pseudallescheria boydii

<400> SEQUENCE: 45

aagtctcttt tgcaagct 18

<210> SEQ ID NO 46  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Penicillium notatum

<400> SEQUENCE: 46

gatcaaccca aattttta 18

<210> SEQ ID NO 47  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Penicillium marneffei

<400> SEQUENCE: 47

gggttggtca ccaccata 18

<210> SEQ ID NO 48  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Penicillium marneffei

<400> SEQUENCE: 48

tggtcaccac catattta 18

<210> SEQ ID NO 49  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Fusarium moniliforme

<400> SEQUENCE: 49

tctagtgacg gtctcgct 18

<210> SEQ ID NO 50  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Fusarium oxysporum

<400> SEQUENCE: 50

cgttaattcg cgttcctc 18

<210> SEQ ID NO 51  
<211> LENGTH: 21  
<212> TYPE: DNA  
<213> ORGANISM: Fusarium solani

<400> SEQUENCE: 51

ctaacacctc gcaactggag a 21

<210> SEQ ID NO 52  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Cunninghamella elegans

<400> SEQUENCE: 52

tagtcggctt taatagat 18

<210> SEQ ID NO 53  
<211> LENGTH: 18

-continued

---

<212> TYPE: DNA  
<213> ORGANISM: Cunninghamella elegans

<400> SEQUENCE: 53

tattaagttt atagccat 18

<210> SEQ ID NO 54  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Cunninghamella elegans

<400> SEQUENCE: 54

taagtttata gccattct 18

<210> SEQ ID NO 55  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Absidia corymbifera

<400> SEQUENCE: 55

gttgctgtca tggccta 18

<210> SEQ ID NO 56  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Sporothrix schenckii

<400> SEQUENCE: 56

gacgcgcagc tctttta 18

<210> SEQ ID NO 57  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Rhizopus microsporus

<400> SEQUENCE: 57

catataatct aggggttc 18

<210> SEQ ID NO 58  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Aspergillus sp.

<400> SEQUENCE: 58

cctcgagcgt atggggct 18

<210> SEQ ID NO 59  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Fusarium sp.

<400> SEQUENCE: 59

cccaacttct gaatgttg 18

<210> SEQ ID NO 60  
<211> LENGTH: 18  
<212> TYPE: DNA  
<213> ORGANISM: Mucor sp.

<400> SEQUENCE: 60

mtggggcctc ccaaataa 18

<210> SEQ ID NO 61

-continued

```

<211> LENGTH: 21
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Description of Artificial Sequence: B-58 biotin
probe

<400> SEQUENCE: 61
gaatcatcga rtctttgaac g

```

21

We claim:

1. An isolated nucleic acid probe that consists essentially of 10 to 50 consecutive nucleotides for species-specific identification of *Aspergillus*, wherein the probe selectively hybridizes under stringent conditions to the internal transcribed spacer 2 nucleic acid sequence of one of *Aspergillus flavus* (SEQ ID NO:1), *Aspergillus fumigatus* (SEQ ID NO:2), *Aspergillus niger* (SEQ ID NO:3), *Aspergillus terreus* (SEQ ID NO:4), or *Aspergillus nidulans* (SEQ ID NO:5), but does not selectively hybridize under stringent conditions to the internal transcribed spacer 2 region of any other *Aspergillus* species, nor does it hybridize to the internal transcribed spacer 2 nucleic acid sequence of *Fusarium solani* (SEQ ID NO:6), *Fusarium moniliforme* (SEQ ID NO:7), *Mucor rouxii* (SEQ ID NO:8), *Mucor racemosus* (SEQ ID NO:9), *Mucor plumbeus* (SEQ ID NO:10), *Mucor indicus* (SEQ ID NO:11), *Mucor circinelloides* f. *circinelloides* (SEQ ID NO:12), *Rhizopus oryzae* (SEQ ID NO:13 and NO:14), *Rhizopus microsporus* (SEQ ID NO:15 and 16), *Rhizopus circinans* (SEQ ID NO:17 and 18), *Rhizopus stolonifer* (SEQ ID NO:19), *Rhizomucor pusillus* (SEQ ID NO:20), *Absidia corymbifera* (SEQ ID NO:21 and 22), *Cunninghamella elegans* (SEQ ID NO:23), *Pseudallescheria boydii* (teleomorph of *Scedosporium apiospermum*) (SEQ ID NO:24, 25, 26, and 27), *Penicillium notatum* (SEQ ID NO:28), or *Sporothrix schenckii* (SEQ ID NO:29).

2. The isolated nucleic acid probe of claim 1 wherein the probe selectively hybridizes with an *Aspergillus flavus* nucleic acid of SEQ ID NO:1, or a complementary sequence thereof.

3. The isolated nucleic acid probe of claim 1 wherein the probe selectively hybridizes with an *Aspergillus fumigatus* nucleic acid of SEQ ID NO:2, or a complementary sequence thereof.

4. The isolated nucleic acid probe of claim 1 wherein the probe selectively hybridizes with an *Aspergillus niger* nucleic acid of SEQ ID NO:3, or a complementary sequence thereof.

5. The isolated nucleic acid probe of claim 1 wherein the probe selectively hybridizes with an *Aspergillus terreus* nucleic acid of SEQ ID NO:4, or a complementary sequence thereof.

6. The isolated nucleic acid probe of claim 1 wherein the probe selectively hybridizes with an *Aspergillus nidulans* nucleic acid of SEQ ID NO:5, or a complementary sequence thereof.

7. A method of detecting a species of *Aspergillus flavus* (SEQ ID NO:1), *Aspergillus fumigatus* (SEQ ID NO:2), *Aspergillus niger* SEQ ID NO:3), *Aspergillus terreus* (SEQ ID NO:4), or *Aspergillus nidulans* (SEQ ID NO:5) in a sample comprising

contacting the sample with a nucleic acid probe consisting essentially of 10 to 50 consecutive nucleotides that

selectively hybridizes with a nucleic acid having a sequence as set forth as SEQ ID NO:1, SEQ ID NO:2, SEQ ID NO:3, SEQ ID NO:4, or SEQ ID NO:5, or a complementary sequence thereof;

wherein hybridization of the nucleic acid probe with the sample indicates the detection of the *Aspergillus* species in the sample.

8. The method of claim 7, wherein the probe selectively hybridizes with an *Aspergillus flavus* nucleic acid of SEQ ID NO:1, or a complementary sequence thereof.

9. The method of claim 7, wherein the probe selectively hybridizes with an *Aspergillus fumigatus* nucleic acid of SEQ ID NO:2, or a complementary sequence thereof.

10. The method of claim 7, wherein the probe selectively hybridizes with an *Aspergillus niger* nucleic acid of SEQ ID NO:3, or a complementary sequence thereof.

11. The method of claim 7, wherein the probe selectively hybridizes with an *Aspergillus terreus* nucleic acid of SEQ ID NO:4, or a complementary sequence thereof.

12. The method of claim 7, wherein the probe selectively hybridizes with an *Aspergillus nidulans* nucleic acid of SEQ ID NO:5, or a complementary sequence thereof.

13. An isolated nucleic acid probe for identifying a filamentous fungus wherein the probe consists essentially of a nucleic acid having a sequence as set forth as SEQ ID NO:61, or a complementary sequence thereof, respectively.

14. The isolated nucleic acid probe of claim 1, wherein the probe consists essentially of a nucleotide sequence having a sequence as set forth as SEQ ID NO:30, SEQ ID NO:31, SEQ ID NO:32, SEQ ID NO:33, SEQ ID NO:34, or SEQ ID NO:35.

15. The isolated nucleic acid probe of claim 1, wherein the probe consists essentially of the nucleotide sequence set forth as SEQ ID NO:30.

16. The isolated nucleic acid probe of claim 1, wherein the probe consists essentially of the nucleotide sequence set forth as SEQ ID NO:31.

17. The isolated nucleic acid probe of claim 1, wherein the probe consists essentially of the nucleotide sequence set forth as SEQ ID NO:32.

18. The isolated nucleic acid probe of claim 1, wherein the probe consists essentially of the nucleotide sequence set forth as SEQ ID NO:33.

19. The isolated nucleic acid probe of claim 1, wherein the probe consists essentially of the nucleotide sequence set forth as SEQ ID NO:34.

20. The isolated nucleic acid probe of claim 1, wherein the probe consists essentially of the nucleotide sequence set forth as SEQ ID NO:35.

21. The method of claim 7, wherein the probe consists essentially of a nucleotide sequence as set forth as SEQ ID NO:30, SEQ ID NO:31, SEQ ID NO:32, SEQ ID NO:33, and SEQ ID NO:34.

22. The method of claim 7, wherein the probe consists essentially of the nucleotide sequence set forth as SEQ ID NO:30.

**53**

**23.** The method of claim 7, wherein the probe consists essentially of the nucleotide sequence set forth as SEQ ID NO:31.

**24.** The method of claim 7, wherein the probe consists essentially of the nucleotide sequence set forth as SEQ ID NO:32.

**25.** The method of claim 7, wherein the probe consists essentially of the nucleotide sequence set forth as SEQ ID NO:33.

**26.** The method of claim 7, wherein the probe consists essentially of the nucleotide sequence set forth as SEQ ID NO:34.

5

&lt;/div